

Dec 03, 2019

Version 2

🌐 Loop L1 (odd level) BsaI type IIS cloning into pCk vectors V.2

🔗 Version 1 is forked from [Loop L1 \(odd level\) type IIS cloning - pCk-ye vectors](#)

DOI

dx.doi.org/10.17504/protocols.io.92gh8bw

Eftychis Frangedakis¹, Susana Sauret-Gueto², Anthony West³, Nicola Patron⁴, Marta Marta Tomaselli², Marius Rebmann², Jim Haseloff²

¹University of Cambridge; ²Plant Sciences, University of Cambridge, OpenPlant;


³previously at Earlham Institute, Norwich; ⁴Earlham Institute, Norwich

OpenPlant Project



Susana Sauret-Gueto

Plant Sciences, University of Cambridge, OpenPlant



Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.92gh8bw>

Protocol Citation: Eftychis Frangedakis, Susana Sauret-Gueto, Anthony West, Nicola Patron, Marta Marta Tomaselli, Marius Rebmann, Jim Haseloff 2019. Loop L1 (odd level) Bsal type IIS cloning into pCk vectors. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.92gh8bw>

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: December 03, 2019

Last Modified: December 03, 2019

Protocol Integer ID: 30504

Keywords: iis cloning into pck vectors protocol, dna circuit, pck vectors protocol, loop assembly, cloning, iis cloning, loop l1, bsai, recursive fabrication

Abstract

Protocol based on:

[Pollak B, Cerda A, Delmans M, et al \(2019\) Loop assembly: a simple and open system for recursive fabrication of DNA circuits. New Phytol 222:628–640](#)

<https://doi.org/10.1111/nph.15625>


Materials


MATERIALS

 Bsal - 5,000 units **New England Biolabs Catalog #R0535L**

 T4 DNA Ligase - 20,000 units **New England Biolabs Catalog #M0202S**

 Sterile water


 BSA, molecular biology grade, 20 mg/ml **New England Biolabs Catalog # B9000S**

 10X NEB T4 DNA ligase buffer **New England Biolabs**

Troubleshooting

Loop pCk vectors

1



Loop vectors for nuclear transformation: pCks (A) and for chloroplast transformation pCkchlo (B).


Loop fusion sites in the pCk vectors to assemble different L0 parts into a L1 construct using a pCk vector and BsaI are: A (GGAG) and F (CGCT).

Loop fusion sites in the pCk vectors to assemble different L1 constructs into a L2 construct using a pCs vector and SapI are: a (ATC), b (GCA), d (TAC), e (CAG) and o (GCT).

Left (LB) and right border (RB) repeats from nopaline C58 T-DNA for Agrobacterium-mediated nuclear transformation. KanR: kanamycin bacterial resistance cassette. LacZ: lacZa cassette for blue-white screening of colonies.

Example of assembly of L0 parts into a transcription unit (L1)

2



Loop assembly of multiple L0 parts into a transcription unit (L1) using a pCk plasmid and BsaI.

Protocol for assembly of L0 parts into a transcription unit (L1)

- Determine the concentrations of each DNA plasmid needed (L0 plasmids and pCk acceptor plasmid) by spectrophotometry (Nanodrop).
In the example in step 2, determine concentration of plasmids L0_PROM5, L0_CDS, L0_3TERM and pCk1.
- Prepare aliquots for each plasmid at a concentration of 15 nM for the L0 plasmids and of 7.5 nM for the acceptor pCk vector. With this final concentration, 1 μ L of each plasmid is added to the plasmids mix (see step 6).

To calculate the concentration in ng/ μ L:

- For a final concentration of 15 nM, the concentration in [ng/ μ L] equals N (the length in bp of the plasmid) divided by 110. This is an approximation of the formula:

$$15 \cdot 10^{-9} \text{ mol/L} \times ((607.4 \times N) + 157.9) \text{ g/mol} \times 10^{-6} \text{ L}/\mu\text{L} \times 10^9 \text{ ng/g} = \text{concentration (ng}/\mu\text{L)}$$

- For a final concentration of 7.5 nM, the concentration in [ng/ μ L] equals N divided by 220.

- 5 Prepare Loop assembly Level 1 reaction master mix (MM) according to Table , if four or less number of L0 parts are assembled into a pCk vector (otherwise see step 8)

| Com pon ents | Volu me (μ L) |
|---|--------------------------|
| Steri le wate r | 3 |
| 10x T4 ligas e buff er (NE B) | 1 |
| 1 mg/ mL bovi ne seru m albu min (NE B) | 0.5 |
| T4 DNA ligas e at 400 U/ μ L (NE B) | 0.25 |
| 10 U/ μ L Bsal (NE B) | 0.25 |

| | | |
|--|---------------------|---|
| | Final volu me | 5 |
|--|---------------------|---|

- 6 Prepare plasmids mix by adding in a 0.2 mL tube: 1 μ L of each L0 plasmid , 1 μ L of the pCk vector (see step 4), and sterile water up to 5 μ L. Mix well.

When 4 L0 parts are assembled into a pCk plasmid, the volume of the plasmid mix is 5 μ L, and thus no volume of water is added.
- 7 Add 5 μ L of MM (step 5) to the 5 μ L of plasmids mix (step 6), to a final volume of 10 μ L. Mix well.
- 8 If more than 4 L0 parts are to be assembled into a pCk vector, reduce the water volume in the MM by 1 μ L (step 5) for each extra 1 μ L of DNA part added in the plasmids mix (step 6).
- 9 Place samples in a thermocycler and use the following program:
Assembly: 26 cycles of 37 °C for 3 min and 16 °C for 4 min.
Termination and enzyme denaturation: 50 °C for 5 min and 80 °C for 10 min.
- 10 Transform 20 μ L of chemically competent E. coli cells (transformation efficiency of 1×10^7 transformants/ μ g plasmid DNA) using 2 μ L of the Loop assembly reaction and then plate on LB agar plates containing 50 μ g/mL kanamycin and 40 μ g/mL of X-gal for blue-white screening.
- 11 Incubate overnight at 37 °C.
- 12 Colonies with white color are likely to contain an L1 insert cloned into the pCk vector (In the example in step 2: PROM5:CDS-3TERM)
Blue color colonies will contain undigested pCk vector with LacZ
- 13 Confirm the presence of the correct insert with Sanger sequencing using the primers pC_F (GCAACGCTCTGTCATCGTTAC) and pC_R (GTAACCTAGGACTTGTGCGACATGTC) for pCk vectors, and pC_F and pC_R2 (CAATCTGCTCTGATGCCGCATAGTTAAG) for pCkchlo vectors.