

Jan 29, 2020

Low Quantity single strand CAGE protocol

In 1 collection

DOI

dx.doi.org/10.17504/protocols.io.bbwwkipcw


Hazuki Takahashi¹, Hiromi Nishiyori-Sueki¹, Piero Carninci¹

¹Laboratory for Transcriptome Technology, RIKEN Center for Integrative Medical Sciences, Yokohama, Kanagawa 230-0045, Japan



Hazuki Takahashi

Laboratory for Transcriptome Technology, RIKEN Center for In...



Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.bbwwkipcw>

Protocol Citation: Hazuki Takahashi, Hiromi Nishiyori-Sueki, Piero Carninci 2020. Low Quantity single strand CAGE protocol. [protocols.io https://dx.doi.org/10.17504/protocols.io.bbwwkipcw](https://dx.doi.org/10.17504/protocols.io.bbwwkipcw)

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this protocol in our group and it is working.

Created: January 29, 2020

Last Modified: January 29, 2020

Protocol Integer ID: 32428

Keywords: CAGE (Cap Analysis of Gene Expression), Transcription Starting Sites, Gene promoter, enhancer RNA, capped rna, functional annotation of the mammalian genome, low quantity of rna, mapping at single nucleotide resolution, rna amount, gene expression analysis, single strand cage protocol, single nucleotide resolution, mammalian genome, gene expression, rna, measuring transcript level, nuclear rna, large cage dataset, transcription, single strand, cellular compartments like nuclear rna, transcript level, genome, sequencing platform, avoiding pcr amplification, dna, encyclopedia of dna element, active promoter, pcr amplification, dna element,

Abstract

The Cap Analysis of Gene Expression (CAGE) is a powerful method to identify the Transcription Starting Sites (TSSs) of capped RNAs while simultaneously measuring transcript level. CAGE allows mapping at single nucleotide resolution all active promoters and enhancers. Large CAGE datasets have been over the years by individual labs and consortia, including the Encyclopedia of DNA Elements (ENCODE) and Functional Annotation Of the Mammalian Genome (FANTOM). These dataset constitute open resource for TSS annotations and gene expression analysis.

Here we provide a protocol for the most recent CAGE version of the method, which is called Low Quantity (LQ) single strand (ss) CAGE (LQ-)ssCAGE. The methods allow to profile samples starting from low quantity of RNAs. Samples may typically include cells cultured in small volumes to maximizing large scale perturbations, cellular compartments like nuclear RNAs or samples from developmental stages. The advantage of (LQ-)ssCAGE is to allow working with small samples and simultaneously avoiding PCR amplification: this step help to preserve exquisite quantification of expression. The protocol makes use of simplified steps, by loading single stranded cap selected cDNA. Further, it allows reducing the RNA amount down to 25 ng and to shorten preparation time (less than 3 days) and to prepare libraries parallel, increasing the number of samples that one operator can reasonably handle (up to 96 samples). By proper organization of the flow of libraries preparation, an operator could further scale up the throughput. The obtained libraries are suitable for sequencing by various Illumina sequencing platforms using indexed reads and original barcode identifiers.

Materials

<Materials and Reagents>

Product name	Maker	Cat. No.
Agencourt AMPure XP	BECKMAN COULTER	A63881
Agencourt RNAClean XP	BECKMAN COULTER	A63987
KAPA Library Quantification Kits (500 reactions)- Illumina / ABI Prism®	KAPA BIOSYSTEMS	KK4835
Sodium periodate ACS Reagent Grade, 25 g	MP Biomedicals	0215257725-
5M Sodium Chloride, Molecular Biology Grade, RNase One Ribonuclease,	Promega	V4221
Lithium chloride solution 8 M, for molecular biology, ≥99%	Promega	M4261
Sodium acetate buffer solution, BioXtra, pH 7.0±0.05 (25 °C), for molecular biology, 3 M, 100ML	SIGMA-ALDRICH	L7026-100ML
non-sterile; 0.2 µm filtered,	SIGMA-ALDRICH	S2404-
Sodium periodate ACS Reagent Grade, DNA Ligation Kit < Mighty Mix > ,	SIGMA-ALDRICH	311448-5G
Ribonuclease H (RNase H) (20 ~ 60 U/µL),	TAKARA BIO INC.	6023
10 mM dNTP Mix,	TAKARA BIO INC.	2150A
Dynabeads M-270 Streptavidin,	ThermoFisher SCIENTIFIC	18427013
Ribonuclease H (2 U/µL),	ThermoFisher SCIENTIFIC	65305
RNaseZap RNase Decontamination Solution,	ThermoFisher SCIENTIFIC	18021014
SuperScript III Reverse Transcriptase,	ThermoFisher SCIENTIFIC	AM9780
UltraPure 0.5 M EDTA, pH 8.0,	ThermoFisher SCIENTIFIC	18080044
UltraPure DNase/RNase-Free Distilled Water,	ThermoFisher SCIENTIFIC	15575020
Biotin (Long Arm) Hydrazide,	ThermoFisher SCIENTIFIC	10977015
0.5M EDTA (pH 8.0),	Vector Laboratories	SP-1100
10w/v% Polyoxyethylene(20) Sorbitan Monolaurate Solution,	Wako Pure Chemical Industries, 311-90075	
1 M Tris-HCl (pH 7.0),	Wako Pure Chemical Industries, 161-24801	
1M Tris-HCl (pH 7.5),	Wako Pure Chemical Industries, 311-90411	
1M Tris-HCl (pH 8.5),	Wako Pure Chemical Industries, 316-90221	
3M Sodium Acetate (pH 5.2),	Wako Pure Chemical Industries, 316-90405	
Dimethyl Sulfoxide,	Wako Pure Chemical Industries, 316-90081	
Ethanol (99.5),	Wako Pure Chemical Industries, 046-21981	
2M NaOH,	Wako Pure Chemical Industries, 057-00456	
Hybridization buffer HT1,	Wako Pure Chemical Industries, 194-05631	
	Illumina	

<Equipment>

Product name,	Maker	Cat. No.
Axygen 0.2 mL Polypropylene PCR Tube Strips and Domed Cap Strips, 8 Tubes/Strip, 8 Domed Caps/Strip,	CORNING, AXYGEN	PCR-0208-CP-C

Clear, Nonsterile		
Axygen 1.5 mL Maxymum Recovery Snaplock	CORNING, AXYGEN	MCT-150-L-C
Microcentrifuge Tube, Polypropylene, Clear, Nonsterile,		
250 Tubes/Pack, 10 Packs/Case		
Axygen 16 Well Polypropylene PCR Microplate,	CORNING, AXYGEN	PCR-16-C
Clear, Nonsterile		
Axygen PCR 1 × 8 Strip Domed Caps,	CORNING, AXYGEN	PCR-02CP-C
Fit 0.2 mL PCR Tube Strips, Clear, Nonsterile,		
X-Pierce Sealing Films, Sterile,	EXCEL Scientific, Inc.	XPS-25
10 µL (0.1 - 10 µL) LongReach Barrier tip,	Sorenson BioScience, Inc.	38000T
10 µL (0.1 - 10 µL) Ultra G Barrier tip,	Sorenson BioScience, Inc.	28200T
1000 µL (100-1000 µl) Barrier tip,	Sorenson BioScience, Inc.	14200T
20 µL (1-20 µL) Barrier tip,	Sorenson BioScience, Inc.	14210T
200 µL (1-200 µL) NX Barrier tip,	Sorenson BioScience, Inc.	30550T
MicroAmp Optical 384-Well Reaction Plate with Barcode,	ThermoFisher SCIENTIFIC	4309849
MicroAmp Optical Adhesive Film,	ThermoFisher SCIENTIFIC	4360954
PIPETMAN Classic, P1000,	GILSON	F123602
PIPETMAN Classic, P2,	GILSON	F144801
PIPETMAN Classic, P20,	GILSON	F123600
PIPETMAN Classic, P200,	GILSON	F123601
MJ Research Tetrad PTC-225 Thermal Cycler,	GMI, MJ Research	8252-30-1004
MINI CENTRIFUGE,	HARMONY (LMS)	CFM-1000
PlateSpin II,	KUBOTA	
Proline Plus Mechanical Pipettes, 8-ch, 0.5 - 10 µL,	sartorius	728120
Proline Plus Mechanical Pipettes, 8-ch, 10 - 100 µL,	sartorius	728130
Proline Plus Mechanical Pipettes, 8-ch, 30 - 300 µL,	sartorius	728140
Vortex-Genie 2,	Scientific Industries, Inc.	SI-0286
miVAC DNA,	SP Scientific, Genevac	DNA-10000-G00
miVac rotor for micro plate,	SP Scientific, Genevac	DRS-00000-200
Dynabeads MPC-S (Magnetic Particle Concentrator),	ThermoFisher SCIENTIFIC	A13346
DynaMag-96 Side Skirted Magnet,	ThermoFisher SCIENTIFIC	12027

<RT primer (for ssCAGE)>

name	5' modification	Sequence	3' modification	base
N6+TCT	P	TCTNNNNNN		9

<RT primers containing barcode (for LQ-ssCAGE)>

name	5' modification	Sequence	Base
TCT_7nt_N15_#001P		TCTGATGCTCNNNNNNNNNNNNNNNNNN	25
TCT_7nt_N15_#002	P	TCTTATGAGCNNNNNNNNNNNNNNNNNN	25
TCT_7nt_N15_#003	P	TCTCGACGATNNNNNNNNNNNNNNNNNN	25

TCT_7nt_N15_#004	P	TCTTAGTCACNNNNNNNNNNNNNNNNNN	25
TCT_7nt_N15_#005	P	TCTCGTACTGNNNNNNNNNNNNNNNNNN	25
TCT_7nt_N15_#006	P	TCTTACGTAGNNNNNNNNNNNNNNNNNN	25
TCT_7nt_N15_#007	P	TCTAGACTCGNNNNNNNNNNNNNNNNNN	25
TCT_7nt_N15_#008	P	TCTTCGCTGANNNNNNNNNNNNNNNNNN	25
TCT_7nt_N15_#009	P	TCTCGATCTGNNNNNNNNNNNNNNNNNN	25
TCT_7nt_N15_#010 P		TCTTGTCACGNNNNNNNNNNNNNNNNNN	25
TCT_7nt_N15_#011 P		TCTATGCACTNNNNNNNNNNNNNNNNNN	25
TCT_7nt_N15_#012 P		TCTGACATCGNNNNNNNNNNNNNNNNNN	25
TCT_7nt_N15_#013 P		TCTATCTAGCNNNNNNNNNNNNNNNNNN	25
TCT_7nt_N15_#014 P		TCTTGTCGCANNNNNNNNNNNNNNNNNNN	25
TCT_7nt_N15_#015 P		TCTCAGATCTNNNNNNNNNNNNNNNNNN	25
TCT_7nt_N15_#016 P		TCTTCAGAGCNNNNNNNNNNNNNNNNNN	25
TCT_7nt_N15_#017 P		TCTATACTGCNNNNNNNNNNNNNNNNNN	25
TCT_7nt_N15_#018 P		TCTATCGTGANNNNNNNNNNNNNNNNNN	25
TCT_7nt_N15_#019 P		TCTGAGCTCTNNNNNNNNNNNNNNNNNN	25
TCT_7nt_N15_#020	P	TCTTCTCGAGNNNNNNNNNNNNNNNNNN	25
TCT_7nt_N15_#021 P		TCTCAGTGTCNNNNNNNNNNNNNNNNNN	25
TCT_7nt_N15_#022	P	TCTGCTACGANNNNNNNNNNNNNNNNNNN	25
TCT_7nt_N15_#023	P	TCTAGTACGCNNNNNNNNNNNNNNNNNN	25
TCT_7nt_N15_#024	P	TCTGACTAGCNNNNNNNNNNNNNNNNNN	25
TCT_7nt_N15_#025	P	TCTTACGCTANNNNNNNNNNNNNNNNNNN	25
TCT_7nt_N15_#026	P	TCTGCTATAGNNNNNNNNNNNNNNNNNN	25
TCT_7nt_N15_#027	P	TCTACACTGTNNNNNNNNNNNNNNNNNN	25
TCT_7nt_N15_#028	P	TCTACATGTCNNNNNNNNNNNNNNNNNN	25
TCT_7nt_N15_#029	P	TCTGCAGACTNNNNNNNNNNNNNNNNNN	25
TCT_7nt_N15_#030	P	TCTTAGCGACNNNNNNNNNNNNNNNNNN	25
TCT_7nt_N15_#031 P		TCTACAGTCGNNNNNNNNNNNNNNNNNN	25
TCT_7nt_N15_#032	P	TCTGTATGACNNNNNNNNNNNNNNNNNN	25
TCT_7nt_N15_#033	P	TCTCAGCAGTNNNNNNNNNNNNNNNNNN	25
TCT_7nt_N15_#034	P	TCTACAGCTGNNNNNNNNNNNNNNNNNN	25
TCT_7nt_N15_#035	P	TCTGCTGATANNNNNNNNNNNNNNNNNNN	25
TCT_7nt_N15_#036	P	TCTTGTCGACNNNNNNNNNNNNNNNNNN	25
TCT_7nt_N15_#037	P	TCTTGCGATANNNNNNNNNNNNNNNNNNN	25
TCT_7nt_N15_#038	P	TCTAGTGTCANNNNNNNNNNNNNNNNNNN	25
TCT_7nt_N15_#039	P	TCTATACGTCNNNNNNNNNNNNNNNNNN	25
TCT_7nt_N15_#040	P	TCTCACTATGNNNNNNNNNNNNNNNNNN	25
TCT_7nt_N15_#041 P		TCTTCACAGTNNNNNNNNNNNNNNNNNN	25
TCT_7nt_N15_#042	P	TCTTACGCGTNNNNNNNNNNNNNNNNNN	25
TCT_7nt_N15_#043	P	TCTCGACGTANNNNNNNNNNNNNNNNNNN	25
TCT_7nt_N15_#044	P	TCTAGCATGCNNNNNNNNNNNNNNNNNN	25
TCT_7nt_N15_#045	P	TCTATCGCGTNNNNNNNNNNNNNNNNNN	25
TCT_7nt_N15_#046	P	TCTGATATGCNNNNNNNNNNNNNNNNNN	25

TCT_7nt_N15_#047	P	TCTTAGCATCNNNNNNNNNNNNNNNNNN	25
TCT_7nt_N15_#048	P	TCTGCTCATGNNNNNNNNNNNNNNNNNN	25
TCT_7nt_N15_#049	P	TCTCTGATGCNNNNNNNNNNNNNNNNNN	25
TCT_7nt_N15_#050	P	TCTGCATATGNNNNNNNNNNNNNNNNNN	25
TCT_7nt_N15_#051 P		TCTGTACATNNNNNNNNNNNNNNNNNN	25
TCT_7nt_N15_#052	P	TCTCTAGCATNNNNNNNNNNNNNNNNNN	25
TCT_7nt_N15_#053	P	TCTCACGTGTNNNNNNNNNNNNNNNNNN	25
TCT_7nt_N15_#054	P	TCTCTAGTACNNNNNNNNNNNNNNNNNN	25
TCT_7nt_N15_#055	P	TCTGATCGTANNNNNNNNNNNNNNNNNNN	25
TCT_7nt_N15_#056	P	TCTCTAGTCGNNNNNNNNNNNNNNNNNN	25
TCT_7nt_N15_#057	P	TCTGACACTGNNNNNNNNNNNNNNNNNN	25
TCT_7nt_N15_#058	P	TCTAGCTCGANNNNNNNNNNNNNNNNNNN	25
TCT_7nt_N15_#059	P	TCTGACTGCANNNNNNNNNNNNNNNNNNN	25
TCT_7nt_N15_#060	P	TCTAGACGCTNNNNNNNNNNNNNNNNNN	25
TCT_7nt_N15_#061 P		TCTCAGTCTGNNNNNNNNNNNNNNNNNN	25
TCT_7nt_N15_#062	P	TCTTGACTGCNNNNNNNNNNNNNNNNNN	25
TCT_7nt_N15_#063	P	TCTATCAGCGNNNNNNNNNNNNNNNNNN	25
TCT_7nt_N15_#064	P	TCTCTCAGAGNNNNNNNNNNNNNNNNNN	25
TCT_7nt_N15_#065	P	TCTTATGCACNNNNNNNNNNNNNNNNNN	25
TCT_7nt_N15_#066	P	TCTTGCTACANNNNNNNNNNNNNNNNNNN	25
TCT_7nt_N15_#067	P	TCTATGCGCANNNNNNNNNNNNNNNNNNN	25
TCT_7nt_N15_#068	P	TCTTCGTGACNNNNNNNNNNNNNNNNNN	25
TCT_7nt_N15_#069	P	TCTCGATACGNNNNNNNNNNNNNNNNNN	25
TCT_7nt_N15_#070	P	TCTGCTAGCANNNNNNNNNNNNNNNNNNN	25
TCT_7nt_N15_#071 P		TCTTAGCGTANNNNNNNNNNNNNNNNNNN	25
TCT_7nt_N15_#072	P	TCTCTGACGANNNNNNNNNNNNNNNNNNN	25
TCT_7nt_N15_#073	P	TCTTACGTCANNNNNNNNNNNNNNNNNNN	25
TCT_7nt_N15_#074	P	TCTCATAGCTNNNNNNNNNNNNNNNNNN	25
TCT_7nt_N15_#075	P	TCTTCGATCANNNNNNNNNNNNNNNNNNN	25
TCT_7nt_N15_#076	P	TCTCATATCGNNNNNNNNNNNNNNNNNN	25
TCT_7nt_N15_#077	P	TCTCTCGATANNNNNNNNNNNNNNNNNNN	25
TCT_7nt_N15_#078	P	TCTAGTCAGTNNNNNNNNNNNNNNNNNN	25
TCT_7nt_N15_#079	P	TCTTCGATAGNNNNNNNNNNNNNNNNNN	25
TCT_7nt_N15_#080	P	TCTGCGTATANNNNNNNNNNNNNNNNNNN	25
TCT_7nt_N15_#081 P		TCTATGTCAGNNNNNNNNNNNNNNNNNN	25
TCT_7nt_N15_#082	P	TCTACGTGCGANNNNNNNNNNNNNNNNNNN	25
TCT_7nt_N15_#083	P	TCTCACTAGTNNNNNNNNNNNNNNNNNN	25
TCT_7nt_N15_#084	P	TCTTCACGTGNNNNNNNNNNNNNNNNNN	25
TCT_7nt_N15_#085	P	TCTGCTGTACNNNNNNNNNNNNNNNNNN	25
TCT_7nt_N15_#086	P	TCTTACTGAGNNNNNNNNNNNNNNNNNN	25
TCT_7nt_N15_#087	P	TCTCGATGCANNNNNNNNNNNNNNNNNNN	25
TCT_7nt_N15_#088	P	TCTTGTCAGANNNNNNNNNNNNNNNNNNN	25
TCT_7nt_N15_#089	P	TCTCTCTAGANNNNNNNNNNNNNNNNNNN	25

TCT_7nt_N15_#090	P	TCTCGTACATNNNNNNNNNNNNNNNNNN	25
TCT_7nt_N15_#091 P		TCTCTAGCTGNNNNNNNNNNNNNNNNNN	25
TCT_7nt_N15_#092	P	TCTCGCTATANNNNNNNNNNNNNNNNNN	25
TCT_7nt_N15_#093	P	TCTTAGCTCANNNNNNNNNNNNNNNNNN	25
TCT_7nt_N15_#094	P	TCTGTATAGCNNNNNNNNNNNNNNNNNN	25
TCT_7nt_N15_#095	P	TCTGTGACTCNNNNNNNNNNNNNNNNNN	25
TCT_7nt_N15_#096	P	TCTAGCATCGNNNNNNNNNNNNNNNNNN	25
TCT_7nt_N15_#097	P	TCTTGAGCATNNNNNNNNNNNNNNNNNN	25
TCT_7nt_N15_#098	P	TCTCGTATAGNNNNNNNNNNNNNNNNNN	25
TCT_7nt_N15_#099	P	TCTACTAGCGNNNNNNNNNNNNNNNNNN	25
TCT_7nt_N15_#100 P		TCTACTGATGNNNNNNNNNNNNNNNNNN	25
TCT_7nt_N15_#101 P		TCTGTACATGNNNNNNNNNNNNNNNNNN	25
TCT_7nt_N15_#102 P		TCTTGACTCANNNNNNNNNNNNNNNNNN	25
TCT_7nt_N15_#103 P		TCTGATCTCANNNNNNNNNNNNNNNNNN	25
TCT_7nt_N15_#104 P		TCTCGACTGTNNNNNNNNNNNNNNNNNN	25
TCT_7nt_N15_#105 P		TCTCAGTAGCNNNNNNNNNNNNNNNNNN	25
TCT_7nt_N15_#106 P		TCTACTATGCNNNNNNNNNNNNNNNNNN	25
TCT_7nt_N15_#107 P		TCTACGATCGNNNNNNNNNNNNNNNNNN	25
TCT_7nt_N15_#108 P		TCTGATCACGNNNNNNNNNNNNNNNNNN	25
TCT_7nt_N15_#109 P		TCTTGACATAGNNNNNNNNNNNNNNNNNN	25
TCT_7nt_N15_#110 P		TCTCGTCATANNNNNNNNNNNNNNNNNN	25
TCT_7nt_N15_#111 P		TCTAGCTCAGNNNNNNNNNNNNNNNNNN	25
TCT_7nt_N15_#112 P		TCTGTAGCACNNNNNNNNNNNNNNNNNN	25
TCT_7nt_N15_#113 P		TCTCTATGCGNNNNNNNNNNNNNNNNNN	25
TCT_7nt_N15_#114 P		TCTGACTACTNNNNNNNNNNNNNNNNNN	25
TCT_7nt_N15_#115 P		TCTGATCGCTNNNNNNNNNNNNNNNNNN	25
TCT_7nt_N15_#116 P		TCTCTGATAGNNNNNNNNNNNNNNNNNN	25
TCT_7nt_N15_#117 P		TCTGTGCATANNNNNNNNNNNNNNNNNN	25
TCT_7nt_N15_#118 P		TCTGATGCATNNNNNNNNNNNNNNNNNN	25
TCT_7nt_N15_#119 P		TCTACTGCAGNNNNNNNNNNNNNNNNNN	25
TCT_7nt_N15_#120 P		TCTTGATCACNNNNNNNNNNNNNNNNNN	25

<5' linkers>

name, 5' modification, sequence, 3' modification, base

5'adaptor N6, -, AATGATACGGCGACCACCGAGATCTACACTCTTTCCCTACACGACGCTCTTCCGATCTNNNNNN, P, 64

5'adaptor GN5, -, AATGATACGGCGACCACCGAGATCTACACTCTTTCCCTACACGACGCTCTTCCGATCTGNNNNNN, P, 64

5'adaptor down, P, AGATCGGAAGAGCGTCGTGTAGGGAAAGAGTGTAGATCTCGGTGGTCGCCGTATCATT, P, 58

<3' linkers >

name

sequence, 3' modification, base

3'end adaptor up 01

NNNAGATCGGAAGAGCACACGTCTGAACTCCAGTCACATCACGATCTCGTATGCCGTCTTCTGCTTG, P, 67

3'end adaptor down 01

CAAGCAGAAGACGGCATACGAGATCGTGATGTGACTGGAGTTCAGACGTGTGCTCTTCCGA, - , 61

3'end adaptor up 02

NNNAGATCGGAAGAGCACACGTCTGAACTCCAGTCACCGATGTATCTCGTATGCCGTCTTCTGCTTG, P, 67

3'end adaptor down 02

CAAGCAGAAGACGGCATACGAGATACATCGGTGACTGGAGTTCAGACGTGTGCTCTTCCGA, - , 61

3'end adaptor up 03

NNNAGATCGGAAGAGCACACGTCTGAACTCCAGTCACTTAGGCATCTCGTATGCCGTCTTCTGCTTG, P, 67

3'end adaptor down 03

CAAGCAGAAGACGGCATACGAGATGCCTAAGTGACTGGAGTTCAGACGTGTGCTCTTCCGA, - , 61

3'end adaptor up 04

NNNAGATCGGAAGAGCACACGTCTGAACTCCAGTCACTGACCAATCTCGTATGCCGTCTTCTGCTTG, P, 67

3'end adaptor down 04

CAAGCAGAAGACGGCATACGAGATTGGTCAGTGACTGGAGTTCAGACGTGTGCTCTTCCGA, - , 61

3'end adaptor up 05

NNNAGATCGGAAGAGCACACGTCTGAACTCCAGTCACACAGTGATCTCGTATGCCGTCTTCTGCTTG, P, 67

3'end adaptor down 05

CAAGCAGAAGACGGCATACGAGATCACTGTGTGACTGGAGTTCAGACGTGTGCTCTTCCGA, - , 61

3'end adaptor up 06

NNNAGATCGGAAGAGCACACGTCTGAACTCCAGTCACGCCAATATCTCGTATGCCGTCTTCTGCTTG, P, 67

3'end adaptor down 06

CAAGCAGAAGACGGCATACGAGATATTGGCGTGACTGGAGTTCAGACGTGTGCTCTTCCGA, - , 61

3'end adaptor up 07

NNNAGATCGGAAGAGCACACGTCTGAACTCCAGTCACCAGATCATCTCGTATGCCGTCTTCTGCTTG, P, 67

3'end adaptor down 07

CAAGCAGAAGACGGCATACGAGATGATCTGGTGACTGGAGTTCAGACGTGTGCTCTTCCGA, - , 61

3'end adaptor up 08

NNNAGATCGGAAGAGCACACGTCTGAACTCCAGTCACACTTGAATCTCGTATGCCGTCTTCTGCTTG, P, 67

3'end adaptor down 08

CAAGCAGAAGACGGCATACGAGATTCAAGTGACTGGAGTTCAGACGTGTGCTCTTCCGA, - , 61

3'end adaptor up 09

NNNAGATCGGAAGAGCACACGTCTGAACTCCAGTCACGATCAGATCTCGTATGCCGTCTTCTGCTTG, P, 67

3'end adaptor down 09

CAAGCAGAAGACGGCATACGAGATCTGATCGTGACTGGAGTTCAGACGTGTGCTCTTCCGA, - , 61

3'end adaptor up 10

NNNAGATCGGAAGAGCACACGTCTGAACTCCAGTCACTAGCTTATCTCGTATGCCGTCTTCTGCTTG, P, 67

3'end adaptor down 10

CAAGCAGAAGACGGCATACGAGATAAGCTAGTGACTGGAGTTCAGACGTGTGCTCTTCCGA, - , 61

3'end adaptor up 11

NNNAGATCGGAAGAGCACACGTCTGAACTCCAGTCACGGCTACATCTCGTATGCCGTCTTCTGCTTG, P, 67

3'end adaptor down 11

CAAGCAGAAGACGGCATAACGAGATGTAGCCGTGACTGGAGTTCAGACGTGTGCTCTTCCGA, - , 61

3'end adaptor up 12

NNNAGATCGGAAGAGCACACGTCTGAACTCCAGTCACCTTGTAATCTCGTATGCCGTCTTCTGCTTG, P, 67

3'end adaptor down 12

CAAGCAGAAGACGGCATAACGAGATTACAAGGTGACTGGAGTTCAGACGTGTGCTCTTCCGA, - , 61

3'end adaptor up 13

NNNAGATCGGAAGAGCACACGTCTGAACTCCAGTCACAGTCAACAATCTCGTATGCCGTCTTCTGCTTG, P, 69

3'end adaptor down 13

CAAGCAGAAGACGGCATAACGAGATTGTTGACTGTGACTGGAGTTCAGACGTGTGCTCTTCCGA, - , 63

3'end adaptor up 14

NNNAGATCGGAAGAGCACACGTCTGAACTCCAGTCACAGTTCCGTATCTCGTATGCCGTCTTCTGCTTG, P, 69

3'end adaptor down 14

CAAGCAGAAGACGGCATAACGAGATACGGAAGTGTGACTGGAGTTCAGACGTGTGCTCTTCCGA, - , 63

3'end adaptor up 15

NNNAGATCGGAAGAGCACACGTCTGAACTCCAGTCACATGTCAGAATCTCGTATGCCGTCTTCTGCTTG, P, 69

3'end adaptor down 15

CAAGCAGAAGACGGCATAACGAGATTCTGACATGTGACTGGAGTTCAGACGTGTGCTCTTCCGA, - , 63

3'end adaptor up 16

NNNAGATCGGAAGAGCACACGTCTGAACTCCAGTCACCCGTCCCGATCTCGTATGCCGTCTTCTGCTTG, P, 69

3'end adaptor down 16

CAAGCAGAAGACGGCATAACGAGATCGGGACGGGTGACTGGAGTTCAGACGTGTGCTCTTCCGA, - , 63

3'end adaptor up 18

NNNAGATCGGAAGAGCACACGTCTGAACTCCAGTCACGTCCGCACATCTCGTATGCCGTCTTCTGCTTG, P, 69

3'end adaptor down 18

CAAGCAGAAGACGGCATAACGAGATGTGCGGACGTGACTGGAGTTCAGACGTGTGCTCTTCCGA, - , 63

3'end adaptor up 19

NNNAGATCGGAAGAGCACACGTCTGAACTCCAGTCACGTGAAACGATCTCGTATGCCGTCTTCTGCTTG, P, 69

3'end adaptor down 19

CAAGCAGAAGACGGCATAACGAGATCGTTTCACGTGACTGGAGTTCAGACGTGTGCTCTTCCGA, - , 63

3'end adaptor up 20

NNNAGATCGGAAGAGCACACGTCTGAACTCCAGTCACGTGGCCTTATCTCGTATGCCGTCTTCTGCTTG, P, 69

3'end adaptor down 20

CAAGCAGAAGACGGCATAACGAGATAAGGCCACGTGACTGGAGTTCAGACGTGTGCTCTTCCGA, - , 63

3'end adaptor up 21

NNNAGATCGGAAGAGCACACGTCTGAACTCCAGTCACGTTTCGGAATCTCGTATGCCGTCTTCTGCTTG, P, 69

3'end adaptor down 21

CAAGCAGAAGACGGCATAACGAGATTCCGAAACGTGACTGGAGTTCAGACGTGTGCTCTTCCGA, - , 63

3'end adaptor up 22

NNNAGATCGGAAGAGCACACGTCTGAACTCCAGTCACCGTACGTAATCTCGTATGCCGTCTTCTGCTTG, P, 69

3'end adaptor down 22

CAAGCAGAAGACGGCATAACGAGATTACGTACGGTGACTGGAGTTCAGACGTGTGCTCTTCCGA, - , 63

3'end adaptor up 23

NNNAGATCGGAAGAGCACACGTCTGAACTCCAGTCACGAGTGGATATCTCGTATGCCGTCTTCTGCTTG, P, 69

3'end adaptor down 23

CAAGCAGAAGACGGCATACGAGATATCCACTCGTGACTGGAGTTCAGACGTGTGCTCTTCCGA, - , 63

3'end adaptor up 25

NNNAGATCGGAAGAGCACACGTCTGAACTCCAGTCACACTGATATATCTCGTATGCCGTCTTCTGCTTG, P, 69

3'end adaptor down 25

CAAGCAGAAGACGGCATACGAGATATATCAGTGTGACTGGAGTTCAGACGTGTGCTCTTCCGA, - , 63

3'end adaptor up 27

NNNAGATCGGAAGAGCACACGTCTGAACTCCAGTCACATTCTTTATCTCGTATGCCGTCTTCTGCTTG, P, 69

3'end adaptor down 27

CAAGCAGAAGACGGCATACGAGATAAAGGAATGTGACTGGAGTTCAGACGTGTGCTCTTCCGA, - , 63

<Recipes>

1. 250 mM NaIO₄

Dissolve 1 mg of Sodium periodate (ACS Reagent Grade) in 18.7 µL of RNase/DNase free water and keep in the dark. NOTE: The solution can be aliquoted to 50 µL and stored at 80 °C.

2. 100 mM Biotin (long arm) Hydrazide

Dissolve 50 mg of Biotin (long arm) Hydrazide in 1.345 mL of DMSO 1.345ml.

Biotin is packed in 50mg/tube.

NOTE: The solution can be aliquoted to 50 µL and stored at 80 °C.

3. LiCl buffer

Mix following materials and store at room temperature.

8M Lithium chloride solution	3.64 mL
1 M Tris-HCl (pH 7.5)	0.80 mL
10w/v% Polyoxyethylene (20) Sorbitan Monolaurate Solution	0.40 mL
0.5 M EDTA (pH 8.0)	0.16 mL
RNase, DNase –free water	3.64 mL

4. TE wash buffer

Mix following materials and store at room temperature.

RNase, DNase –free water	9.12 mL
1 M Tris-HCl (pH 7.5)	0.40 mL
10w/v% Polyoxyethylene (20) Sorbitan Monolaurate Solution	0.40 mL
0.5 M EDTA (pH 8.0)	0.08 mL

5. Release buffer

Mix following materials and store at room temperature.

RNase ONE 10x Reaction Buffer	100 μ L
10w/v% Polyoxyethylene (20) Sorbitan Monolaurate Solution	1.0 μ L
RNase, DNase –free water	899 μ L

6. 1 x TE buffer

Mix following materials and store at room temperature.

1 M Tris-HCl (pH 8.0)	500 μ L
0.5 M EDTA (pH 8.0)	100 μ L
RNase, DNase-free water	49.4 mL

7. 0.1 M NaCl/TE buffer

1 M NaCl	500 μ L
1 x TE buffer	4.5 mL

8. 2.5 μ M 5'linker

a. Mix following solution and carry out the annealing reaction to generate 100 μ M GN5 linker

1 mM 5'adaptor GN5	4 μ L
1 mM 5'adaptor down	4 μ L
1 M NaCl	4 μ L
1 x TE buffer	28 μ L

<Annealing condition>

95 °C, 5 min: – 0.1 °C sec⁻¹ down to 83 °C; 5 min at 83 °C; – 0.1 °C sec⁻¹ down to 71 °C; 5 min at 71 °C; – 0.1 °C sec⁻¹ down to 59 °C; 5 min at 59 °C; –0.1 °C sec⁻¹, to 47 °C; 5 min at 47 °C; –0.1 °C sec⁻¹, to 35 °C; 5 min at 35 °C; –0.1 °C s⁻¹ to 23 °C; 5 min at 23 °C; –0.1 °C sec⁻¹ to 11 °C, and then hold at 11 °C.

b. Mix following solutions and repeat the annealing reaction at step **a.** to generate 100 μ M N6 linker

1 mM 5'adaptor N6	1 μ L
1 mM 5'adaptor down	1 μ L
1M NaCl	1 μ L
1 x TE buffer	7 μ L

c. Mix 40 μ L of 100 μ M GN5 linker from step **a.** and 10 μ L of 100 μ M N6 linker from step **b.**

d. Dilute 100 μ M of mixed 5'linkers to 2.5 μ M in 0.1 M NaCl/TE buffer. For instance, add 2.5 μ L of 100 μ M mixed linkers from step **c** to 97.5 μ L of 0.1 M NaCl/TE buffer for 25 samples.

NOTE: 2.5 μ M and 100 M mixed 5'linkers can be stored at -20°C.

9. 2.5 μ M 3' linker



a. Mix following solutions and carry out the annealing reaction at the step **8.a.** to generate 100 μM of 3' linker.

1 mM 3'end adaptor up 1 μL

1 mM 3'end adaptor down 1 μL

1 M NaCl 1 μL

1 x TE buffer 7 μL

b. Dilute 100 μM mixed linker from step **a.** to 2.5 μM in 0.1 M NaCl/TE buffer (see details at the step **8.d.**).


Troubleshooting

Before start

RNAs should be treated in the RNase free solutions, benches and equipments. We highly recommend to perform RNA quality check before starting (see detail at previously published paper [**RNA stability check**](#)).




RT reaction in different well using RT primer with barcode




Step 1.2

Sample mixture



Oxidation/Biotinylation




Step 2-4

RNase ONE




Step 5

Cap-trap



Step 6-8

Release



Step 9-10

5' end adaptor Ligation



Step 11

3' end adaptor Ligation



Step 12

purification



Figure.1 Protocol Flow



1 Reverse Transcription (RT) to generate RNA/cDNA hybrids

1.1 single strand CAGE (ssCAGE)

When the amount of RNAs is more than 1 µg, choose this RT protocol. Here is an example of 5 µg starting RNAs.

When the amount of RNAs is less than 1 µg, go to step **1.2**.

a. Mix 5 µL of RNAs (1 µg/µL), 1 µL of RT primer (N6+TCT) and 4 µL of RNase free water in a tube.

Note

NOTE: Use 1.5 mL tube, 8 strip PCR tube or 16-96 wells plate depending on numbers of samples.

b. Incubate the RNA/primer mix from step **a.** at 65 °C for 5 min and immediately put on ice for 2 min.

c. Add the following components to the solution from step **b.** and carefully mix by pipetting on ice.

5 x First-Strand buffer	4 µL
10 mM dNTPs	1 µL
0.1 M DTT	1 µL
RNase free water	2 µL
SuperScript III Reverse Transcriptase	2 µL

d. Incubate the mixed solution (total 20 µL) at 25 °C for 30 sec, 55 °C for 30 min and keep 4 °C to generate RNA/cDNA hybrids.

e. Add 36 µL (1.8-folds) of RNACleanXP beads to 20 µL of RNA/cDNA hybrids from step **d.** and mix well by pipetting.

①. Incubate at room temperature for 5 min. Spin down and set the tube on magnetic bar for 5 min.

②. Discard the supernatant by pipette aspiration.

③. Wash the beads with 200 µL of 70% EtOH.

④. Set the tube on magnetic bar for 5 min.

⑤. Discard the supernatant by pipette aspiration.

⑥. Repeat step ③ ~ ⑤.

⑦. Discard the 70% EtOH completely by pipette aspiration.

⑧. Add 42 µL of RNase free water and mix by pipetting extensively (more than 60 times) to elute RNA/cDNA hybrids.



- ⑨. Incubate at room temperature for 5 min.
- ⑩. Spin down and set the tube on magnetic bar for 5 min.
- ⑪. Transfer supernatant (~ 40 µL of the cDNA/RNA hybrids) to new tube.

1.2 Low Quantity single strand CAGE (LQ-ssCAGE)



When the amount of starting RNAs is less than 1 µg, choose this protocol.

Note

NOTE: The amount of RNAs in one tube should be 25 ng ~ 100 ng for the RT reaction. After RT, samples have each barcode and can be mixed up to 5 µg. For instance, if the starting RNA amount is 50 ng, 100 samples can be mixed. The number of mixed samples is depending on how much samples you have and how deep the samples should be sequenced. Here is an example of 48 mix samples from each 50 ng starting RNA. **CRITICAL:** The mixed solution should be less than 5 µg. Otherwise the following linkers at step **9** and **11** are not enough for the linker ligation.

- a.** Mix 4 µL of 50 ng RNAs (12.5 ng/µL) and 1 µL of 1.25 mM RT primer containing barcode in a 96 wells plate by pipetting.
- b.** Incubate the RNA/primer mix from step **a.** at 65 °C for 5 min and immediately put on ice for 2 min.
- c.** Mix the following components (enzyme mix) on ice:

5 x First Strand buffer	2.0 µL
0.1 M DTT	0.5 µL
10 mM dNTPs	0.5 µL
RNase free water	1.0 µL
SuperScript III Reverse Transcriptase	1.0 µL
Total volume	5.0 µL

Note

NOTE: Make premix solution for 48 samples x 1.1 samples.

- d.** Add 5 µL of enzyme mix from premix solution at step **c.** to RNA/primer mix from step **b.** and carefully mix 10 times by pipetting on ice.
- e.** Incubate at 25 °C for 30 sec, 50 °C for 30 min and keep 4 °C to generate RNA/cDNA hybrids.
- f.** Mix samples by following steps.

- ①. Transfer each 10 μL of RNA/cDNA hybrids from step **e.** to one of new 1.5 mL tube (total volume is 480 μL).
- ②. Add 15 μL of RNase free water to the first 8 wells at the 96 wells plate from step ①, wash wells by pipetting, and transfer the 15 μL of the solutions to the next 8 wells.
- ③. Wash the 8 wells by pipetting and transfer 15 μL of the solutions to the next 8 wells.
- ④. Repeat step ② three times until the end of 8 wells and transfer all solutions (total volume is 120 μL (15 μL x 8 wells)) to the 1.5 mL tube from step ① (final volume is 600 μL).
- ⑤. Mix 600 μL of the solution from step ④ by vortex, spin down and aliquot 200 μL in new two tubes of 1.5 mL tube (total 200 μL in 3 tubes of 1.5 mL tube).

Note

NOTE: Step **f.** is for the collection of remaining molecules to avoid losing any RNA/cDNA hybrids in the wells. We recommend to aliquot less than 200 μL per each 1.5 mL tube at step **f.** ⑤ due to the volume limitations of next RNAClean XP purification step.

g. Add 300 μL (1.8 -folds) of RNAClean XP beads to the 48 mixed RNA/cDNA hybrids in a 1.5 mL tube from step **f.** ⑤, mix well by pipetting and then elute the mixed RNA/cDNA hybrids by following steps.

- ①. Incubate at room temperature for 10 min, spin down and set the 1.5 mL tube on magnetic bar for 5 min.
- ②. Discard the supernatant by pipette aspiration.
- ③. Wash the beads with 1.2 mL of 70% EtOH.
- ④. Set the 1.5 mL tube on magnetic bar for 5 min.
- ⑤. Discard the supernatant by pipette aspiration.
- ⑥. Repeat step ③ ~ ⑤ twice.
- ⑦. Discard the 70% EtOH completely by pipette aspiration.
- ⑧. Add 100 μL RNase free water and mix by pipetting extensively (more than 60 times) to elute RNA/cDNA hybrids.
- ⑨. Incubate at room temperature for 5 min.
- ⑩. Spin down and set the tube on magnetic bar for 5 min.
- ⑪. Transfer 100 μL of the cDNA/RNA hybrids to new 1.5 mL tube.
- ⑫. Repeat step ⑧ ~ ⑪ twice (final volume is 200 μL).
- ⑬. Concentrate 200 μL of the cDNA/RNA hybrids solutions from step ⑫ to around 40 μL by SpeedVac vacuum concentrator at 37 $^{\circ}\text{C}$, and collect the solutions from 3 of 1.5 mL tubes to 1 of 1.5 mL tube (total volume is around 120 μL) and concentrate to 40 μL in the 1.5 mL tube at 37 $^{\circ}\text{C}$. The timing is around 2 h.

Note

CRITICAL: DO NOT DRY UP the cDNA/RNA hybrids solutions by the concentrating step ⑬. Check the sample volume during the concentration several times and adjust the final volume to 40 μL with RNase/DNase-free when the volume become less than 40 μL .

2 Oxidation to modify diol group of cap structure

a. Mix 40 μL of RNA/cDNA hybrid from step **1.1e.** or **1.2g.**, 2 μL of 1 M NaOAc (pH 4.5) and 2 μL of 250 mM NaIO_4 in a tube by 10 times pipetting.

Note

NOTE: 1.5 mL tube, 8-Strip PCR tube or 16 wells plate can be used depending on how many samples should be prepared.

b. Incubate for 5 min on ice in dark by aluminum foil wrapping.

c. Add 16 μL of 1 M Tris-HCl (pH 8.5) to neutralize the solution and mix well by pipetting.

3 Purification

Add 108 μL (1.8 -folds) of RNACleanXP beads to 60 μL of oxidated RNA/cDNA hybrid from step **2**, mix well by pipetting and then elute the 48 samples mixed RNA/cDNA hybrids by following steps.

①. Incubate at room temperature for 5min.

②. Spin down and set the tube on magnetic bar for 5 min.

③. Discard the supernatant by pipette aspiration.

④. Wash the beads with 200 μL of 70% EtOH

⑤. Discard the 70% EtOH.

⑥. Repeat step ④~ ⑤ twice and discarded 70% EtOH completely.

⑦. Add 42 μL of RNase free water and mix by pipetting extensively (more than 60 times) to elute supernatant.

⑧. Incubate at room temperature for 5min.

⑨. Spin down and set the tube on magnetic bar for 5 min.

⑩. Collect 40 μL of the supernatant to new tube.

4 Biotinylation by the coupling reaction to the oxidized RNA/cDNA hybrids

a. Mix 40 μL of purified oxidized RNA/cDNA hybrid from step **3**, 4 μL of 1M NaOAc (pH 6.3) and 4 μL of 100 mM Biotin (long arm) hydrazide by 10 times pipetting.

b. Incubate for 30 min at 40 °C.

c. Add 86.4 μL of RNACleanXP (1.8 -folds) to 44 μL of solution from step b, and perform purification by following step **3** ① ~ ⑩.

Note

NOTE: Biotinylation happens both 5'end of capped RNAs and 3'end of RNA ([see Figure 2 at published protocol](#)).

5 RNaseONE treatment to digest RNA of RNA/cDNA hybrids

- a.** Add 4.5 μL of 10 x RNaseONE buffer and 0.5 μL of RNaseONE to 40 μL of purified biotinylated RNA/cDNA hybrids from step **4** and mix 10 times by pipetting.
- b.** Incubate for 30 min at 37 °C.
- c.** Add 81 μL of RNACleanXP (1.8 -folds) to the solution from step **b.** and purify by the following step **3** ①~ ⑩.

Note

CRITICAL: This step is critical to digest uncompleted cDNA synthesis to the 5' end of capped RNAs ([See Figure 2 at published protocol](#)).

6 Dynabeads M-270 Streptavidin beads preparation

- a.** Add 30 μL of Dynabeads M-270 Streptavidin to new 1.5ml tube, set on the magnetic bar for 5 min, and discard supernatant.
- b.** Wash the beads with 30 μL of LiCl buffer and discard supernatant.
- c.** Repeat step **b.** for 2 times.
- d.** Resuspend the beads in 35 μL of LiCl buffer.

7 CapTrap reaction

- a.** Add 35 μL of beads from step **6** to 40 μL of RNA/cDNA hybrids from step **5** and mix well by pipetting.
- b.** Incubate for 15 minutes at 37 °C and mix by pipetting at the time of 7 min.
- c.** Spin down and set the plate on the magnetic bar for 2 min and discard the supernatant by the pipette aspiration.
- d.** Add 150 μL of TE wash buffer and mix 60 times by pipetting. Spin down and set the plate on the magnetic bar for 2 min and discard the supernatant.
- e.** Repeat step **d.** for 3 times.
- f.** Add 35 μL of Release buffer to the beads and mix 60 times by pipetting.
- g.** Incubate for 5 min at 95°C for 5min and on ice for 1 min.
- h.** Spin down and set the tube on the magnetic bar for 2 min.
- i.** Transfer the supernatant containing CapTrapped cDNA from step **h.** to new 16 well plate.
- j.** Add 30 μL of Release buffer to the beads and mix well by pipetting, spin down, and set the tube on the magnetic stand for 2 min.
- k.** Transfer 30 μL of supernatant containing CapTrapped cDNA to new tube.

8 RNaseONE and RNaseH reaction to remaining RNAs

- a. Add 2.9 μL of release buffer, 2 μL of RNaseONE and 0.1 μL of RNase H to 30 μL of CapTrapped cDNA from step **7** and mix 10 times by pipetting.
- b. Incubate at 37°C for 30 min.
- c. Add 126 μL of AMPureXP (1.8 -folds) to 33 μL of cDNA from step **b.** and perform purification by following **step 3 ① ~ ⑩.**
- d. Dry up 40 μL of purified CapTrapped cDNA by the SpeedVac concentrator at 37°C for around 75 min.
- e. Add 4 μL of RNase free water to the dried pellet.

9 **5' Single Strand Linker Ligation**

- a. Incubate 4 μL of cDNA at 95 °C for 5 min and put on ice for 2 min.
- b. Incubate 4 μL of 5' linker at 55 °C for 5 min and put on ice for 2 min.
- c. Mix 4 μL of 5' linker, 4 μL of CapTrapped cDNA from step **a.** and 16 μL of Mighty Mix by pipetting (total volume is 24 μL).
- d. Incubate at 16 °C for 16 h (overnight).

10 **Purification to remove excess 5' linkers and linker dimers**



- a. Add 46 μL of RNase free water and 126 μL of AMPureXP (1.8 -folds) to 24 μL of cDNA from step **9**, mix well by pipetting.
- b. Purify the solution by following step **3 ① ~ ⑩.**
- c. Incubate at 95 °C for 5 min and immediately put on ice for 2 min.
- d. Add 48 μL of AMPureXP beads (1.2 -folds) to 40 μL of cDNA and mix well by pipetting.
- e. Repeat purification by following step **3 ① ~ ⑩.**
- f. Dry up 40 μL of cDNA by SpeedVac concentrator at 80°C for around 35 min.
- g. Add 4 μL of RNase free water to the dried pellet.

Note

CRITICAL: It is important to perform heat inactivation and repeat purification step to remove non-ligated excess linkers during the purification step. Otherwise, these excess linkers can disturb library sequencing quality; excess linkers compete with the library to anneal Illumina flow-cells. In addition, 5'linker-3'linker dimers can be created at the next step **11**.

11 **3' Single Strand Linker Ligation**

- a. Incubate 4 μL of cDNA at 95 °C for 5 min and put on ice for 2 min.
- b. Incubate 4 μL of 3' linker at 65 °C for 5 min and put on ice for 2 min.
- c. Mix 4 μL of 3'linker, 4 μL of cDNA from step **a.** and 16 μL of Mighty Mix by pipetting.
- d. Incubate at 30 °C for 4 h.

**Note**

NOTE: When you have more than 2 libraries with same barcode at RT primer and plan to mix libraries at the sequencing step **14**, use different Index at the 3'linker.

12 Purification to remove excess 3' linkers and linker dimers (final library) II

- a.** Purify the solution from step **11** by following step **11a.-f.**
- b.** Dissolve the pellet with 10 μL of RNase, DNase-free water.

13 Library QC II

Quantify the concentration of library by KAPA Library Quantification Kit.

14 Sequencing

- a.** Mix the following components

Library (2250 amol)	X μL
RNase, DNase-free water (19-X)	μL
2N NaOH	1 μL
Total	20 μL

- b.** Incubate at room temperature for 5 min.
- c.** Put the tube on ice and add 20 μL of 1 M Tris-HCl (pH 7.0, pre-chilled) to neutralize.
- d.** Add 110 μL of HT1 buffer.
- e.** Transfer 120 μL of the denatured and diluted library to HiSeq2500.

Sequencing kit

- ssCAGE: 50 bp Single Read sequencing.
- LQ-ssCAGE: 50 bp Paired-End sequencing (see Figure 1 for the library structure).
 - Read 1: to read sequence information of cDNA
 - Read 2: to read sequence information of barcode in RT primer
 - Index 1: to read sequence information of index at the 3'linker