



Jul 13, 2023

LRRK1 expression and purification

DOI

dx.doi.org/10.17504/protocols.io.rm7vzx3b5gx1/v1

Yu Xuan Lin^{1,2}, janice reimer^{1,2}

¹Department of Cellular and Molecular Medicine, University of California, San Diego, La Jolla, CA 92093, USA;

²Aligning Science Across Parkinson's (ASAP) Collaborative Research Network, Chevy Chase, MD, 20815, USA



Robert Fagiewicz

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.rm7vzx3b5gx1/v1>

Protocol Citation: Yu Xuan Lin, janice reimer 2023. LRRK1 expression and purification. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.rm7vzx3b5gx1/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working!

Created: July 10, 2023

Last Modified: July 13, 2023

Protocol Integer ID: 84802

Keywords: purification protein purification protocol, lrrk1 expression, length lrrk1, purification

Funders Acknowledgements:

Aligning Science Across Parkinson's: ASAP

Abstract

Protein purification protocol for full-length LRRK1 as done by Leschziner and Reck-Peterson Labs.

Original protocol by Janice M. Reimer and Yu Xuan Lin.

Materials

Equipment:

Ultracentrifuge, we used one with a Ti70 rotor

Dounce homogenizer

FPLC with SP, His, and S200I columns.

Also make sure Ni-NTA beads and TEV protease are available.

Buffers (suggested volumes listed):

Day 1:

300 mL Lysis Buffer :

50 millimolar (mM) HEPES 7.4

0.5 Molarity (M) NaCl

20 millimolar (mM) Imidazole 8.0

0.5 millimolar (mM) TCEP

5 % volume Glycerol

5 millimolar (mM) MgCl₂

20 micromolar (μM) GDP

Plus protease inhibitors (Protein inhibitor tablets and 0.5 millimolar (mM) Pefabloc)

300 mL Wash Buffer :

50 millimolar (mM) HEPES 7.4

0.5 Molarity (M) NaCl

0.5 millimolar (mM) TCEP

5 % volume Glycerol

5 millimolar (mM) MgCl₂

20 micromolar (μM) GDP

100 mL Elution Buffer :

50 millimolar (mM) HEPES 7.4

0.5 Molarity (M) NaCl

300 millimolar (mM) Imidazole 8.0

0.5 millimolar (mM) TCEP

5 % volume Glycerol

5 millimolar (mM) MgCl₂

20 micromolar (μM) GDP

300 mL Dilution Buffer :



50 millimolar (mM) HEPES 7.4
0.5 millimolar (mM) TCEP
5 % volume Glycerol
5 millimolar (mM) MgCl₂
20 micromolar (μM) GDP

500 mL Buffer A (250 mM NaCl) :
50 millimolar (mM) HEPES 7.4
250 millimolar (mM) NaCl
0.5 millimolar (mM) TCEP
5 % volume Glycerol
5 millimolar (mM) MgCl₂
20 micromolar (μM) GDP

300 mL Buffer B (2.5 M NaCl) :
50 millimolar (mM) HEPES 7.4
2.5 Molarity (M) NaCl
0.5 millimolar (mM) TCEP
5 % volume Glycerol
5 millimolar (mM) MgCl₂
20 micromolar (μM) GDP

Day 2:

500 mL Buffer 0 :
50 millimolar (mM) HEPES 7.4
0.5 Molarity (M) NaCl
0.5 millimolar (mM) TCEP
5 % volume Glycerol
5 millimolar (mM) MgCl₂
20 micromolar (μM) GDP

300 mL Buffer 300 :
50 millimolar (mM) HEPES 7.4
0.5 Molarity (M) NaCl
300 millimolar (mM) Imidazole 8.0
0.5 millimolar (mM) TCEP
5 % volume Glycerol
5 millimolar (mM) MgCl₂
20 micromolar (μM) GDP



300 mL Storage Buffer :
20 millimolar (mM) HEPES 7.4
0.7 Molarity (M) NaCl
0.5 millimolar (mM) TCEP
5 % volume Glycerol
2.5 millimolar (mM) MgCl₂
20 micromolar (μM) GDP

Troubleshooting



His6-Z-TEV-LRRK1 Purification

4d 6h

- 1 N-terminally tagged His6-Z-TEV-LRRK1(FL) was expressed in Sf9 insect cells. Insect cells were infected with baculovirus and grown at 27°C for 3 days. 3d
- 2 Cells were harvested and then resuspended in lysis buffer (50mM HEPES pH 7.4, 500 mM NaCl, 20 mM imidazole, 0.5 mM TCEP, 5% glycerol, 5 mM MgCl₂, 20 μM GDP, 0.5 mM Pefabloc and cOmplete EDTA-free protease inhibitor cocktail (Roche). 1h
- 3 Cells were lysed using a Dounce homogenizer and clarified by centrifugation at 50000 rpm, 4°C, 01:28:00, and Ti70 rotor. 2h
- 4 The supernatant was incubated for 1 hr with Ni-NTA agarose beads (Qiagen) equilibrated in lysis buffer. Beads were applied to a gravity column where they were extensively washed with lysis buffer, wash buffer (50mM HEPES pH 7.4, 500mM NaCl, 0.5mM TCEP, 5% glycerol, 5mM MgCl₂, 20uM GDP), followed by elution in lysis buffer containing 300 mM imidazole (50mM HEPES pH 7.4, 500mM NaCl, 0.5mM TCEP, 5% glycerol, 5mM MgCl₂, 20uM GDP, 300mM imidazole) 2h
- 5 The eluted protein was diluted to 250 mM NaCl using dilution buffer (50mM HEPES pH 7.4, 0.5mM TCEP, 5% glycerol, 5mM MgCl₂, 20uM GDP) and loaded onto a SP Sepharose column (Cytiva) equilibrated in buffer (50 mM HEPES pH 7.4, 250 mM NaCl, 0.5 mM TCEP, 5% glycerol, 5 mM MgCl₂, 20 μM GDP) 1h
- 6 The protein was eluted using a 250 mM to 2.5 M NaCl gradient. Fractions containing LRRK1(FL) were pooled, diluted to ~500 mM NaCl, and incubated with TEV protease overnight at 4°C. 18h
- 7 The protein was concentrated and put directly over a S200 size exclusion column (Cytiva) equilibrated in storage buffer (50 mM HEPES pH 7.4, 150 mM NaCl, 0.5 mM TCEP, 5% glycerol, 5 mM MgCl₂ and 20 μM GDP). 5h
- 8 The protein was concentrated to ~5-6 μM and flash frozen in liquid nitrogen for storage. N-terminally tagged His6-Z-TEV-LRRK120-2015 was purified using the same protocol. 1h

Protocol references

<https://www.biorxiv.org/content/10.1101/2022.11.22.517582v2>