


Aug 04, 2023

🌐 LRRK2 and LAMP1 immunofluorescence staining in various cell lines

DOI

dx.doi.org/10.17504/protocols.io.bp2l6x991lqe/v1



Siyu Chen¹, eva karasmanis¹

¹UC San Diego, HHMI



Siyu Chen

UC San Diego, HHMI

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.bp2l6x991lqe/v1>

Protocol Citation: Siyu Chen, eva karasmanis 2023. LRRK2 and LAMP1 immunofluorescence staining in various cell lines. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.bp2l6x991lqe/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Other

After multiple rounds of optimization, anti-LRRK2 antibody generates significant non-specific signals in the LRRK2 KO cell lines so this antibody may not be reliable for Immunofluorescence staining experiments. The protocol can be used for general IF experiments.

Created: August 03, 2023

Last Modified: May 31, 2024

Protocol Integer ID: 85948

Keywords: IF, LRRK2, RAW264.7, LAMP1, ASAPCRN, Irrk2 ko cell line, lamp1 immunofluorescence, reliable for immunofluorescence, immunofluorescence, various cell lines this protocol, antibody, staining experiment, Irrk2, various cell line, cell line

Disclaimer

DISCLAIMER – FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to **protocols.io** is not peer reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with **protocols.io**, can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.

Abstract

This protocol is being used to test the antibody

⊗ Recombinant Anti-LRRK2 antibody [MJFF2 (c41-2)] (ab133474) **Abcam** , as well as

⊗ Anti-LAMP1 antibody [1D4B] (ab25245) **Abcam** . Please note that after multiple rounds of optimization, anti-LRRK2 antibody generates significant non-specific signals in the LRRK2 KO cell lines so this antibody may not be reliable for Immunofluorescence staining experiments. The protocol can be used for general IF experiments.

Guidelines

GDB buffer was used as blocking buffer in the optimized protocol.

This protocol uses ibidi 8 Well Chamber μ -Slides and will coat them with Fibronectin. Standard volume would be 300 μ L for each well. Volumes will need altering for wells, other plates and slides.

All steps are performed at room temperature (RT) on the lab bench, except for methanol permeabilization.



Materials

- (0.1M) NaPi pH 7.4
3.1 g of NaH₂PO₄•H₂O
0.9 g of Na₂HPO₄ (anhydrous)
distilled H₂O to make a volume of 1 L
The pH of the final solution will be 7.4. This buffer can be stored for up to 1 mo at 4°C
- GDB buffer
30mM NaPi (sodium phosphate) pH 7.4
0.45mM NaCl
0.2% porcine (or fish) gelatin
In ddH₂O
- Kim wipes
- Ethanol (100%, stored in dark chemicals cupboard)
- Water (double-deionised H₂O from Milli-Q, "MQ-H₂O")

Protocol materials

☒ Recombinant Anti-LRRK2 antibody [MJFF2 (c41-2)] (ab133474) **Abcam**

☒ Anti-LAMP1 antibody [1D4B] (ab25245) **Abcam**

☒ Chloroquine diphosphate salt **Merck MilliporeSigma (Sigma-Aldrich) Catalog #C6628-25G**

☒ Leu-Leu methyl ester hydrobromide **Merck MilliporeSigma (Sigma-Aldrich) Catalog #L7393-500MG**

☒ Ammonium chloride **Merck MilliporeSigma (Sigma-Aldrich) Catalog #254134**

☒ Triton X-100 **Merck MilliporeSigma (Sigma-Aldrich) Catalog #X100**

☒ Recombinant Anti-LRRK2 antibody [MJFF2 (c41-2)] (ab133474) **Abcam**

☒ Anti-LAMP1 antibody [1D4B] (ab25245) **Abcam**

☒ Gelatin **Merck MilliporeSigma (Sigma-Aldrich) Catalog #G2500**

☒ Goat Anti-Rabbit IgG H&L (Alexa Fluor® 647) (ab150079) **Abcam**

☒ Goat Anti-Rat IgG H&L (Alexa Fluor® 488) (ab150157) **Abcam**

☒ DAPI **Thermo Fisher Scientific Catalog #D1306**

☒ FluorSave™ Reagent **Merck MilliporeSigma (Sigma-Aldrich) Catalog #345789**

Troubleshooting

Safety warnings

! N/A



Before start

Important note:




While LAMP1 immunofluorescence staining robustly gives expected results, after multiple rounds of optimization, anti-LRRK2 antibody generates significant non-specific signals in the LRRK2 KO cell lines in our hand. This is a warning that this antibody may not be reliable for Immunofluorescence staining experiments. However, this protocol can be used widely for IF experiments.

Please check [This google sheet](#) to learn more about the cell lines and conditions tested for the antibody



Day 0: : Seed cells

1h

- 1 Add  300 μ L of  11 μ g/mL fibronectin into each ibidi well. Incubate at RT for  01:00:00 .



1h

- 2 Rinse the wells with PBS for 3 times

- 3 Make GDB buffer if necessary

	Stock	Amount needed	C	Final conc.	E
	(0.1M) NaPi pH 7.4	15	mL	30	mM
	(5M) NaCl	0.0045	mL	0.45	mM
	Gelatin	0.1	g	0.2	%
	H2O	34.9955	mL		
	Total	50	mL		

Recipe to make GDB buffer

- 4 Seed adherent cells to 40-80% confluency in each well. Incubate at  37 °C  Overnight to get optimal seeding.

Note

For RAW264.7 cells, 6×10^4 in 300 μ L is a good starting point. Less would be needed for other typical cell lines since macrophage cells are smaller.
It is suggested to start with two different cell concentrations for the first time. Please refer to [this page](#) for more information

Day 1: Drug treatment

3h

- 5 Apply any drug treatments or controls and note time of additions before proceeding with fixing. As an example,

3h



Chloroquine diphosphate salt **Merck MilliporeSigma (Sigma-Aldrich) Catalog #C6628-25G**



and



Leu-Leu methyl ester hydrobromide **Merck MilliporeSigma (Sigma-Aldrich) Catalog #L7393-500MG**

can be added at desired concentrations for ⌚ 03:00:00

	Drugs	MW - g/mol	mM	weight in 1 mL	E
	LLOME	339.27	1000	339.27	mg/1mL
	CQ	515.86	100	51.586	mg/1mL

Drug stock recipe and concentrations

In each well, add 1 in 1000 ([M] 1 millimolar (mM) for LLOME and

[M] 0.1 millimolar (mM) for CQ)

Staining

25m

- Put 100% ethanol on ice before proceeding with next steps.
- Bring GDB buffer to 🌡 Room temperature . Prepare and prewarm fixation buffer (4% sucrose, 3% PFA in 1xPBS) at 🌡 37 °C

Note

3% PFA is preferred from 4% as it reduces autofluorescence

Need 🧪 375 µL



Pierce™ 16% Formaldehyde (w/v) Methanol-free **Thermo Fisher Scientific Catalog #28906**

and 🧪 0.08 g



Sucrose **Merck MilliporeSigma (Sigma-Aldrich) Catalog #S0389** , dissolve with



2 mL PBS.

- Get cells, aspirate media and immediately add prewarmed fixation buffer. Incubate for ⌚ 00:10:00 at 🌡 37 °C
- Aspirate PFA, rinse 2x with PBS and wash two times with PBS for 5 minutes each, ⌚ 00:10:00 in total

10m

10m



- 10 Quenching: Only necessary when staining the day of fixation. Incubate 3 times of 10 minutes (⌚ 00:30:00 in total) using 0.4%
⊗ Ammonium chloride Merck MilliporeSigma (Sigma-Aldrich) Catalog #254134 ([M] 75 millimolar (mM)) 30m
- 11 Optional step: For half of the samples, choose to add another permeabilisation step with ⌚ 300 µL 100% ethanol at -20°C for ⌚ 00:20:00 . Leave the rest at ⌚ Room temperature in GDB buffer. 20m
- Note
- This step helps to increase contrast when imaging
⊗ Recombinant Anti-LRRK2 antibody [MJFF2 (c41-2)] (ab133474) Abcam , but doesn't help with the issue of non-specific signals when LRRK2 KO cell lines are used.
- Note
- ethanol/methanol permeabilisation at ⌚ -20 °C will disrupt the microtubule staining. Apply this step with caution when other antibodies are used.
- 12 Aspirate ethanol, wash 2x with PBS.
- 13 Prepare GDB + 0.05%
⊗ Triton X-100 Merck MilliporeSigma (Sigma-Aldrich) Catalog #X100 freshly. Add ⌚ 300 µL for ⌚ 00:10:00 .
⌚ 5 µL in ⌚ 10 mL GDB buffer. Can be stored at ⌚ 4 °C for a couple of days. 10m
- 14 Aspirate Triton X-100, wash 2x with GDB.
- 15 Block cells with ⌚ 300 µL GDB for ⌚ 00:30:00 at ⌚ Room temperature 30m
- 16 During blocking, prepare final concentration of [M] 1 µg/mL for both
⊗ Recombinant Anti-LRRK2 antibody [MJFF2 (c41-2)] (ab133474) Abcam and 10m



🧬 Anti-LAMP1 antibody [1D4B] (ab25245) **Abcam** solution with GDB and

🌀 50000 rpm, 4°C, 00:10:00 - each ibidi dish well requires 🧪 150 µL

Note

🧪 150 µL is the minimum required amount. 🧪 200 µL would be sufficient and optimal to cover the entire well

17 Take most of the supernatant and place in new tube. Mix to get even concentration (due to 🧬 Gelatin **Merck MilliporeSigma (Sigma-Aldrich) Catalog #G2500** in there, there will be a small clear precipitate)

18 Aspirate blocking solution and incubate cells with 🧪 150 µL of primary antibody solution ⌚ Overnight at 🌡 4 °C on a table-top shaker

Note

Primary antibody incubation can be as short as 2 hours without affecting the final outcome

Note

When incubating overnight, consider wrapping up the dish with parafilm and/or put the dish in a humidity chamber to prevent the well from drying up

Day 2

15m

19 Bring GDB to 🌡 Room temperature . Aspirate antibody solution and rinse 2x with GDB.

20 Wash 3× 5min 🌡 Room temperature RT with GDB, ⌚ 00:15:00 in total

15m

21 Prepare 1:500 Alexa-fluor conjugated secondary antibody solution with GDB - each ibidi well requires 🧪 200 µL .



When LRRK2 and LAMP1 are co-stained,

Goat Anti-Rabbit IgG H&L (Alexa Fluor® 647) (ab150079) **Abcam** and

Goat Anti-Rat IgG H&L (Alexa Fluor® 488) (ab150157) **Abcam** are used at final

concentration of **[M] 4 ug/mL** . Lower concentration did not help with the issue of non-specific LRRK2 signal.

22 Incubate cells with 200 µL of secondary antibody solution for 01:30:00 at Room temperature and protect from light with an ice box.

1h 30m

23 15 minutes before the incubation is finished, add DAPI **Thermo Fisher Scientific Catalog #D1306** to a final concentration of **[M] 1 ug/mL** and incubate until the last step finishes.

24 Rinse cells with 1xPBS for 5 times

25 Apply 2-4 drops of FluorSave™ Reagent **Merck MilliporeSigma (Sigma-Aldrich) Catalog #345789** hard mounting media in each well and swirl to make sure the bottom is fully covered.

26 Allow to air-dry for 00:10:00 at Room temperature .

10m

27 Image within 48 h of mounting or the sample will begin to deteriorate (bright debris impeding imaging) and visibly autofluoresce in red.