

Oct 26, 2023

## LRRK2 cloning, plasmid construction, and mutagenesis

DOI

[dx.doi.org/10.17504/protocols.io.kxygx35ddg8j/v1](https://dx.doi.org/10.17504/protocols.io.kxygx35ddg8j/v1)

David Snead<sup>1,2</sup>, Yu Xuan Lin<sup>1,2</sup>

<sup>1</sup>Department of Cellular and Molecular Medicine, University of California, San Diego, La Jolla, CA 92093;

<sup>2</sup>Aligning Science Across Parkinson's (ASAP) Collaborative Research Network, Chevy Chase, MD, 20815, USA



**Amalia Villagran Suarez**

Department of Cellular and Molecular Medicine, University of...

### Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.kxygx35ddg8j/v1>

**Protocol Citation:** David Snead, Yu Xuan Lin 2023. LRRK2 cloning, plasmid construction, and mutagenesis. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.kxygx35ddg8j/v1>

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

We use this protocol and it's working

**Created:** October 24, 2023

**Last Modified:** October 26, 2023

**Protocol Integer ID:** 89780

**Keywords:** mutagenesis of lrrk2, lrrk2 cloning, protocol for cloning, mutagenesi, cloning, lrrk2, plasmid construction, rckw

**Funders Acknowledgements:**

**ASAP**

Grant ID: ASAP-000519

**MJFF**

Grant ID: 18321

## Abstract

Protocol for Cloning, plasmid construction, and mutagenesis of LRRK2 and LRRK2-RCKW as done by Leschziner and Reck-Peterson Labs.

Original protocol by David Snead and Yu Xuan Lin.

## Materials

### Materials












- Q5 Site-Directed Mutagenesis (NEB)
- DH5 $\alpha$  competent cells.
- QIAprep Spin Miniprep Kit (Qiagen)

### Equipment

- Thermocycler

## Troubleshooting

## Cloning, plasmid construction, and mutagenesis

- 1 The DNA coding for LRRK2-RCKW residues 1327 to 2527 (taken from Mammalian Gene Collection) was PCR-amplified using the forward primer TACTTCCAATCCATGAAAA491AGGCTGTGCCTTATAACCGA and the reverse primer TATCCACCTTTACTGTCACTCAACAGATGTTTCGTCTCATTTTTTCA.  
The DNA coding for LRRK2 was codon-optimized for *Spodoptera frugiper* (Sf9) cells and synthesized by Epoch Life Science.
  - 2 The DNA for either LRRK2 and LRRK2-RCKW, containing a N-terminal His6-Z-tag and TEV protease cleavage site, was cloned into a pFB-6HZB vector (**SGC**) by ligation-independent cloning, RRID:Addgene\_53641
  - 3 LRRK2 variants were generated using Q5 Site-Directed Mutagenesis Kit (NEB).
- 3.1 Primers for G2019S mutant:  
Forward: gccaaagatcgctgactacagcattgccagctactgttgc  
Reversed: gcaacagctactgggcaatgctgtagtcagcgatcttggc
- Primers for I2020T mutant:  
Forward: ccaagatcgctgactacggaactgccagctact  
Reversed: agtactgggcagttccgtagtcagcgatcttgg
- Using the following thermocycler conditions.
- |  |   |
|--|---|
|  98 °C    |  00:00:30              |
| (  98 °C  |  00:00:10              |
|  50-72 °C |  00:00:30              |
|  72 °C    |  00:13:00              |
|  72 °C    |  00:02:00 ) X25 cycles |
|  4 °C     | Hold.   |
- 4 PCR products were subject to a KLD reaction (NEB), followed by a transformation into DH5α competent cells and plated on LB with antibiotics.
  - 5 Colonies were picked and grown overnight for DNA extraction using Qiagen miniprep kit.
  - 6 Extracted DNA plasmids were submitted for sequencing.



The resulting plasmids were utilized for the generation of recombinant Baculoviruses according to the Bac-to-Bac expression system protocol (Invitrogen).

## Protocol references

<https://www.neb.com/en-us/protocols/2013/01/26/q5-site-directed-mutagenesis-kit-quick-protocol-e0554>