

Oct 23, 2023

LRRK2-RCKW: MLI-2: E11 DARPin cryo-EM sample preparation

DOI

dx.doi.org/10.17504/protocols.io.q26g7p17qgwz/v1

Marta Sanz Murillo¹, Andres Leschziner¹

¹University of California, San Diego



Marta Sanz Murillo

University of California, San Diego

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.q26g7p17qgwz/v1>

Protocol Citation: Marta Sanz Murillo, Andres Leschziner 2023. LRRK2-RCKW: MLI-2: E11 DARPin cryo-EM sample preparation. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.q26g7p17qgwz/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: September 05, 2023

Last Modified: May 31, 2024

Protocol Integer ID: 87407

Keywords: ASAPCRN, em sample preparation protocol, lrrk2, em grid preparation, sample preparation, darpin cryo, cryo, rckw, preparation

Funders Acknowledgements:

Aligning Science Across Parkinson's: ASAP

Grant ID: ASAP-000519

Abstract

Protocol used to prepare LRRK2-RCKW: DARPin:MLi-2 complex and cryo-EM grid preparation.

Troubleshooting

Before start

Make



Protein purification and buffer exchange

- 1 His6-Z-TEV-LRRK2-RCKW was expressed and purified as described in a previous protocol.

Protocol

NAME

LRRK1 expression and purification

CREATED BY

Robert Fagiewicz

Preview

- 2 Prepare LRRK2 buffer exchange. Keep it at 4°C.
20 millimolar (mM) HEPES pH=7.4
150 millimolar (mM) NaCl
2.5 millimolar (mM) MgCl₂
20 micromolar (μM) GDP
0.5 millimolar (mM) TCEP
- 3 Spin down purified LRRK2-RCKW (10000 rcf, 4°C, 10 minutes). Leave protein on ice afterward.
For the best result, keep protein on ice and reduce the amount of time between spinning and freezing cryo-EM samples.
- 4 Exchange buffer using a spin desalting column (Zeba™ Spin Desalting Columns, 7K MWCO (Catalog number: 89877)).
- 5 Spin down again the exchange buffer LRRK2-RCKW (10000 rcf, 4°C, 10 minutes) and measure the concentration. Leave protein on ice afterward

Expected result

The initial concentration range was 20-40 μM. The final concentration might be half of the initial one. Final volume might be 13-16 μL.

- 5.1 Thaw E11 DARPIn and spin it down. Measure its concentration.



- 6 Dilute MLI-2 stock (diluted in 100% DMSO) to a desired concentration
- 7 Based on LRRK2-RCKW concentration, add the necessary volume to get a proportional ratio
LRRK2:DARPin:MLi-2 1:1.25:3 and dilute to a final 10 micromolar (μM) LRRK2-RCKW concentration using exchange buffer (150 mM NaCl).
- 8 Incubate 10 minutes at RT. Afterward, keep it on ice until grid preparation.

E11 DARPin purification

- 9 E11 DARPin purified as described in the next protocol

cryo-EM sample preparation

- 10 We used UltraAuFoil Holey Gold 2/2 200 mesh grids and plasma cleaning them in the Solarus II (Gatan) using the QuantiFoil Au preset.
- 11 Dilute the sample to the desired concentration using the LRRK2 exchange buffer. We used 6 micromolar (μM).
- 12 Apply 3 to 3.5 microliters (μl) of sample and plunge freeze. We used a Vitrobot (FEI) to blot away excess sample and plunge freeze in ethane liquid. (In our case, we use 4 seconds as a time blot as 20 sec as a wait time and 4 as a blot force, but these parameters are slightly different from one Vitrobot to another. I would try with the Vitrobot parameters already tested in your machine first).
- 13 Store grids in liquid nitrogen until ready for imaging