

Jun 10, 2024

## LungMAP2 URM C cellular indexing of transcriptomes and epitopes by sequencing (CITE-seq) protocol

 Forked from [SenNet URM C 10X Genomics Single-Nucleus RNA-Sequencing for Transcriptomic Profiling](#)

DOI

[dx.doi.org/10.17504/protocols.io.81wgbzzyqgpk/v1](https://dx.doi.org/10.17504/protocols.io.81wgbzzyqgpk/v1)

Jeffrey Malik<sup>1</sup>, Gloria S Pryhuber<sup>1</sup>, John Ashton<sup>1</sup>, Ravi Misra<sup>1</sup>, Cameron Baker<sup>1</sup>, Gautam Bandyopadhyay<sup>1</sup>,  
Heidie Huyck<sup>1</sup>

<sup>1</sup>University of Rochester Medical Center

Jeffrey Malik: Genomics Research Center University of Rochester

Gloria S Pryhuber: Department of Pediatrics

John Ashton: Genomics Research Center University of Rochester

Ravi Misra: Department of Pediatrics

Cameron Baker: University of Rochester Medical Center

Gautam Bandyopadhyay: Department of Pediatrics

Heidie Huyck: Department of Pediatrics

### LungMap2 Consortium

Tech. support email: [lungmap2dcc@gmail.com](mailto:lungmap2dcc@gmail.com)

[Click here to message tech. support](#)



Ravi Misra

URMC

## Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account



**DOI:** <https://dx.doi.org/10.17504/protocols.io.81wgbzzyqgpk/v1>

**Protocol Citation:** Jeffrey Malik, Gloria S Pryhuber, John Ashton, Ravi Misra, Cameron Baker, Gautam Bandyopadhyay, Heidie Huyck 2024. LungMAP2 URM C cellular indexing of transcriptomes and epitopes by sequencing (CITE-seq) protocol.

**protocols.io** <https://dx.doi.org/10.17504/protocols.io.81wgbzzyqgpk/v1>

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** May 30, 2024

**Last Modified:** June 10, 2024

**Protocol Integer ID:** 100933

**Keywords:** RNA Sequencing, 10X Genomics, Cellular Indexing of Transcriptomes and Epitopes by Sequencing, LungMAP, lungmap2 urmc cellular indexing of transcriptome, dissociated lung cells for cellular indexing, 10x genomics single cell, epitopes by sequencing, single cells with adt, single cell citeseq data, lungmap2 urmc, protocol the lungmap program, lungmap program, generation of single cell suspension, sequencing data, single cell suspension, dissociated lung cell, lung research community, staining single cell, urmc htc lung tissue digestion sop, single cell, transcriptome, utilized single cells suspension, single cells suspension, sequencing platform, sequencing, cellular indexing, sequencing step, tagged antibody, 10x single cell, lung community, lung, rna

**Funders Acknowledgements:**

**National Institute of Health**

**Grant ID:** U01HL14886



## Abstract

The LungMAP program has been working for over a decade to provide the lung community with a comprehensive resource for the lung research community, including reference data sets and protocols [1]. The University of Rochester has served as the Human Tissue Core and worked extensively with Research Center lab groups. We have established protocols for procuring human tissue including, but not limited to lung, and processing them for a variety of downstream analysis modalities including the generation of single cell suspensions [2, 3, 4]. We have successfully utilized single cell suspensions to generate single cell CITESeq data utilizing the 10X Genomics Single Cell 3' (v3.1) RNA sequencing platform [5-7]. This protocol describes the steps necessary for preparing generating sequencing data following staining single cells with ADT-tagged antibodies (including library preparation, sequencing steps, and mapping steps).

## References

1. <https://www.lungmap.net/about-lungmap/lungmap-objectives/>
2. 602.2 Donor Acceptance Criteria for URMHC HTC SenNet Inclusion [dx.doi.org/10.17504/protocols.io.biz7kf9n](https://doi.org/10.17504/protocols.io.biz7kf9n).
3. 603.3 & 604.5\_URMC\_HTC\_Lung and Lobe Processing for SenNet  
[dx.doi.org/10.17504/protocols.io.n2bvj395plk5/v1](https://doi.org/10.17504/protocols.io.n2bvj395plk5/v1)
4. 702.B.3 URMHC HTC Lung Tissue Digestion SOP + Worksheet 052820  
[dx.doi.org/10.17504/protocols.io.biz5kf86](https://doi.org/10.17504/protocols.io.biz5kf86)
5. 709.1 Staining of Dissociated Lung Cells for Cellular Indexing of Transcriptomes and Epitopes by Sequencing (CITE sequencing) V.2 [dx.doi.org/10.17504/protocols.io.e6nvw1jk7lmk/v2](https://doi.org/10.17504/protocols.io.e6nvw1jk7lmk/v2)
6. TotalSeq<sup>TM</sup>-A Antibodies and Cell Hashing with 10x Single Cell  
3' Reagent Kit v3 3.1 Protocol [dx.doi.org/10.17504/protocols.io.8aahsae](https://doi.org/10.17504/protocols.io.8aahsae)
7. Chromium Next GEM Single Cell 3' Reagent Kits v3.1(Dual Index) User Guide. Document Number: CG000315. October 2022.



## Attachments



CG000315\_ChromiumNe

X...

4.3MB

## Guidelines

Full protocol is from 10X Genomics. All modifications are to the original protocol (Chromium Next GEM Single Cell 3' Reagent Kits v3.1 User Guide (Rev E) CG000315, support.10xgenomics.com) and TotalSeq<sup>TM</sup>-A Antibodies and Cell Hashing with 10x Single Cell 3' Reagent Kit v3.1 (Dual Index) Protocol (<https://www.biolegend.com/en-us/protocols/totalseq-a-dual-index-protocol>).

## Troubleshooting



## GEM Generation and Barcoding

### 1 Prepare Reagents for use

1. Equilibrate gel beads to room temperature for 00:30:00 .
2. Thaw, vortex, and centrifuge RT Reagent, Template Switch Oligo, and Reducing Agent B. Verify no precipitate in Reducing Agent B.
3. Centrifuge RT Enzyme C before adding to Master Mix.

### 2 Prepare Master Mix

1. Prepare on ice. Pipette mix 15x and centrifuge briefly.
2. Add 31.9  $\mu\text{L}$  Master Mix to each tube of a PCR 8-tube strip on ice.

### 3 Load Chromium Single Cell G Chip

1. Assemble Next GEM Chromium Chip G in a 10X Chip Holder.
2. Must load chip in order according to row label for microfluidic channels to work properly:  
Row 1  $\rightarrow$  Row 2  $\rightarrow$  Row 3
3. Make sure NO bubbles are introduced while loading chip.
4. Dispense 50% glycerol solution into unused Chip Wells (if <8 samples used per chip)
  - a. 75  $\mu\text{L}$  into unused wells in Row 1
  - b. 40  $\mu\text{L}$  into unused wells in Row 2
  - c. 45  $\mu\text{L}$  into unused wells in Row 3

DO NOT add 50% glycerol solution to the bottom row of NO FILL wells. DO NOT use any substitute for 50% glycerol solution.




5. Use Cell Suspension Volume Calculator Table (see step 1.2 of Chromium Next GEM Single Cell 3' v3.1 protocol) to add the appropriate volume of nuclease-free water to Master Mix already in PCR 8-tube strip. Add corresponding volume of gently pipette-mixed single nuclei suspension to Master Mix. Total of 75  $\mu\text{L}$  in each tube. Gently pipette mix, avoid introducing bubbles. DO NOT add nuclease free-water directly to single cell suspension, add instead to Master Mix.

### 6. Cell Loading

- a. load 32,000 nuclei  $\rightarrow$  target recovery 20,000


7. Load 70  $\mu\text{L}$  Master Mix + Nuclei Suspension into the bottom center of each well in Row 1 without introducing bubbles.
8. Snap the Gel Bead strip into a 10X Vortex Adapter. Vortex 00:00:30 . Remove the Gel Bead strip and flick sharply downward to ensure maximum recovery. Confirm there are no bubbles at the bottom of the tubes and liquid levels look even.



10. Load  50  $\mu\text{L}$  gel beads into wells of Row 2. Only puncture the foil seal for gel bead tubes being used. Dispense slowly and without introducing bubbles.
11. Load  45  $\mu\text{L}$  partitioning oil into each Row 3 by pipetting two aliquots of  140  $\mu\text{L}$  .
12. Attach 10X gasket. Align the top-notch. Ensure gasket holes are aligned with the wells. Avoid touching the smooth gasket surface. Do not press down on the gasket. Run chip immediately after loading the partitioning oil.

#### 4 **Run the Chromium Controller**

18m

1. Press the eject button on the controller to eject tray.
2. Place the assembled chip with the gasket in the tray. Press the button to retract the tray.
3. Confirm the program on screen. Press the play button.
4. At completion of the run (~  00:18:00 ), the controller will chime. Proceed immediately to the next step.


Firmware Version 4.0 or higher is required in the Chromium Controller or the Chromium Single Cell Controller used for this protocol.

#### 5 **Transfer GEMs**

1. Chill strip tubes on ice.
2. Press the eject button to remove the chip.
3. Discard the gasket. Open the chip holder. Fold the lid back until it clicks to expose the wells at 45 degrees.
4. Check the volume in rows 1-2. Abnormally high volume in any well indicates a clog.
5. Slowly aspirate 100  $\mu\text{L}$  GEMs from the lowest points of the Recovery Wells without creating a seal between the pipette tips and the wells.

##### Note

Inspect the GEMs in the pipette tip. GEMS should appear opaque and uniform across all channels. Excess partitioning oil (clear) in the pipette tips indicates a potential clog.

6. Slowly dispense (~  00:00:20 ) GEMs into the tube strip on ice with the pipette tips against the sidewalls of the wells.  
If multiple chips are run back-to-back, cap/cover the GEM-containing tube strip and place on ice for no more than 1 h.

#### 6 **GEM-RT Incubation**

3d 0h 50m

1. Set up the cDNA reaction.



2. Incubate in a thermocycler to complete reverse transcription ( 53 °C

00:45:00 , 85 °C 00:05:00 , 4 °C HOLD).

3. Store at 4 °C for up to 72:00:00 or at -20 °C for up to a week. Or proceed to the next step.

## 7 **Post GEM-RT Cleanup & cDNA Amplification - Prepare Reagents for use (See 10X User Guide for details)**

1. Thaw, vortex, and centrifuge a tube of Reducing Agent B and cDNA primers.
2. Maintain Amp Mix on ice after vortex and centrifuge
3. Thaw Cleanup Buffer at 65 °C for 00:10:00 with shaking at max rpm then cool to room temperature. Verify no visible crystals. Cool to RT.

## 8 **Post GEM RT-Cleanup --Dynabead**

17m



1. Add 125 µL Recovery Agent to each sample at room temperature. DO NOT MIX. Wait 00:02:00 .
2. Carefully aspirate 125 µL of pink oil phase from the bottom of the tube. DO NOT aspirate any aqueous sample.
3. Prepare Dynabeads Cleanup Mix (See 10X User Guide for detail)
4. Vortex and briefly centrifuge mixture. Add 200 µL to each sample. Pipette 10x to mix.
5. Incubate at room temperature for 00:10:00 ; pipette mix again 00:05:00 into incubation.
6. Prepare Elution Solution I (See User Guide for Detail). Vortex and centrifuge briefly.
7. Place sample on magnetic separator-HIGH until solution clears. Remove the supernatant (aqueous phase and Recovery Agent). Add 300 µL freshly prepared 80% ethanol. Wait 30 sec. Remove Ethanol. Repeat. Remove the ethanol.
8. Centrifuge briefly. Place on the magnet•Low. Remove remaining ethanol. Air dry for 1 min.
9. Remove from the magnet. Immediately add 35.5 µl Elution Solution I (prepared in step above). Pipette mix (pipette set to 30 µl) without introducing bubbles.
10. Incubate 2 min at room temperature.
11. Place on the magnet•Low until the solution clears. Transfer 35 µl sample to a new tube strip.

## 9 **cDNA Amplification with the addition of ADT and HTO primers**

1. Prepare cDNA Amplification Mix on ice (see 10x User Guide for details and Step 16 of TotalSeq™-A Antibodies and Cell Hashing with 10x Single Cell 3' Reagent Kit v3.1



(Dual Index) Protocol [dx.doi.org/10.17504/protocols.io.8aahsae](https://dx.doi.org/10.17504/protocols.io.8aahsae)). Vortex and centrifuge briefly.








2. Add  65  $\mu\text{L}$  cDNA Amplification Reaction Mix including ADT and HTO primers to  35  $\mu\text{L}$  sample.
3. Pipette mix 15x to mix and centrifuge briefly.
4. Incubate in a thermocycler to amplify cDNA.

#### Note

- Use cycle number optimization table in 10x User Guide) for total number of cycles
- When using nuclei, increase the optimized number of cycles by one cycle
- Maximum: load 1600 nuclei - perform 12 cycles

## 10 cDNA Cleanup --SPRIselect

3d 0h 2m

1. Perform a 0.6X beads purification with SPRIselect reagent (  60  $\mu\text{L}$  )
2. Incubate x 5 minutes at RT.
3. Place on magnet-HIGH until solution clears. **Transfer 70  $\mu\text{L}$  of supernatant to a new tube. Note: the ADT/HTO material will be in the supernatant and the cellular RNA expression pool will be on the SPRI beads. Library construction will be completed in parallel with 3' Gene Expression Library Construction**
4. Add  200  $\mu\text{L}$  freshly prepared 80% ethanol, wait 30 sec and remove ethanol. Repeat.
5. Place on magnet-HIGH until solution clears. Transfer EB. Pipette mix x 15x. Incubate  00:02:00 RT.
6. Place on magnet-HIGH until solution clears. Transfer  40  $\mu\text{L}$  to a new strip tube.
7. Store at  4  $^{\circ}\text{C}$  for up to  72:00:00 or at  -20  $^{\circ}\text{C}$  for up to 4 weeks. Or proceed to the next step.

## 11 cDNA QC & Quantification

1. Run TapeStation to obtain concentration and size
2. Expected size range for amplified cDNA: 400 to 2500 bp
3. cDNA total yield range: 80 to 1000 ng





## Note


**QC Cutoff: Minimum cDNA total yield of 80 ng**

### 3' Gene Expression Library Construction

**12 Prepare Reagents for use**



1. Thaw, vortex, and centrifuge Fragmentation Buffer, Adaptor Oligos, Ligation Buffer, SI Primer. Verify no precipitate.
2. Maintain on ice Fragmentation Enzyme, DNA Ligase, and Amp Mix.

**13 Fragmentation, End Repair & A-tailing**

1. Prepare Fragmentation Mix on ice (See User Manual for volume details). Pipette mix and centrifuge briefly.
2. Transfer  10  $\mu\text{L}$  purified cDNA to a new strip tube





## Note

Note that 10  $\mu\text{L}$  (25%) cDNA sample is sufficient for generating 3' Gene Expression Library. The remaining 30  $\mu\text{L}$  (75%) cDNA sample can be stored at 4 C for up to 72 hours or at -20 C for up to 4 weeks for generating additional 3' Gene Expression Libraries.

3. Add  25  $\mu\text{L}$  Buffer EB to each sample
4. Add  15  $\mu\text{L}$  Fragmentation Mix to each sample
5. Pipette mix, centrifuge briefly, and transfer to pre-cooled block. "SKIP" hold step to initiate the protocol

**14 Post Fragmentation, End Repair & A-tailing Double Sided Size Selection --SPRIselect**

12m

1. Vortex to resuspend the SPRIselect reagent. Add  30  $\mu\text{L}$  SPRIselect reagent (0.6X) to each sample and pipette mix 15x.
2. Incubate  00:05:00 at room temperature.
3. Centrifuge briefly.
4. Place sample on magnetic separator-HIGH until solution clears. DO NOT discard supernatant.
5. Transfer  75  $\mu\text{L}$  supernatant to a new tube strip tube.
6. Vortex to resuspend the SPRIselect reagent. Add  10  $\mu\text{L}$  SPRIselect reagent (0.8X) to each transferred supernatant and pipette mix 15x.



7. Incubate 00:05:00 RT
8. Place sample on magnetic separator-HIGH until solution clears.
9. Remove 80  $\mu$ L supernatant. DO NOT discard any beads.
10. Wash twice with 125  $\mu$ L freshly prepared 80% ethanol.
11. Centrifuge briefly. Place on magnet-LOW until solution clears. Remove remaining ethanol, do not overdry.
12. Remove from magnet. Add 50.5  $\mu$ L Qiagen Buffer EB, pipette mix x 15. Incubate x 00:02:00 RT.
13. Place on magnet-HIGH until clear. Transfer 50  $\mu$ L to a new strip tube.

## 15 Adaptor Ligation

15m

1. Prepare Adaptor Ligation Mix on ice. Pipette mix and centrifuge briefly.
2. Add 50  $\mu$ L Adaptor Ligation Mix to sample. Pipette mix 15x. Centrifuge briefly.
3. Incubate in a thermocycler to ligate adaptor. ( 20  $^{\circ}$ C x 00:15:00 , 4  $^{\circ}$ C Hold)

## 16 Post Ligation Cleanup --SPRIselect



1. Perform a 0.8X beads purification with SPRIselect reagent (add 80  $\mu$ L reagent to each sample, pipette-mix x 15x, incubate 5 min at RT)
2. Place on magnet-HIGH until solution clears, remove supernatant. Wash twice with freshly prepared 80% ethanol (add 200  $\mu$ L , wait 30 sec, remove ethanol and repeat)
3. Centrifuge briefly. Place on magnet-LOW. Remove remaining ethanol and air dry up to 2 min.
4. Remove from magnet. To elute, add 30.5  $\mu$ L Qiagen Buffer EB. Pipette mix x 15. Incubate 2 min at RT. Place on magnet-LOW until solution clears.
5. Transfer 30  $\mu$ L to a new strip tube.

## 17 Sample Index PCR













1. Choose the appropriate sample index sets to ensure that no sample indices overlap in a multiplexed sequencing run.
2. Prepare Sample Index PCR Mix on ice. Pipette mix and centrifuge briefly.
3. Add 50  $\mu$ L Amp Mix to the 30  $\mu$ L sample already in the new strip tube.
4. Add 20  $\mu$ L of an individual Dual Index TT Set A to each sample. Record well assignment. Pipette-mix x5 and centrifuge briefly.
5. Incubate in a thermocycler to incorporate sample indices.

**Note**

Use cycle number recommendation table (see Chromium Next GEM Single Cell 3' v3.1(Dual Index) User Guide, Rev E) for total number of cycles; recommended number of cycles is based on cDNA input

6. Store at  4 °C for up to  72:00:00 . Or proceed to the next step.

**18 Post Sample Index PCR Double Sided Size Selection --SPRIselect**

1. Vortex to resuspend the SPRIselect reagent. Add  60 µL SPRIselect reagent (0.6X) to each sample and pipette mix 15x.
2. Incubate  00:05:00 at room temperature.
3. Centrifuge briefly
4. Place sample on magnetic separator. DO NOT discard supernatant.
5. Transfer  150 µL supernatant to a new tube strip tube.
6. Vortex to resuspend the SPRIselect reagent. Add  20 µL SPRIselect reagent (0.8X) to each sample and pipette mix 15x.
7. Incubate  00:05:00 at room temperature.
8. Centrifuge briefly
9. Place sample on magnetic separator-HIGH until solution clears.
10. Remove  165 µL supernatant. DO NOT discard any beads.
11. With tube still on magnet, add  200 µL freshly prepared 80% ethanol to pellet.  
Wait 30 sec and then remove ethanol.
12. Repeat Step 11 for total of 2 washes then centrifuge briefly. Place on magnet-LOW.  
Remove remaining ethanol.
13. Remove from magnet. Add  35.5 µL Qiagen Buffer EB. Pipette-mix x 15)
14. Incubate 2 min at RT.
15. Place on magnet-LOW until solution clears then transfer  35 µL to a new strip tube.
16. Store at  4 °C for up to  72:00:00 or at  -20 °C for long-term storage.

**Note****QC Post Library Construction:**

- **Quantify Library (e.g. using Qubit dsDNA HS Assay)**
- **Estimate Library size range (e.g. using TapeStation or BioAnalyzer) - expected size range of 300 to 800 bp, average: 475 bp**
- If additional peaks below 200 bp are observed, see User Guide v 3.1 for recommendations

## Amplify Dual-Indexed ADT Sequencing Library

- 19 Aliquot 22.5 uL of each sample (ADT cDNA) to a new, labeled 0.5-mL strip tube.
- 20 Add 25 uL of 2X Q5 HotStart HiFi Master Mix to each sample.
- 21 Add 2.5 uL of appropriate indexed primer-pair (10uM each primer) to each sample. Mix by pipetting and spin briefly.
- 22
- | Step | Temperature        | Time |
|------|--------------------|------|
| 1    | 95°C               | Hold |
| 2    | 95°C               | 3:00 |
| 3    | 95°C               | 0:20 |
| 4    | 60°C               | 0:30 |
| 5    | 72°C               | 0:20 |
| 6    | Complete 12 cycles |      |
| 7    | 72°C               | 5:00 |
| 8    | 4°C                | Hold |
- 23 Add 60 ul AMPure beads (1.2X ratio) to each sample and pipette to mix 10X
- 24 Incubate 5 minutes at room temperature.
- 25 Place tube on magnet and wait 1 minute until solution is clear.



- 26 Carefully remove and discard the supernatant.
- 27 Add 200  $\mu$ l 80% Ethanol to the tube without disturbing the pellet and stand for 30 seconds (first Ethanol wash).
- 28 Carefully remove and discard the ethanol wash.
- 29 Add 200  $\mu$ l 80% Ethanol to the tube without disturbing the pellet and stand for 30 seconds (second Ethanol wash).
- 30 Carefully remove and discard the ethanol wash.
- 31 Centrifuge tube briefly and return it to magnet.
- 32 Remove and discard any remaining ethanol and allow the beads to air dry for 2 minutes.
- 33 Resuspend beads in 32  $\mu$ l EB.
- 34 Pipette mix vigorously and incubate at room temperature for 5 minutes.
- 35 Place tube on magnet and transfer 30  $\mu$ l supernatant to low-bind tube.
- 36 Quantify library by standard methods (Qubit & BioAnalyzer). ADT library will be around 210 bp.

## Sequencing

- 37 NovaSeq Sequencing (target - 100,000 reads per nucleus)
  1. Paired End, Dual Indexing
    - a. Read 1: 28 cycles
    - b. i7 Index: 10 cycles
    - c. i5 Index: 10 cycles



- d. Read 2: 91 cycles
- 2. Library Loading
  - a. 10X recommended Loading concentration: 300 pM
  - b. Optional: 1% PhiX
- 3. Output
  - a. SP: 650–800 million reads
  - b. S1: 1.3–1.6 billion reads
  - c. S2: 3.3 –4.1 billion reads
  - d. S4: 8-10 billion reads

## Cell Ranger Mapping and Analysis Pipeline

- 38 Create a feature reference file corresponding to the antibody panel and library mapping file corresponding to the file paths for the ADT and GEX libraries. Run cellranger count, incorporating both of these files and the combined Human+COVID reference
- 39 Create fastq from bcl files using cellranger mkfastq with default parameters
- 40 Insert 10X provided well indices into Clarity LIMS generated sample sheet.
- 41 Create combined human and COVID reference by downloading the Human 2020-A reference from <https://www.10xgenomics.com/support/software/cell-ranger/latest/release-notes/cr-reference-release-notes#2020-a>. Download the full COVID sequence from <https://www.ncbi.nlm.nih.gov/nuccore/MT547814>. Append the COVID fasta onto the human genome fasta. Create a record within the human GTF file corresponding to the full COVID sequence.



## Protocol references

### References

1. <https://www.lungmap.net/about-lungmap/lungmap-objectives/>
2. 602.2 Donor Acceptance Criteria for URMHC HTC SenNet Inclusion [dx.doi.org/10.17504/protocols.io.biz7kf9n](https://doi.org/10.17504/protocols.io.biz7kf9n).
3. 603.3 & 604.5\_URMC\_HTC\_Lung and Lobe Processing for SenNet  
[dx.doi.org/10.17504/protocols.io.n2bvj395plk5/v1](https://doi.org/10.17504/protocols.io.n2bvj395plk5/v1)
4. 702.B.3 URMHC HTC Lung Tissue Digestion SOP + Worksheet 052820  
[dx.doi.org/10.17504/protocols.io.biz5kf86](https://doi.org/10.17504/protocols.io.biz5kf86)
5. 709.1 Staining of Dissociated Lung Cells for Cellular Indexing of Transcriptomes and Epitopes by Sequencing (CITE sequencing) V.2 [dx.doi.org/10.17504/protocols.io.e6nvw1jk7lmk/v2](https://doi.org/10.17504/protocols.io.e6nvw1jk7lmk/v2)
6. TotalSeq<sup>TM</sup>-A Antibodies and Cell Hashing with 10x Single Cell  
3' Reagent Kit v3 3.1 Protocol [dx.doi.org/10.17504/protocols.io.8aahsae](https://doi.org/10.17504/protocols.io.8aahsae)
7. Chromium Next GEM Single Cell 3' Reagent Kits v3.1(Dual Index) User Guide. Document Number: CG000315. October 2022.

<https://www.10xgenomics.com/support/software/cell-ranger/latest/release-notes/cr-reference-release-notes#2020-a>.

<https://www.ncbi.nlm.nih.gov/nuccore/MT547814>