



Aug 29, 2022

Lysis and transduction of E. coli with P1 phage - creation of double-knockout mutants

DOI

dx.doi.org/10.17504/protocols.io.n92ldpbmxl5b/v1

Saul Moore¹

¹Imperial College London, MRC London Institute of Medical Sciences

Saul Moore: This protocol was carried out by Cassandra Backes of the Host-Microbe Co-Metabolism laboratory, MRC-LMS

Behavioural Genomics



Saul Moore

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.n92ldpbmxl5b/v1>

Protocol Citation: Saul Moore 2022. Lysis and transduction of E. coli with P1 phage - creation of double-knockout mutants. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.n92ldpbmxl5b/v1>



License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: August 28, 2022

Last Modified: August 29, 2022

Protocol Integer ID: 69283

Keywords: knockout mutants from single gene deletion mutant, single gene deletion mutant, knockout mutants the creation, knockout mutant, p1 phage, cassandra backes of the host, metabolism, transduction, cassandra backe

Disclaimer

DISCLAIMER – FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to [protocols.io](#) is not peer reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with [protocols.io](#), can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.

Abstract

The creation of double-knockout mutants from single gene deletion mutants of the Keio Collection was performed by Cassandra Backes of the Host-Microbe Co-Metabolism laboratory, MRC-LMS.

Troubleshooting

Lysis

- 1 Dilute an overnight culture of donor strain grown with selection for the marker to be transduced 1:100 in fresh LB supplemented with 5 mM CaCl₂, and 0.2% (or 10 mM) D-Glucose. DO NOT ADD ANTIBIOTIC TO THIS CULTURE. Also prepare an extra tube (for no-phage control).
E. g. For each 5 mL culture add:
 - 25 ul of 1 M CaCl₂,
 - 50 uL of 20%(or 1M) D-Glucose,
 - 50 uL O/N culture.
- 2 Incubate for 30 to 45 mins at 37C with shaking (220rpm) (or until you get cells to early log phase – tube should be slightly turbid, but noticeable growth)
- 3 Add P1 phage – 1 to 2 drops – using sterile/discardable Pasteur pipette to cultures (except to no-phage control tube)
- 4 Continue to incubate, with shaking. Monitor for 1–3 hr until the culture has lysed (you'll see cellular debris in the tube and the culture will have significantly lessened in its turbidity – compare with no-phage control tube).
 - a. Might take longer than 3 hours in some cases.
 - b. Sometimes the culture does not lessen in turbidity but if you can see some cellular debris when you shake the tube then you can proceed to the next step.
- 5 Transfer lysed culture into 15 mL blue lid falcon tubes* and in the chemical fume hood add 125 uL of Chloroform and vortex for 1 min.
 - a. The chloroform ensures complete cell lysis and kills bacteria. It does not harm phage, but care should be taken to avoid transferring the chloroform to the storage tubes during step 8, because storing of phage in presence of chloroform can result in decreased viability of the stock overtime.
 - b. *Chloroform can dissolve the plastic on the 30 mL sterillin tubes – so make sure to transfer lysed culture to 15 mL falcon tubes first!
- 6 Spin down 10 min at 2000 RPM.
- 7 Pour supernatant into 5 mL syringe and sterile filter directly into 2× 2 mL tubes.
 - a. Avoid pouring the chloroform into syringe as it will be very hard to filter through otherwise.



- 8 Store at 4C (write date and donor strain genotype)

Transduction

1d

- 9 Overnight culture of the donor strain, carrying the mutation to be transduced, in LB + 5 mM CaCl₂ at 37°C. 1d
- 10 In a 15mL Falcon, add 50 uL of overnight culture + 5-10 uL phage P1, 20 min at 37°C WITHOUT shaking.
Negative control: tube without P1 phage
- 11 Add 5 mL LB CaCl₂ 5 mM and transfer to an Erlenmeyer flask at 37°C with stirring until OD₆₀₀ ≈ 1
cell lysis and the stock is ready
- 12 Put everything back in a Falcon 15mL and add 500 uL chloroform
Store at 4°C. (CaCl₂ concentration can be increased to 10 mM)
- 13 Overnight culture of recipient strain in LB
- 14 Dilute the overnight culture to 1/100th in 5 mL CaCl₂
- At OD₆₀₀ ≈ 1, take 5 Eppendorf tubes
1. Negative control tube – P1 lysate + 750 uL of recipient cells
 2. 5 uL tube of P1 lysate + 750 uL of recipient cells
 3. 10 uL tube of P1 lysate + 750 uL of recipient cells
 4. 15 uL tube of P1 lysate + 750 uL of recipient cells
 5. 20 uL tube of P1 lysate + 750 uL of recipient cells
- Incubate for 20 min at 37°C WITHOUT stirring
- 15 Invert, centrifuge for 5 min at 5000 rpm and empty the supernatant
- 16 Resuspend the pellet in 1 mL of LB + 7.5 mM citrate and leave for 30 to 45 min at RT
- 17 Centrifuge again, resuspend in 100 µL and plate on LB + AB + 7.5 mM Citrate
- 18 The next day, restrict the clones on LB + citrate dishes to have isolated clones



19 Verification of clones by PCR on colonies