



Jan 05, 2023

Version 1

# Lysogeny Broth (LB) medium V.1

DOI

[dx.doi.org/10.17504/protocols.io.ewov1o262lr2/v1](https://dx.doi.org/10.17504/protocols.io.ewov1o262lr2/v1)



Andreas Sagen<sup>1</sup>

<sup>1</sup>University of Oslo



Andreas Sagen

University of Oslo

## Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.ewov1o262lr2/v1>

**Protocol Citation:** Andreas Sagen 2023. Lysogeny Broth (LB) medium. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.ewov1o262lr2/v1>

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



**Protocol status:** Working

**I use this protocol and it's working**

**Created:** January 05, 2023

**Last Modified:** January 06, 2023

**Protocol Integer ID:** 74793

**Keywords:** LB, bacteria, e coli, medium lysogeny broth, lysogeny broth, lb broth, broth, lb medium, lb, cultivating escherichia, rich medium, bacteria

## Abstract

Lysogeny broth (LB) is a nutritionally rich medium which is primarily used for the growth of bacteria<sup>[1]</sup>. LB broth is commonly used when cultivating *Escherichia coli*. There exist different formulations of LB and lead to the development of derivations for specialized use.

## Guidelines

Follow step by step, unless stated otherwise. Equipment needed should be standard to a microbiology lab.

## Materials

Analytical scale, autoclave, bottle, weight vessel, LAF bench

## Protocol materials

☒ Sodium chloride Merck MilliporeSigma (Sigma-Aldrich) Catalog #S9625

☒ Tryptone Merck Millipore (EMD Millipore) Catalog #T9410

☒ Yeast Extract Merck MilliporeSigma (Sigma-Aldrich) Catalog #Y0875

☒ Tryptone Merck Millipore (EMD Millipore) Catalog #T9410

☒ Yeast Extract Merck MilliporeSigma (Sigma-Aldrich) Catalog #Y0875

☒ Sodium chloride Merck MilliporeSigma (Sigma-Aldrich) Catalog #S9625

☒ Agar Merck MilliporeSigma (Sigma-Aldrich) Catalog #A1296

## Troubleshooting

## Safety warnings

⚠ When removing autoclaved components, be sure to take care as this can be very hot. If using antibiotics, use sufficient PPE to protect yourself, as some can be toxic to humans.

## Before start

Prepare glassware by cleaning it, and ensure that scale is sufficiently calibrated




## 500 mL LB-Lennox (broth) medium

1 All compounds are measured using a high precision analytical scale from powdered compounds. Each compound is measured to within 1% of the target weight. All compounds are mixed in a Duran bottle

1.1 Fill the bottle with  400 mL double-distilled water


1.2 Measure  5000 mg Tryptone,  2500 mg Yeast extract and  2500 mg Sodium chloride

Powdered compounds:


 Tryptone **Merck Millipore (EMD Millipore) Catalog #T9410**


 Yeast Extract **Merck MilliporeSigma (Sigma-Aldrich) Catalog #Y0875**

 Sodium chloride **Merck MilliporeSigma (Sigma-Aldrich) Catalog #S9625**

1.3 Add powdered solids into bottle, and use a magnetic mixer with a stir bar to mix for  00:05:00

5m

1.4 Adjust pH while mixing to  6.7 using concentrated sodium hydroxide

1.5 Add distilled water to a total of  500 mL

1.6 Autoclave liquid at  121 °C for  00:15:00

15m

### Note

Cool to 50°C and supplement with antibiotics as appropriate





## 500 mL LB-Lennox (agar) medium

20m


2 All compounds are measured using a high precision analytical scale from powdered compounds. Each compound is measured to within 1% of the target weight. All compounds are mixed in a Duran bottle

2.1 Fill the bottle with  400 mL double-distilled water



2.2 Measure  5000 mg Tryptone,  2500 mg Yeast extract,  2500 mg Sodium chloride and  7500 mg agar

Powdered compounds:


 Tryptone **Merck Millipore (EMD Millipore) Catalog #T9410**

 Yeast Extract **Merck MilliporeSigma (Sigma-Aldrich) Catalog #Y0875**


 Sodium chloride **Merck MilliporeSigma (Sigma-Aldrich) Catalog #S9625**

 Agar **Merck MilliporeSigma (Sigma-Aldrich) Catalog #A1296**

2.3 Add powdered solids into bottle, and use a magnetic mixer with a stir bar to mix for

 00:05:00

5m

2.4 Adjust pH while mixing to  6.7 using concentrated sodium hydroxide

2.5 Add distilled water to a total of 500 mL

2.6 Autoclave liquid at  121 °C for  00:15:00

15m

#### Note

Cool to 50°C and supplement with antibiotics as appropriate

Agar can be stored, then reheated to 50°C to be poured