



Jan 06, 2023

Version 2

Lysogeny Broth (LB) medium V.2

DOI

<https://dx.doi.org/10.17504/protocols.io.ewov1o262lr2/v2>



Andreas Sagen¹

¹University of Oslo



Andreas Sagen

University of Oslo

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.ewov1o262lr2/v2>

Protocol Citation: Andreas Sagen 2023. Lysogeny Broth (LB) medium. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.ewov1o262lr2/v2> Version created by **Andreas Sagen**

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this protocol and it's working

Created: January 06, 2023

Last Modified: January 06, 2023

Protocol Integer ID: 74863

Keywords: LB, bacteria, Lennox broth, E. coli, Escherichia coli, medium lysogeny broth, lysogeny broth, several common formulations of lb, formula of the lb medium, abbreviation lb, lennox broth, helpful for bacterial growth, broth, luria broth, lb medium, bacterial growth, essential amino acid, lb, ingredient, essential amino acids such as peptide, similar composition of ingredient, magnesium, compound, vitamin, growing bacteria, yeast extract, peptide, sodium chloride, peptone, plethora of organic compound, mineral, osmotic balance, bertani medium, sodium, nitrogen, sodium ions for transport, several common formulation, bacteria, organic compound, yeast, casein peptone

Abstract

Lysogeny broth (LB) is a nutritionally rich medium which is primarily used for the growth of bacteria^[1]. The initialism is also commonly, albeit incorrectly, taken to mean **Luria broth**, **Lennox broth**, or **Luria-Bertani medium**. According to its creator Giuseppe Bertani, the abbreviation **LB** was actually intended to stand for **lysogeny broth**. The formula of the **LB medium** was published in 1951 in the first paper of Bertani on lysogeny^[2]. There are several common formulations of **LB**. Although they are different, they generally share a somewhat similar composition of ingredients used to promote growth, including the following: Peptides and casein peptones, Vitamins (including B vitamins), Trace elements (e.g. nitrogen, sulfur, magnesium) and Minerals. Sodium ions for transport and osmotic balance are provided by sodium chloride (NaCl). Tryptone is used to provide *essential amino acids* such as peptides and peptones to the growing bacteria, while the yeast extract is used to provide a plethora of organic compounds helpful for bacterial growth. These compounds include vitamins and certain trace elements.

Guidelines

Follow step by step, unless stated otherwise. Equipment needed should be standard to a microbiology lab.

Materials

Analytical scale, autoclave, bottle, weight vessel, LAF bench

Protocol materials

⊗ Tryptone Merck Millipore (EMD Millipore) Catalog #T9410

⊗ Yeast Extract Merck MilliporeSigma (Sigma-Aldrich) Catalog #Y0875

⊗ Sodium chloride Merck MilliporeSigma (Sigma-Aldrich) Catalog #S9625

⊗ Agar Merck MilliporeSigma (Sigma-Aldrich) Catalog #A1296



Safety warnings

- ⚠ When removing autoclaved components, be sure to take care as this can be very hot. If using antibiotics, use sufficient PPE to protect yourself, as some can be toxic to humans.

Before start

Prepare glassware by cleaning it, and ensure that scale is sufficiently calibrated



500 mL LB-Lennox (broth) medium

- 1 All compounds are measured using a high precision analytical scale from powdered compounds. Each compound is measured to within 1% of the target weight. All compounds are mixed in a Duran bottle
- 1.1 Fill the bottle with  400 mL double-distilled water
- 1.2 Measure  5000 mg Tryptone,  2500 mg Yeast extract and  2500 mg Sodium chloride

Powdered compounds:
 -  Tryptone **Merck Millipore (EMD Millipore) Catalog #T9410**
 -  Yeast Extract **Merck MilliporeSigma (Sigma-Aldrich) Catalog #Y0875**
 -  Sodium chloride **Merck MilliporeSigma (Sigma-Aldrich) Catalog #S9625**
- 1.3 Add powdered solids into bottle, and use a magnetic mixer with a stir bar to mix for  00:05:00
- 1.4 Adjust pH while mixing to  6.7 using concentrated sodium hydroxide
- 1.5 Add distilled water to a total of  500 mL

5m

- 1.6 Autoclave liquid at  121 °C for  00:15:00

15m

Note

Cool to 50°C and supplement with antibiotics as appropriate

500 mL LB-Lennox (agar) medium

20m

- 2 All compounds are measured using a high precision analytical scale from powdered compounds. Each compound is measured to within 1% of the target weight. All compounds are mixed in a Duran bottle
- 2.1 Fill the bottle with  400 mL double-distilled water



2.2 Measure  5000 mg Tryptone,  2500 mg Yeast extract,  2500 mg Sodium chloride and  7500 mg agar

Powdered compounds:

 Tryptone **Merck Millipore (EMD Millipore) Catalog #T9410**

 Yeast Extract **Merck MilliporeSigma (Sigma-Aldrich) Catalog #Y0875**

 Sodium chloride **Merck MilliporeSigma (Sigma-Aldrich) Catalog #S9625**

 Agar **Merck MilliporeSigma (Sigma-Aldrich) Catalog #A1296**

2.3 Add powdered solids into bottle, and use a magnetic mixer with a stir bar to mix for

 00:05:00

5m

2.4 Adjust pH while mixing to  6.7 using concentrated sodium hydroxide

2.5 Add distilled water to a total of 500 mL

2.6 Autoclave liquid at  121 °C for  00:15:00

15m

Note

Cool to 50°C and supplement with antibiotics as appropriate

Agar can be stored, then reheated to 50°C to be poured