

Jan 11, 2026

Magic Red® Cathepsin B Assay Protocol

DOI

<https://dx.doi.org/10.17504/protocols.io.kqdg32mb1v25/v1>

Ali Ghoochani^{1,2,3,4}, Monther Abu-Remaileh^{1,2,3,4}

¹Department of Chemical Engineering, Stanford University, Stanford, CA 94305, USA;

²Department of Genetics, Stanford University, Stanford, CA 94305, USA;

³The Institute for Chemistry, Engineering and Medicine for Human Health (Sarafan ChEM-H), Stanford University, Stanford, CA 94305, USA;

⁴Aligning Science Across Parkinson's (ASAP) Collaborative Research Network, Chevy Chase, MD, 20815, USA

asap



Ali Ghoochani

Stanford University

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.kqdg32mb1v25/v1>

Protocol Citation: Ali Ghoochani, Monther Abu-Remaileh 2026. Magic Red® Cathepsin B Assay Protocol. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.kqdg32mb1v25/v1>



License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: October 13, 2024

Last Modified: January 11, 2026

Protocol Integer ID: 109768

Keywords: red fluorescent signal, magic red substrate, based assay, fluorophore cresyl violet, fluorescence microscope, fluorescence, cresyl violet fluorophore, immunochemistry technology, assay

Abstract

The Magic Red Cathepsin B Assay Kit from ImmunoChemistry Technologies is a fluorescence-based assay designed to detect and measure intracellular cathepsin B activity in live cells. This kit utilizes the Magic Red substrate (MR-RR2), which contains the fluorophore cresyl violet conjugated to a cathepsin B-specific dipeptide sequence (arginine-arginine). Upon cleavage by active cathepsin B, the substrate releases the cresyl violet fluorophore, resulting in a red fluorescent signal that can be detected using a fluorescence microscope or plate reader.

Materials

Magic Red Cathepsin B Assay Kit (ImmunoChemistry Technologies, Cat# 937.



Prepare Magic Red Staining Solution

- 1 Reconstitute Magic Red Substrate (MR-RR2) with 50 μ L DMSO.

Add Staining Solution to Cells

- 2 Once cells are ready for the experiment:
 - 2.1 Dilute the stock 1:10 with diH₂O to form the staining solution.
 - 2.2 Dilute the Magic Red staining solution 1:25. Add 20 μ L staining solution to 480 μ L media.
 - 2.3 Incubate the cells with the Magic Red staining solution for  00:30:00 at  37 °C 

Imaging

- 3 Immediately image the cells using a fluorescent or confocal microscope.
 - 3.1 Use filters appropriate for Magic Red detection (excitation: 550-590 nm, emission: >610 nm)

Note: If the background fluorescence appears high or if there is excess staining solution, a gentle PBS wash can be considered to reduce background.