

Jun 25, 2024

Maintenance & Differentiation: Embryonic Mouse Hippocampal Cells (CLU198)

DOI

<https://dx.doi.org/10.17504/protocols.io.bp2l628wrgqe/v1>

Md Razaul Karim¹

¹University of Minnesota

Md Razaul Karim: Lee Lab

Team Lee



Jane Balster

ASAP - Team Lee

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.bp2l628wrgqe/v1>

Protocol Citation: Md Razaul Karim 2024. Maintenance & Differentiation: Embryonic Mouse Hippocampal Cells (CLU198). protocols.io <https://dx.doi.org/10.17504/protocols.io.bp2l628wrgqe/v1>



License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: June 07, 2024

Last Modified: June 25, 2024

Protocol Integer ID: 102396

Keywords: Maintenance, Differentiation, Embryonic Mouse Hippocampal Cells (CLU198), differentiation of embryonic mouse hippocampal cell, embryonic mouse hippocampal cell, cell, clu198, mouse

Funders Acknowledgements:

ASAP

Grant ID: 000592

Abstract

This protocol details the maintenance & differentiation of Embryonic Mouse Hippocampal Cells (CLU198).

Materials

Media:

1. **Maintenance Media:** DMEM full media containing (DMEM/10% FBS/1% Pen-Strep).

2. **Differentiation Media:** Complete Neurobasal media containing

- I. Neurobasal-A (1x) Media Gibco
- II. Glutamax (100x) Gibco-(1x final)
- III. Pen-Strep (1% final)
- IV. B27-Supple (50x) Gibco-(**0.5x** final)

For differentiation:

- 12 Well plate (1 ml/well)
- 6 Well plate/35 mm dish (9.5 cm²-2 ml/well)
- 60 mm dish (21 cm²-3.5 ml/dish)
- 100 mm dish (56 cm²-10 ml/dish)
- 250 ml Flask for sub-culture/maintenance (10 ml) [Tissue culture flask-Greiner bio-one- Cat.No.-658 170]



CELL CULTURE FLASK, 250 ML, 75 CM², PS, RED STANDARD SCREW CAP, CLEAR, CELLSTAR® TC, STERILE, 5 PCS. **greiner bio-one Catalog #658170**



Maintenance

- 1 For regular maintenance of CLU cells, use DMEM full media.

Differentiation (CLU 198): Takes about a week of differentiation.

33m

2

30m

Note

No need to add Retinoic acid.



Day-01 (Mon): Plating with DMEM full media.

- Warm DMEM full media, PBS, and Trypsin in the 37 °C bead bath for 00:30:00 . Clean the working area by using 70% ethanol.

- 2.1 Sup out old media without touching cells.

- 2.2 Wash by adding 5 mL PBS slowly, rinse, and rock back and forth.



- 2.3 Add 2 mL - 3 mL trypsin (0.25%); keep in incubator for 00:03:00 .

3m



- 2.4 Check under microscope if cells are detached, add 5 mL media and transfer to a tube.

- 2.5 Spin 300 x g, 00:03:00 .



- 2.6 Sup out and add 10 mL fresh media & re-suspend cells gently and carefully.

- 2.7 Count cells density and split accordingly. 15,000 cells/ml for maintenance



- (i) Usually 1.5×10^4 /ml cells for Biochem, and
- (ii) 0.5×10^4 /ml cells for IF.

3 **Day-02 (Tue):**

Replace with Complete Neurobasal Media.

4 **Day-03 (Wed):**

Rest.

5 **Day-04 (Thu):**

Rest.

6 **Day-05 (Fri):**

Replace with Complete Neurobasal Media/ (Start drug treat if necessary).

7 **Day-06 (Sat):**

Rest.

8 **Day-07 (Sun):**

Rest.

9 **Day-08 (Mon):**

Replace with Neurobasal Media/Drug treat.

10 **Day-09 (Tue):**

Drug treat if necessary /Harvesting.

11 **Day-10 (Wed):**

Drug treat if necessary /Harvesting.

Cells harvesting:

20m 2s

12 Wash once with cold PBS.



13 Add cold lysis buffer.



14 Keep  On ice & scrap immediately in Eppendorf tube.



15 Sonicate (10 S on  00:00:02 off 20% Amplitude, 2 Pulses).

2s

16 Boil ( 100 °C ,  00:10:00).

10m



17 Centrifuge  13.000 rpm, 4°C, 00:10:00 / Collect sup.

10m



18 Keep in  -80 °C Freezer.



19 BCA to measure protein concentration.

20 Prepare with sample buffer and run WB analysis.