



Jun 26, 2024

Maintenance & Differentiation: SHSY-5Y Neuroblastoma Cells

DOI

<https://dx.doi.org/10.17504/protocols.io.rm7vzj6r2lx1/v1>

Md Razaul Karim¹

¹University of Minnesota

Md Razaul Karim: Lee Lab

Team Lee



Jane Balster

ASAP - Team Lee

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.rm7vzj6r2lx1/v1>

Protocol Citation: Md Razaul Karim 2024. Maintenance & Differentiation: SHSY-5Y Neuroblastoma Cells. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.rm7vzj6r2lx1/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this protocol and it's working

Created: June 07, 2024

Last Modified: June 26, 2024

Protocol Integer ID: 102398

Keywords: Maintenance, Differentiation, SHSY-5Y Neuroblastoma Cells, Cell harvesting, neuroblastoma cells this protocol detail, neuroblastoma cell, cells this protocol detail, cell

Funders Acknowledgements:

asap

Grant ID: 000592

Abstract

This protocol details the maintenance & differentiation of SHSY-5Y Neuroblastoma Cells.



Materials

Media:

1. Maintenance Media: DMEM full media containing (DMEM/10% FBS/1% Pen-Strep).
2. Differentiation Media: Complete Neurobasal media containing
 - I. Neurobasal-A (1x) Media Gibco
 - II. Pen-Strep (1% final)
 - III. B27-Supple (50x) Gibco-(**1.0x** final)
 - IV. Retinoic Acid in DMSO (Stock: 10 mM); use 10 μ M (Final). 1 μ l of stock per 1ml media.

Note

Add **Ingredient IV Retinoic acid** freshly each time

For differentiation:

- 12 Well plate (1 ml/well)
- 6 Well plate/35 mm dish (9.5 cm²-2 ml/well)
- 60 mm dish (21 cm²-3.5 ml/dish)
- 100 mm dish (56 cm²-10 ml/dish)
- 250 ml Flask for sub-culture/maintenance (10 ml) [Tissue culture flask-Greiner bio-one- Cat.No.-658 170]



CELL CULTURE FLASK, 250 ML, 75 CM², PS, RED STANDARD SCREW CAP, CLEAR, CELLSTAR® TC, STERILE, 5 PCS. **greiner bio-one Catalog #658170**



Maintenance

- 1 For regular maintenance of SHSY-5Y cells, use DMEM full media.

Differentiation (SHSY-5Y) Cells: About 5-7 days differentiation. No need Glutamax

3m

2

30m

Note

No need Glutamax.



Day-01 (Mon): Plating with DMEM full media

- Warm DMEM full media, PBS, and Trypsin in the 37 °C bead bath for 00:30:00 . Clean the working area by using 70% ethanol.

- 2.1 Sup out old media without touching cells.

- 2.2 Wash by adding 5 mL PBS slowly, rinse, and rock back and forth.



- 2.3 Add 2 mL - 3 mL trypsin (0.25%); keep in incubator for 00:03:00 .

3m



- 2.4 Check under microscope if cells are detached, add 5 mL media and transfer to a tube.

- 2.5 Spin 300 x g, 00:03:00 .



- 2.6 Sup out and add 10 mL fresh media & re-suspend cells gently and carefully.

- 2.7 Count cells density and split accordingly. 15,000 cells/ml for maintenance



- Usually 1.0×10^4 /ml cells for Biochem, and
- 0.5×10^4 /ml cells for IF.

3 **Day-02 (Tue):**

Replace with Complete Neurobasal Media (Without Glutamax).

Note

#Add Retinoic Acid freshly

4 **Day-03 (Wed):**

Rest.

5 **Day-04 (Thu):**

Rest.

6 **Day-05 (Fri):**

Replace with Complete Neurobasal Media (Without Glutamax) /(Start drug treat if necessary)

Note

Add Retinoic Acid freshly

7 **Day-06 (Sat):**

Rest.

8 **Day-07 (Sun):**

Rest.

9 **Day-08 (Mon):**

Replace with Complete Neurobasal Media (**No Retinoic acid**)/Drug treat.

10 **Day-09 (Tue):**



Drug treat if necessary /Harvesting.

11 Day-10 (Wed):

Drug treat if necessary /Harvesting.

Cells harvesting:

20m 2s

12 Wash once with cold PBS.



13 Add cold lysis buffer.

14 Keep On ice & scrap immediately in Eppendorf tube.

15 Sonicate (10 S on 00:00:02 off 20% Amplitude, 2 Pulses)

2s

16 Boil (100 °C , 00:10:00).

10m



17 Centrifuge 13.000 rpm, 4°C, 00:10:00 / Collect sup.

10m



18 Keep in -80 °C Freezer.



19 BCA to measure protein concentration.

20 Prepare with sample buffer and run WB analysis.