


Feb 06, 2019

Making Differentiation Media for SH-SY5Y

DOI

<https://dx.doi.org/10.17504/protocols.io.xuufnww>



Kenneth Schackart¹

¹University of Arizona

Yoon Lab



Kenneth Schackart

University of Arizona

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.xuufnww>

Protocol Citation: Kenneth Schackart 2019. Making Differentiation Media for SH-SY5Y. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.xuufnww>



License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: In development

We are still developing and optimizing this protocol

Created: February 06, 2019

Last Modified: February 06, 2019

Protocol Integer ID: 20084

Keywords: differentiation media for sh, making differentiation media, differentiation medium for sh, differentiation medium, sy5y

Abstract

This describes how to make differentiation medium for SH-SY5Y using all-*trans*-retinoic acid.

Materials

- all-*trans*-retinoic acid
- B-27 Supplement
- L-glutamine
- Neurobasal Medium (NBM)
- 100% Ethanol
- Autoclaved deionized water (ADIW)



Make ATRA Stock Solution

- 1 In a 15 mL centrifuge tube, add  14.25 mL ethanol and  0.75 mL autoclaved deionized water. This will make a 95% ethanol solution.
- 2 Transfer  3330 μL of 95% EtOH solution to another 15 mL centrifuge tube.
- 3 Weigh out  5 mg ATRA.

Note

ATRA is light sensitive, so do this step with the lights off.
- 4 Add the  5 mg ATRA to the  3330 μL 95% EtOH solution. This makes a 5 mM stock ATRA solution.
- 5 Vortex mix the solution.
- 6 This solution is used at a rate of  100 μL stock solution per  50 mL differentiation media. The final concentration should be 10 μM .

Note

We most recently made 10 aliquots, each with  62 μL . This is to make 10  30 mL media batches, each using  60 μL .



7 Stock solution can be stored for up to 6 weeks at 4°C .

Make I-Glutamine Stock Solution

8 Weigh out 175 mg I-glutamine.

9 Aliquot 15 mL NBM into a 15 mL centrifuge tube.

10 Add the I-glutamine to the NBM.

11 Vortex mix the solution.

12 Aliquot the solution into 1.5 mL aliquots.

This solution is 20x. The concentration is 80 mM and the final desired concentration is 4 mM.

Note

This protocol makes 10 aliquots of 1.5 mL . Each stock aliquot makes 30 mL differentiation medium.

13 Store aliquots at -20°C .

Aliquot B-27 Supplement



- 14 Aliquot  600 μL batches of 50x B-27 supplement.
- 15 Each aliquot will make  30 mL differentiation media.
- 16 Store aliquots at  $-20\text{ }^{\circ}\text{C}$.

Making the media

17

Note

Do not mix the media until you are ready to change the media in a flask. It degrades quickly.

Warm  30 mL NBM in water bath.

- 18 Thaw one 50x B-27 supplement aliquot, and one 20x l-glutamine supplement.
- 19 In biosafety cabinet, combine the following in a 50 mL tube
 -  29.34 mL Nuerobasal Medium
 -  600 μL 50x B-27 Supplement
 -  1.5 mL 20x l-glutamine Stock Solution
 -  60 μL 500x All-*trans*-retinoic Acid Stock Solution
- 20 Sterile filter the differentiation media.
- 21 Use immediately.