



Jun 13, 2020

Making inexpensive light-powered *Chlamydomonas reinhardtii* (a green micro-alga) bead bracelets/necklaces for teaching the interplay of photosynthesis and cellular respiration to K4-K16 students.

 [PeerJ](#)

DOI


dx.doi.org/10.17504/protocols.io.bgpyjvpw

Mautusi Mitra¹, Sara Broom¹, Kysis Pinto¹, Sovi-Mya Doan Wellons², Ariel Dominique Roberts¹

¹University of West Georgia; ²The Heritage School



Mautusi Mitra



Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

[Create free account](#)

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.bgpyjvpw>

External link: <https://www.facebook.com/pages/UWG-Mitra-Lab/1618923318377133>



Protocol Citation: Mautusi Mitra, Sara Broom, Kysis Pinto, Sovi-Mya Doan Wellons, Ariel Dominique Roberts 2020. Making inexpensive light-powered *Chlamydomonas reinhardtii* (a green micro-alga) bead bracelets/necklaces for teaching the interplay of photosynthesis and cellular respiration to K4-K16 students.. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.bgpyjvpw>

Manuscript citation:

Mitra M, Broom SM, Pinto K, Wellons SD, Roberts AD. 2020. Engaging inexpensive hands-on activities using *Chlamydomonas reinhardtii* (a green micro-alga) beads to teach the interplay of photosynthesis and cellular respiration to K4–K16 Biology students. *PeerJ* 8:e9817

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We have used this protocol at ten different schools in Georgia, USA (shared it with over 946 students) and at two universities (30 undergraduates) and it works well, if you follow the detailed instructions.

Created: May 20, 2020

Last Modified: June 13, 2020

Protocol Integer ID: 37336

Keywords: *Chlamydomonas*, micro-alga, algal beads, photosynthesis, cellular respiration, Plant Biology education, K16 students, Biology, algae bracelets, algae necklace, bioenergy, biotechnology, biofuel, algae bead, cellular respiration under light, photosynthesis in the algal cell, algal bead bracelet, alkaline in the light, *chlamydomonas* bead, photosynthesis convert, algal cell, interplay of photosynthesis, water inside the bracelet, algae, pH testing strip, respiration in live cell, cellular respiration to k4, useful for bioenergy researcher, using pH testing strip, acidic in the dark, relative rates of photosynthesis, pH change, powered *chlamydomona*, photosynthetic mutant strain, respiration, live cell, pH of the water, making inexpensive light, college biology student, alkaline, inexpensive light, releases CO₂ into the environment, bioenergy researcher, green micro, bead, lab activity, light condition, CO₂ from the water, bracelet, necklace, pH, access to an oxygen electrode, *chlamydomona*, carbon in the presence, O₂,

Abstract

Scientific concept: Cellular respiration oxidizes organic chemicals and releases CO₂ into the environment and photosynthesis converts CO₂ into fixed carbon in the presence of light. Irrespective of the light conditions, respiration in live cells in the beads will release CO₂ that will dissolve in water inside the bracelet/necklace in which the beads are immersed to generate carbonic acid. Conversely in the light, photosynthesis in the algal cells in the beads will remove CO₂ from the water surrounding the beads. Hence pH of the water inside the algal bead bracelets will be acidic in the dark and alkaline in the light. Students can monitor the dramatic color changes in a gradient fashion in the water inside the bracelet that is induced by the relative rates of photosynthesis/cellular respiration. Students can measure pH of the water inside the bracelets/necklaces using pH testing strips. This protocol might be also useful for bioenergy researchers who are interested in exploring the use of immobilized motile micro-algae for biofuel production.

Applications: This educational protocol can be used to teach the interplay of photosynthesis and cellular respiration to K4 - college Biology students in a fun and engaging way. The lab activities can be customized according to grade levels. Students will make their *Chlamydomonas* beads and use these beads to make algae bead bracelets or necklaces to test the opposing nature of chemical reactions of cellular respiration and photosynthesis. This interplay can be studied by visually monitoring the water color/pH changes induced by relative rates of photosynthesis and cellular respiration under light and in darkness. Alternatively, this protocol can be used by educators at institutions with limited resources who do not have access to an oxygen electrode to comparatively study relative rates of photosynthesis in *Chlamydomonas* wild type and a photosynthetic mutant strain using strain specific-algae beads that approximately the same number of cells per bead.

Attachments



[Article S2.pdf](#)

25KB

Guidelines

Detailed guidelines and protocol related warnings have been included in each step of the described protocol.

Materials

Materials/equipment required:

1. 200-300 ml of dense healthy culture of *Chlamydomonas reinhardtii* wild type strain (4A+ [CC-4051 4A+ mt+] or CC4533 [CC-4533cw15mt-])
2. TAP or HSA growth media for growing *Chlamydomonas*
3. YA slants or TAP agar media plates for maintaining *Chlamydomonas* strain stocks on agar media.
4. 100 mL of 2% sodium alginate solution
5. Pre-chilled 600 mL of 3% CaCl_2
6. Graduated measuring cylinders
7. P1000 micropipette (**Note:** you will need also a P200 micropipette if you conduct the experiment in a glass vial)
8. Micropipette tips (that can be fitted to P1000 and P200 micropipettes)
9. 50 mL Falcon tubes
10. Tea strainer or oil strainer or any strainer with fine mesh or coffee filters
11. Plastic spoons
12. Plastic transfer pipettes
13. Tap water
14. 500 mL or 1000 mL beakers
15. Screw capped 500 ml or 1000 mL bottles for chemical solution storage
16. pH test strips
17. Bicarbonate indicator solution
18. Petri dishes
19. Funnels




20. Magnetic stirrer and stir bars
21. Alumina foil to cover the mouth of the algae culture flasks/dark exposure
22. Flexible tubing
23. Colorful yarn from arts/craft stores like Hobby Lobby
24. 1 L and 50 mL glass Erlenmeyer flask to culture algae
25. Inoculating loops for inoculating algae media
26. Ruler
27. Scissor
28. 1.5 mL Eppendorf tubes
29. 5.5 ml screw capped glass vials (if one wants to perform the same experiment in vials instead of in a bracelet)
30. Hausser Scientific Bright-Line™ Counting Chamber
31. Phenol red indicator solution

Additionally, the lab requires few simple equipment like bench top centrifuge, an orbital open- air shaker for shaking and aerating algae culture flasks, a weighing balance for weighing chemicals, a refrigerator for storing CaCl_2 solution and a Bunsen burner for flaming the mouth of algae culture flask during media inoculation, as part of sterile technique. pH meter would be required to accurately measure pH in glass vials but is not required for algae bead bracelets as pH strips are used to measure pH of the water inside the bracelets. Cool white fluorescent lights are required to grow algae. To simulate darkness, bracelets or vials can be kept inside a cabinet drawer; alternatively, one can keep them covered with alumina foil on the bench) The color change of the algae bead bracelets/necklaces and glass vials can be photographed with a cell phone camera.

Troubleshooting

Safety warnings

 There are no safety warnings associated with this protocol.

Before start

Read every step of the protocol very carefully. Please feel to email me: mmitra@westga.edu for any question that you might have.

1 Algal media and cultures:

Maintain *Chlamydomonas* wild type strain 4A+ strain in the lab on Tris-Acetate Phosphate (TAP) agar media plates [17] in dim light intensities ($15\text{--}20\ \mu\text{mol m}^{-2}\text{s}^{-1}$) at 25°C . Start a starter culture of 4A+ approximately 10–12 days before the lab activity by inoculating 10 mL of liquid TAP media in a 50 mL flask with 4A+ cells from a 5-day old TAP agar media plate of respective strains. After 5 days of growth, add 1 mL of the starter culture to inoculate 300 mL of fresh TAP media in a 1L flask. Grow the TAP liquid 4A+ culture for 6–7 days for dense dark green growth. We have found that the *Chlamydomonas* culture should be grown under low light ($80\text{--}100\ \mu\text{mol photons m}^{-2}\text{s}^{-1}$) to obtain a healthy culture that is not photo-oxidatively stressed to be used for our lab activities. Grow *Chlamydomonas* TAP liquid culture at 25°C under continuous illumination of $100\ \mu\text{mol photons m}^{-2}\text{s}^{-1}$, provided by the combined light intensities of four to six cool white fluorescent lights on an open-air orbital shaker at a speed of 150–180 rpm.

2 Preparation of 2% sodium alginate and 3% Calcium chloride solutions 2 days before bead making:

Dissolve 2 grams of sodium alginate in 100 mL of E-pure water overnight at room temperature by stirring at a speed of 500–600 rpm using a magnetic stirrer. [Note: sodium alginate forms a very viscous solution when dissolved at 1.5%–4%]. 2% sodium alginate solution was stored at room temperature. Dissolve thirty grams of calcium chloride in 1000 mL of E-pure water and store at 4°C in a fridge.

3 Cell counts:

Determine cell density before harvesting *Chlamydomonas* cells from the TAP liquid culture to estimate the volume of culture that needs to be harvested per 50 mL falcon tube. Calculate cell density (number of cells per mL of the culture) by counting the cells using a Hausser Scientific Bright-Line™ Counting Chamber (Hausser Scientific, Philadelphia, PA). Please calculate the cell density of the culture to determine the volume of culture you need to spin down in a benchtop centrifuge to harvest $250\text{--}395 \times 10^6$ cells.

4 Preparation of *Chlamydomonas* 4A+ beads:

Harvest *Chlamydomonas* cells by spinning down the dense TAP liquid 4A+ culture in 50 mL falcon tubes at 1000–1500 g for 3 minutes in a benchtop centrifuge. Do not perform a long spin (10–15 minutes) as it generates a hard cell pellet which takes a long time to resuspend in 2% sodium alginate. Discard the supernatant and collect the cell pellet. Maximum number of total cells harvested per 50 mL falcon tube is 395×10^6 cells/50 mL falcon tube and should not be exceeded (Best result is obtained when the harvested cell

numbers are between 250×10^6 - 395×10^6 cells/ 50 mL falcon tube). Note: Harvesting 100 mL of dense *Chlamydomonas* culture having cell density of 21×10^6 cells/mL should generate 200-300 beads of 4 - 4.5 mm in diameter. The culture volume to spin down for harvesting cells depends on the cell density and the number of beads desired. To a cell pellet volume of 1.25- 1.5 mL, add 4-5 mL of 2% sodium alginate using a P1000 micropipette to get a final volume of 6-7 mL algae cell pellet-sodium alginate suspension. [Crude ratio: 2% well mixed-sodium alginate is added to the cell pellet in a 5:1 or 4:1 ratio depending on the total number of cells harvested]. Mix the algae cell pellet and sodium alginate gently to obtain a uniform mixture by gently swirling the falcon tubes for 5-10 minutes or till the entire cell pellet is completely resuspended without any visible cell clumps. Vortexing of the cell pellet tubes should be avoided during resuspension of cells in sodium alginate.

Add the algae-sodium alginate mix drop wise steadily and quickly with uniform pipetting motion by using a P1000 micropipette into a beaker of pre-chilled 3% CaCl_2 kept on ice. [Note: Bead sizes vary depending on pipetting techniques and type of micropipette or plastic transfer pipette used]. Air bubbles create problems when plastic transfer pipettes are used but they are suitable for use in K4-K12 classrooms. If pipetting is not smooth and regular and algae-sodium alginate mixture is not mixed by swirling in between pipetting, irregular shaped beads or beads with different cell numbers/bead (light and dark green beads) will form. Avoid holding the micropipette pipette tip very close to the CaCl_2 liquid surface to prevent clogging of the pipette tip by calcium alginate because of the instantaneous reaction of sodium alginate with CaCl_2 . (Notes: DO NOT DISTURB the CaCl_2 beaker while you are adding the drops. Make sure that you do not have any air gap or bubbles in your pipette tip while you pipette the algae-sodium alginate mixture). While pipetting dropwise the algae-sodium alginate suspension in CaCl_2 solution, the beaker of chilled CaCl_2 should not be shaken and should be given a very gentle swirl in between one batch of pipetting. This swirling mixes the CaCl_2 solution around the beads that have formed in the CaCl_2 beaker to harden them. Pipetting can be resumed only when CaCl_2 liquid surface is still after the swirling of the CaCl_2 solution. As soon as the algae-sodium alginate mixture drops hit the chilled CaCl_2 liquid surface, they solidify into tiny beads. Continue making beads until all the mixture is used. Keep the CaCl_2 beaker containing the beads on ice for 10-15 minutes to allow complete solidification of the algal beads.

Separate the algal beads from the CaCl_2 solution by filtering through a strainer/coffee filter and collect the filtered CaCl_2 in a bottle using a funnel. This filtered CaCl_2 can be stored in a fridge at 4°C and can be reused three- five times for making algal beads in future. Algal beads should be rinsed gently and thoroughly with tap water for at least five minutes to remove residual sodium chloride which is formed when sodium alginate reacts with calcium chloride to form calcium alginate and sodium chloride. This step is a very important step and must not be skipped as any residual sodium chloride will hinder photosynthesis in the experiment. Use a plastic spoon to scoop out the beads gently

from the strainer to put them in a petri dish containing small amount of tap water to prevent drying of beads while you make the bracelets/necklaces (see below).

Note: Any leftover algae beads can be stored in tap water in the refrigerator for 2-3 days for future experiments. Algae bead making demonstration video clips are available at: <https://youtu.be/u4BbZ29qIWQ> and at <https://youtu.be/elxbzeHW8IM>.

5 Preparation of Chlamydomonas 4A+ bead bracelets and necklaces:

Cut flexible tubing into 10 pieces, each 5 inches long for bracelets and 6 inches long for necklaces. Cut off caps of 1.5 mL Eppendorf tubes with a scissor. Use de-capped Eppendorf tubes to plug the ends of the bracelet tubing (one de-capped tube at each end of the cut tubing). Cut colorful yarns according to wrist width, intertwine them to make braids. Loop one braid tightly onto the mouth of each de-capped Eppendorf tube at each end of the bracelet. Cut plastic transfer pipette tip ends at an angle to form scoops. Next, unplug one end of the bracelet tubing by removing one of the de-capped Eppendorf tube that you used for sealing the tubing end. Add about 3-5 mL of tap water [pH 6.9 -7] (depending on the length of the bracelet tubing ranging from 5 inches to 7 inches) into the bracelet flexible tubing either using a micropipette or a plastic transfer pipette. De-ionized/E-pure water is known to contain less dissolved gases and minerals than tap water. Amounts of dissolved oxygen and carbon-dioxide in the water used for photosynthesis monitoring experiment will affect the results in a photosynthesis lab. We have found that tap water works best for the photosynthesis experiment. It is important to have about 0.5 cm-1 cm air gap at each end inside the tubing to provide enough air for cells to breathe as the cells are oxygen-stressed inside the beads. To prepare 10 algae bead bracelets, each 5 inches long and containing 10-12 beads, you would need to make approximately 100-120 beads. Always make extra beads and store them for future use for 2-5 days in the fridge. Add 8-9 drops of the bicarbonate indicator solution to 3.5 ml- 4 mL of tap water inside the bracelets. For algae bead necklaces each of 6-7 inches long, you will need 20-35 beads per necklace. Add 10-11 drops of the pH indicator bicarbonate indicator solution to 4.5 ml- 5.5 mL of tap water inside the algal bead necklaces. Algae beads should be scooped from the petri dish using the transfer pipettes which have been converted to a scoop. Gently introduce the algal beads using the scoop into the tap water inside the tubing without squeezing the bulb of the transfer pipette and just let the water flow inside the tubing get the beads from the immersed scoop. Plug the end of the tubing back with the de-capped Eppendorf tube. Give a gentle swish to the tubing to mix the water and the bicarbonate indicator solution inside the tubing. Check the sealing at the tubing ends to ensure the tubing ends were sealed tightly and there is no water leakage. Plugging must be tight but not so tight so that it makes it hard for you to take out the Eppendorf tubes to flush the bracelet tubing and reuse it in future. Hence do not push the Eppendorf tubes all the way into the tubing.

6 Dark and light exposures of algae bead bracelets/necklaces:

Image the bracelet and record the pH of the water before shifting it to light. The bracelet can be kept under 150-200 micromole $m^{-2}s^{-1}$ light intensity [equivalent to the combined light intensities of 12 to 14 cool white fluorescent lights] or under cool LED indoor plant grow light. But to our experience, the bracelet takes longer time to change color when kept under one LED indoor plant grow light because the light intensity is not high enough compared to what one gets with cool white fluorescent lights. You might need 2-3 indoor plant grow light to speed up the color change in the bracelets (we have not tried it). After light exposure for 2-4 hours (depending on the length of the bracelet tubing and number of algae beads inside the tubing), the bracelet will change from light red/magenta to a dark blue/purple color via gradient color changes because of an increase in pH induced by photosynthesis in light. When shifted to dark for 3-4 hours the bracelet color will change from dark blue/purple to bright yellow via gradient color changes because of a decrease in pH induced by cellular respiration and absence of photosynthesis. You can switch the bracelet back forth between light and dark for 2-3 times over a period of 2-3 days and it should show color changes. Demonstration of algae bead bracelet making video clips are available at: <https://youtu.be/A7VIjLDGSCc> and https://youtu.be/vh_1ASpQgS8 and <https://youtu.be/enctr0yhWQ8>