


Jun 06, 2019

Mammalian Cell Culture: Freezing

DOI

<https://dx.doi.org/10.17504/protocols.io.3szgnf6>



Kenneth Schackart¹

¹University of Arizona

Yoon Lab



Kenneth Schackart

University of Arizona

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.3szgnf6>

Protocol Citation: Kenneth Schackart 2019. Mammalian Cell Culture: Freezing. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.3szgnf6>



License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: June 06, 2019

Last Modified: June 06, 2019

Protocol Integer ID: 24121

Keywords: mammalian cell culture, cell culture, tissue culture flask, cell

Abstract

This protocol describes how to freeze cells that are being cultured in a tissue culture flask (T-75 or T-25).

Guidelines

- Gloves must be worn at all times.
- Perform all tasks within biosafety cabinet.
- Anything entering biosafety cabinet must be generously sprayed with 70% ethanol (even you).
- When finished, wipe biosafety cabinet with 70% ethanol, and UV for at least 15 minutes.

Materials

- Cell culture media (at least 4 mL for T-75 or 1 mL for T-25)
- Fetal Bovine Serum (FBS)
- DPBS
- Trypsin-EDTA
- Dimethylsulfoxide (DMSO)
- 15 mL centrifuge tube
- Cryovial(s)

Before start

- Warm cell culture media, DPBS, and Trypsin-EDTA (and possibly FBS) in  37 °C water bath.
- Wash waste beaker with soap and warm water, then dry with paper towel.
- Expose serological pipet tips, centrifuge tube, and waste beaker to UV for at least  00:15:00 .
- Make label for cryovial with cell line, passage number, and date; affix to cryovial before UV.



Wash Cells

- 1 Remove media from flask.
- 2 Using serological pipet, add  4 mL DPBS to flask [ 1 mL for T-25].
- 3 Using serological pipet, remove DPBS and dispose into waste beaker.
- 4 Repeat the above 2 steps, for a total of 2 washes.

Trypsinize

- 5 Add  4 mL warmed Trypsin-EDTA to flask [ 1 mL for T-25].
- 6 Wait for cells to detach.
- 7 Add  4 mL warmed cell culture media to flask to neutralize Trypsin-EDTA [ 1 mL for T-25]

Pellet Cells

- 8 Using a serological pipet, transfer cell suspension into 15 mL centrifuge tube.
- 9 Centrifuge cell suspension at  1500 rpm for  00:03:00 .

Resuspend



- 10 Remove bulk of supernatant with serological pipet.
- 11 Remove remaining supernatant with 1000 μL pipette. For small cell pellets, it is best to leave a little supernatant to avoid disturbing the pellet.
- 12 Add  950 μL of either FBS or cell culture media and allow to sit for  00:01:00 to make resuspension easier.

Note

Some researchers prefer to use FBS as a freezing medium while others prefer whole media.

- 13 Gently pipette mix cell pellet until resuspended.

- 14 Add  50 μL DMSO to cell suspension.

Note

5% DMSO solution prevents ice crystals from forming in the liquid, minimizing cell rupture during freezing.

Freeze

- 15 Place cryovial in deep freezer cell vial container.
- 16 Allow at least  04:00:00 in deep freezer before transferring cryovial to liquid nitrogen tank.

Note

This prevents too rapidly freezing, which may cause cell rupture.



Document

- 17 Update Lab Frozen Storage Inventory to reflect the new cell vial.