


Feb 06, 2019

Mammalian Cell Culture: Subculturing

DOI

dx.doi.org/10.17504/protocols.io.xukfnuw



Kenneth Schackart¹, Kattika Kaarj¹

¹University of Arizona

481b Laboratory



Kenneth Schackart

University of Arizona

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.xukfnuw>

Protocol Citation: Kenneth Schackart, Kattika Kaarj 2019. Mammalian Cell Culture: Subculturing. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.xukfnuw>



License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: February 06, 2019

Last Modified: February 06, 2019

Protocol Integer ID: 20076

Keywords: mammalian cell culture, cell culture, confluent mammalian cell, mammalian cell, cell, mammalian, culture

Abstract

This protocol details how to subculture nearly confluent mammalian cells grown in a T-25 flask.

Materials

- Gloves
- 0.05% or 0.25% warmed Trypsin-EDTA
- Warmed cell culture Media (e.g. DMEM:F12, EMEM)
- DPBS
- 15 mL centrifuge tube
- Serological pipette and tips
- 1000 µL pipette and tips

Troubleshooting

Safety warnings


! Lab coat and gloves must be worn at all times.

Assess Cell Confluency


- 1 Under light microscope, look at the cells and assess level of confluency. This is how you will determine the need to subculture.

Note

Confluency can be estimated by evaluating the percentage of surface covered by cells.




Low confluency SH-SY5Y



High confluency SH-SY5Y

Wash Cells



2 Using serological pipette, add  1 mL DPBS to T-25 flask.


3 Using serological pipette, remove DPBS and dispose into waste beaker.

4 Repeat the above 2 steps, so that you will wash the cells twice.

Note

Always use a fresh pipette tip when drawing liquid from a stock solution.

Trypsinize

5 Add  1 mL warmed trypsin-EDTA to T-25 flask.

6 Wait  00:05:00 for trypsin-EDTA to detach the cells.

Note

This time will vary in practice, and depends on cell type and trypsin concentration (i.e. 0.05% vs 0.25%). Some cell types may take up to 15 minutes. In those cases, assess detachment progress using a light microscope.



7 Add  1 mL cell culture media.

Note

Trypsin-EDTA is neutralized by adding a volume of cell culture media equal to that of trypsin-EDTA.


Spin Down

8 Using a serological pipette, transfer the cell suspension (cells, trypsin-EDTA, and cell culture media) into a 15 mL centrifuge tube.

- 9 Add  3.5 mL fresh cell culture media to T-25 flask, this will preserve any remaining cells.
- 10 Centrifuge the cell suspension on 1.5 kRPM for  00:03:00 .

Note

This has been performed by T.A. ahead of time. Bring your cell suspension to trade for a spun down cell pellet.



Centrifuge tubes with large cell pellets at the bottom and 8 mL of supernatant.

Resuspend and Reseed

- 11 Remove supernatant, dispose into waste beaker.



**Note**

You can leave a small amount with the serological pipette, the rest will be taken off in the next step.

- 12 Using a 1000 μL pipette, carefully remove the remaining supernatant, being cautious not to disturb the cell pellet.

Safety information

Always dispose of pipette tips in sharps container. Do not use the same tip twice.

- 13 Add  1000 μL cell culture media to the cell pellet, and allow to sit for  00:01:00 .


- 14 Gently pipette mix the cell pellet until the pellet is resuspended.

Note

Pipette mixing is done by slowly drawing in solution and pushing it out several times, all without removing the pipette tip from the solution. Ask T.A. for help on this if you need some pointers.

- 15 Reseed your T-25 flask by transferring  500 μL cell suspension to the flask.

Note

In practice you would reseed the other  500 μL into another flask, or even split the cell suspension into 3 or 4 flasks. Alternatively, you make take the remaining cell suspension for an experiment.

Incubate

- 16 Incubate at  37 $^{\circ}\text{C}$ in CO_2 incubator.

