

Jun 06, 2019

## Mammalian Cell Staining

 Forked from [Mammalian Cell Staining](#)

DOI

<https://dx.doi.org/10.17504/protocols.io.x97fr9n>

Kenneth Schackart<sup>1</sup>


<sup>1</sup>University of Arizona

Yoon Lab



Kenneth Schackart

University of Arizona



### Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.x97fr9n>

**Protocol Citation:** Kenneth Schackart 2019. Mammalian Cell Staining. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.x97fr9n>



**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** February 18, 2019

**Last Modified:** June 06, 2019

**Protocol Integer ID:** 20511

**Keywords:** mammalian cell, cell, actin filament, focal adhesion site, well plate

## Abstract

This protocol details how to stain mammalian cells cultured on a 96-well plate. Actin filaments, focal adhesion sites (as indicated by the presence of vinculin), and nuclei will be stained.

## Guidelines

- Gloves must be worn at all times.
- Perform all tasks within biosafety cabinet.
- Anything entering biosafety cabinet must be generously sprayed with 70% ethanol (even you).
- When finished, wipe biosafety cabinet with 70% ethanol, and UV for at least 15 minutes.

## Materials

- 4% Paraformaldehyde solution
- 0.1% Triton X-100 solution in PBS
- Blocking buffer (PBS + 1% bovine serum albumin)
- Anti-vinculin solution (1:500 in blocking buffer)
- TRITC and FITC-conjugated secondary antibody solution (1:1:248, TRITC:FITC:blocking buffer) referred to as FITC:TRITC, prepared in opaque centrifuge tube
- DAPI solution (1:249 DAPI:blockingbuffer) prepared in opaque centrifuge tube
- Phosphate buffered saline (PBS)
- Washing buffer (PBS + 0.05% Tween-20)

## Before start

- Warm DPBS in  37 °C water bath.
- Wash waste beaker with soap and warm water, then dry with paper towel.
- Expose waste beaker to UV for at least  00:15:00 .



## Fix the cells

- 1 Remove cell culture media.

### Note

You may wash twice with DPBS at this point.

- 2 Add  100  $\mu$ L of  4 % volume paraformaldehyde solution.

- 3 Incubate for  00:15:00 .

## Perforate cell membrane

- 4 Remove paraformaldehyde solution.
- 5 Wash twice with  100  $\mu$ L washing buffer.
- 6 Add  100  $\mu$ L of  0.1 % volume Triton X-100.
- 7 Incubate for  00:05:00 .

## Block unspecific binding

- 8 Remove Triton X-100.
- 9 Wash twice with washing buffer.



10 Add  100  $\mu$ L blocking buffer.

11 Incubate for  00:30:00 .

## Bind anti-vinculin to vinculin

12 Remove blocking buffer.

13 Wash twice with washing buffer.

14 Add  100  $\mu$ L of Anti-Vinculin and blocking buffer mixture.

15 Incubate for  01:00:00 .

### Note

Increasing this duration can greatly increase FITC-conjugated secondary antibody binding. It is acceptable to place covered plate in refrigerator overnight. While waiting on this incubation, it is a good time to prepare the combined dye solution.

## Stain actin filaments and focal adhesion sites

16 Remove anti-vinculin blocking buffer mixture.

17 Wash twice.

18 Add  100  $\mu$ L FITC:TRITC solution



19 Incubate for  00:30:00 .

## Stain nucleus

20 Remove stains.

21 Wash twice with washing buffer.

22 Add  100  $\mu$ L of DAPI solution

23 Incubate for  00:05:00

24 Remove DAPI.

25 Add  100  $\mu$ L PBS to prevent dessication of cells during imaging/storage.