


Feb 05, 2019 Version 2

Mammalian Cell Staining V.2

DOI

dx.doi.org/10.17504/protocols.io.xs7fnhn



Kenneth Schackart¹, Kattika Kaarj¹

¹University of Arizona

481b Laboratory



Kenneth Schackart

University of Arizona

OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.xs7fnhn

Protocol Citation: Kenneth Schackart, Kattika Kaarj 2019. Mammalian Cell Staining. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.xs7fnhn>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

Created: February 05, 2019

Last Modified: February 05, 2019

Protocol Integer ID: 20031

Abstract



This protocol details how to stain mammalian cells cultured on a 96-well plate. Actin filaments, focal adhesion sites (as indicated by the presence vinculin), and nuclei will be stained.






Materials

- 4% Paraformaldehyde solution
- 0.1% Triton X-100
- Blocking buffer
- Anti-vinculin with blocking buffer
- TRITC and FITC-conjugated secondary antibody solution
- NucBlue™ LiveReady Probes™ Reagent solution
- Phosphate buffered saline (PBS)
- Dulbecco's phosphate buffered saline (DPBS)


Fix the cells

- 1 Remove cell culture media.
- 2 Add  100 μ L of [M] 4 % volume paraformaldehyde solution.
- 3 Incubate for  00:05:00 .


Perforate cell membrane

- 4 Remove paraformaldehyde solution.
- 5 Wash twice with  100 μ L DPBS.
- 6 Add  100 μ L of [M] 0.1 % volume Triton X-100.
- 7 Incubate for  00:05:00 .

Block Unspecific Binding

- 8 Remove Triton X-100.
- 9 Wash twice with DPBS.
- 10 Add  100 μ L blocking buffer.




11 Incubate for  00:10:00 .


Stain for focal adhesion sites and actin filaments


12 Remove blocking buffer.

13 Wash twice with DPBS.

14 Add  250 μL of Anti-Vinculin and blocking buffer mixture.

15 Incubate for  00:20:00 .

16 Add  100 μL FITC-conjugated secondary antibody and TRITC.

17 Incubate for  00:30:00 .

Stain nucleus

18 Remove stains.

19 Add  100 μL of NucBlue™ solution.