


Nov 28, 2019

🌐 Marchantia Cryopreservation of Gemmae

🔗 Forked from [Marchantia cryopreservation of gemmae](#)

DOI

dx.doi.org/10.17504/protocols.io.9vxh67n



Marius Rebmann¹, Marta Marta Tomaselli², Linda Silvestri³, Susana Sauret-Gueto¹, Michelle Lim⁴

¹Plant Sciences, University of Cambridge, OpenPlant; ²University of Cambridge, Open Plant;

³University of Cambridge; ⁴Plant Sciences, University of Cambridge, OpenPlant,

OpenPlant Project



Marius Rebmann

Plant Sciences, University of Cambridge, OpenPlant

OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.9vxh67n

Protocol Citation: Marius Rebmann, Marta Marta Tomaselli, Linda Silvestri, Susana Sauret-Gueto, Michelle Lim 2019. Marchantia Cryopreservation of Gemmae. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.9vxh67n>

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: November 28, 2019

Last Modified: November 28, 2019

Protocol Integer ID: 30359

Keywords: Cryopreservation, Marchantia, Gemmae



Abstract

Simplified cryopreservation protocol for *Marchantia polymorpha* gemmae, based on ([Tanaka et al. 2016, Plant and Cell Physiology](#)). Enables long term storage of gemmae at -80°C.



Materials

MATERIALS

✕ Agar **Melford Catalog #A20021**

✕ Absciscic acid **Merck MilliporeSigma (Sigma-Aldrich) Catalog #A1049**

✕ Sodium Alginate **Duchefa Biochemie Catalog #S1320**

✕ Calcium chloride **VWR International (Avantor) Catalog #22328.262**

✕ Sucrose **Fisher Scientific Catalog #10634932**

✕ Glycerol **Merck MilliporeSigma (Sigma-Aldrich) Catalog #G5516**

✕ Gamborg B5 Medium with Vitamins **Duchefa Biochemie Catalog #G0210**

- Petridish, 9cm/4.5cm or 6/12-well transparent multi-well plates
- Liquid nitrogen (N₂)

Reagent setup:

Preculture plates

1/2 strength Gamborg B5, pH 5.8 + 1.2% agar + 0.3 M sucrose + 10μM ABA

Autoclave media without ABA, add ABA from a filter sterilised stock solution just before pouring plates.

9cm/4.5cm petridishes or 6/12-well transparent multi-well plates can be used for plates, depending on desired throughput.

Alginate Solution

1/2 strength Gamborg B5 + 3% Sodium alginate

CaCl₂ Solution

1/2 strength Gamborg B5 + 0.1M CaCl₂

Loading solution

1/2 strength Gamborg B5 + 2M glycerol, + 1M sucrose

Thawing solution

1/2 strength Gamborg B5 + 1M sucrose

Rinse solution


1/2 strength Gamborg B5 + 1% sucrose

Preculture

- 1 Collect fresh gemmae and plate on **preculture plates**, incubate 1-3 days under normal growth conditions (e.g. 21°C constant light)


Encapsulation

- 2 Place small drops (approx. 40µL) of **alginate solution** on an empty 4.5cm petridish/multi-well plate to form beads.




- 3 Use forceps to transfer gemmae from preculture plates, place 1-5 gemmae inside each bead.

Proceed quickly to step 4 to avoid gemmae sinking to the bottom or edge of the bead.




- 4 Add a drop of **CaCl₂ solution** to each bead, make sure to not displace gemmae from the bead.

- 4.1 Allow beads to solidify for 10min.



- 5 Use a serological pipette to fill the plate with **loading solution**, submerging the beads. Soak beads for 30min.
- 6 Use a serological pipette to remove the loading solution without disrupting the beads.
- 6.1 Remove lid and air dry beads in flowhood for >2h.



- 7 Use forceps or scalpels to transfer beads into 1.5mL Eppendorf tubes (max. 10 beads/tube).
- 7.1 Flash freeze tubes in liquid nitrogen for >2 minutes.
- 7.2 Move tubes to -80°C for long term storage.

Thawing and Recovery

- 8 Remove tubes from freezer and immediately place in a 37°C water bath or heat block for 2 minutes.
- 9 Add 1.5mL **thawing solution** to each tube, incubate at room temperature for 10 minutes.



- 10 Replace thawing solution with **rinse solution**. Incubate at room temperature for 10 minutes.
- 11 Remove rinse solution, transfer beads on 1/2 Gamborg plates, incubate under normal growth conditions (e.g. 21°C constant light) until thalli grow