


Sep 24, 2019

Marchantia genotyping (quick and dirty genomic DNA extraction)

DOI

dx.doi.org/10.17504/protocols.io.4wagxae



Eftychis Frangedakis¹, Marta Marta Tomaselli², Marius Rebmann³, Susana Sauret-Gueto³

¹University of Cambridge; ²University of Cambridge, Open Plant;

³Plant Sciences, University of Cambridge, OpenPlant

OpenPlant Project



Eftychis Frangedakis

University of Cambridge, Plant Sciences

OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.4wagxae

Protocol Citation: Eftychis Frangedakis, Marta Marta Tomaselli, Marius Rebmann, Susana Sauret-Gueto 2019. Marchantia genotyping (quick and dirty genomic DNA extraction). **protocols.io** <https://dx.doi.org/10.17504/protocols.io.4wagxae>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: June 28, 2019

Last Modified: September 24, 2019

Protocol Integer ID: 25250

Abstract

This protocol allows for quick and dirty genomic DNA extraction. It can easily be used for genotyping with PCR. The quality of the genomic DNA extracted is not suitable for any other application.

Materials

MATERIALS

 KOD Hot Start DNA Polymerase **Merck MilliporeSigma (Sigma-Aldrich) Catalog #71086-3**



- 1 Take small pieces (3×3 mm) of thalli from individual plants and place in a 1.5 mL Eppendorf tube (A in Figure).
- 2 Add 100 µl genotyping buffer (B in Figure).
- 3 Crush with an autoclaved micro-pestle (C in Figure).
- 4 Place the tube(s) at 80 °C for 10 min.
- 5 Add 380 µL of sterile water to each tube (D in Figure).
- 6 Use 5 µl aliquot of the extract as a template for PCR using preferably the KOD Hot start polymerase.

Note

We found KOD to be more reliable amplifying fragments from a crude genomic DNA extract such as the one used here.

- 7 Check PCR products on a 1.5% (w/v) agarose gel.

8

