


Sep 24, 2019

Marchantia protoplast

DOI

dx.doi.org/10.17504/protocols.io.5kag4se



Eftychis Frangedakis¹

¹University of Cambridge

OpenPlant Project



Eftychis Frangedakis

University of Cambridge, Plant Sciences

OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.5kag4se

Protocol Citation: Eftychis Frangedakis 2019. Marchantia protoplast. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.5kag4se>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: July 17, 2019

Last Modified: September 24, 2019

Protocol Integer ID: 25954

Materials

8% (w/v) mannitol at pH 5.7

- 1 Grow gemmae on ½ Gamborg B5 plus vitamins media plates (the plates should be covered like Fig. 1A) for 4 days (Fig. 1B).
- 2 Add 10 ml of 8% (w/v) mannitol at pH 5.7 on the plate and incubate at room temperature 30min.
- 3 Prepare 2% Driselase solution: add 0.2 g Driselase in 10 ml 8% mannitol solution and incubate in dark for 20 min at RT. Gently invert to mix at intervals. Spin at 3.3k for 3 min. Filter sterilize supernatant using a 0.2 µm filter attached to a syringe.
- 4 Replace mannitol solution with 5 ml 2% Driselase and incubate at RT for 4 hours with gentle shaking (50 rpm).
- 5 Using a tip gently release protoblasts from gemmae (Fig.1C).
- 6 Filter the solution through a 70 µm cell strainer into a 50 ml falcon tube and then transfer in a round bottom falcon tube (Fig. 1D).
- 7 Incubate for further 10 min.
- 8 Spin at 120 x g for 3min. Remove the supernatant without disturbing the cells.
- 9 Wash the cells by re-suspending them in 6ml 8% mannitol (Mix by swirling not shaking) and spin at 120 x g for 3 min. Repeat the wash 2 times.
- 10 Add 300 µl 50% glucose in the autoclaved liquid Gamborg B5 plus vitamins + mannitol solution

11

