



Jan 08, 2026

Mass spectrometry analyses

DOI

<https://dx.doi.org/10.17504/protocols.io.eq2ly579rvx9/v1>

Junyu Xiang¹

¹Department of Gastroenterology, Chongqing Key Laboratory of Digestive Malignancies, Daping Hospital, Chongqing 400042, China



Junyu Xiang

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.eq2ly579rvx9/v1>

Protocol Citation: Junyu Xiang 2026. Mass spectrometry analyses. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.eq2ly579rvx9/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: January 08, 2026

Last Modified: January 08, 2026

Protocol Integer ID: 238195

Keywords: ms metabolomic, metabolomic, mass spectrometry, detailed description for mass spectrometry, based quantitative proteomic, quantitative proteomic, proteins screening, ms for isotype, interacted proteins screening, m

Abstract

A detailed description for mass spectrometry-based quantitative proteomics, interacted proteins screening, LC-MS metabolomics and GC-MS for isotype-labeled metabolomics.

Quantitative proteomics assay

- 1 Cells are collected and sonicated three times on ice using a high intensity ultrasonic processor (Scientz) in lysis buffer (8 M urea, 1% protease inhibitor cocktail).
- 2 The remaining debris is removed by centrifugation at 12,000g at 4°C for 10 m.
- 3 Supernatant is collected, and the protein concentration is determined with BCA kit according to the manufacturer's instructions.
- 4 The protein solution is reduced with 5 mM dithiothreitol for 30 m at 56°C and alkylated with 11 mM iodoacetamide for 15 m at room temperature in the dark.
- 5 Protein samples are diluted by adding 100 mM TEAB to a urea concentration less than 2 M.
- 6 Trypsin is added at 1:50 trypsin-to-protein mass ratio for the first digestion overnight and 1:100 trypsin-to-protein mass ratio for a second 4 h digestion.
- 7 The peptides are desalted by C18 SPE column and dissolved in solvent A (0.1% formic acid, 2% acetonitrile/in water), directly loaded onto a home-made reversed-phase analytical column (25-cm length, 75/100 µm ID).
- 8 The peptides are separated with a gradient from 6% to 24% solvent B (0.1% formic acid in acetonitrile) over 70 m, 24% to 35% solvent B over 14 m, and 35% to 80% over 3 m then holding at 80% for 3 m, all at a constant flow rate of 450 nL/m on a nanoElute UHPLC system (Bruker Daltonics).
- 9 The peptides are subjected to capillary source followed by the timsTOF Pro (Bruker Daltonics) mass spectrometry.
- 10 The electrospray voltage applied is 1.60 kV.
- 11 Precursors and fragments are analyzed at the TOF detector, with a MS/MS scan range from 100 to 1700 m/z.
- 12 The timsTOF Pro is operated in parallel accumulation serial fragmentation (PASEF) mode.





- 13 Precursors with charge states 0 to 5 are selected for fragmentation, and 10 PASEF-MS/MS scans are acquired per cycle.
- 14 The dynamic exclusion is set to 30 sec.
- 15 The resulting MS/MS data are processed using MaxQuant search engine (v.1.6.15.0). 
- 16 Tandem mass spectra are searched against the human SwissProt database (20422 entries) concatenated with reverse decoy database. 
- 17 Trypsin/P is specified as cleavage enzyme allowing up to 2 missing cleavages. 
- 18 The mass tolerance for precursor ions is set as 20 ppm in first search and 5 ppm in main search, and the mass tolerance for fragment ions is set as 0.02 Da. 
- 19 Carbamidomethyl on Cys is specified as fixed modification, and acetylation on protein N-terminal and oxidation on Met are specified as variable modifications. FDR is adjusted to < 1%. 

Interacted proteins screening

- 20 Cells were lysed with Triton X-100 lysis buffer containing EDTA-free protease inhibitors (Roche).
- 21 Cell lysates are centrifuged at 12,000 rpm at 4°C for 10 min.
- 22 The indicated antibody is added to supernatants for immunoprecipitation based on Protein A/G Magnetic Beads.
- 23 Immunoprecipitates are denatured by boiling in SDS sample buffer and resolved by SDS-PAGE.
- 24 Protein bands in the whole lane are visualized via Coomassie blue staining (Beyotime).

- 25 Protein bands are subjected to in-gel digestion in 50 mM ammonium bicarbonate buffer overnight at 37°C with modified sequencing-grade trypsin (0.01 µg/µl).
- 26 The digested protein samples are separated by liquid phase chromatography and ionized by a nanoESI source, then passed to a tandem mass spectrometer Q-Exactive HF X (Thermo Fisher) for DDA (Data Dependent Acquisition) mode detection. 
- 27 Proteins are identified by searching the fragment spectra against the UniProt protein database (EMBL-EBI) using the Mascot search engine (v.2.3; Matrix Science) with the Proteome Discoverer software program (v.1.4; Thermo Fisher). 

LC-MS metabolomics

- 28 Samples are thawed at 4°C and mixed with cold methanol/acetonitrile (1:1, v/v) to remove protein.
- 29 The mixture is centrifuged at 14,000g at 4°C for 15 m and dried in a vacuum centrifuge.
- 30 The samples are re-dissolved in 100 µl acetonitrile/water (1:1, v/v) solvent.
- 31 To monitor the stability and repeatability of instrument analysis, quality control (QC) samples are prepared by pooling 10 µl of each sample and analyzed together with the other samples.
- 32 The LC-MS are performed using a UHPLC system (1290 series, Agilent Technologies) coupled to a triple quadrupole mass spectrometer (5500 QTRAP, AB SCIEX) in the multiple reaction monitoring (MRM) mode. 
- 33 Energy metabolites are monitored in electrospray negative ionization multiple reaction monitoring (MRM) mode only.
- 34 Amino column (particle size, 1.7 µm, 100 mm length × 2.1 mm ID) is used at a temperature of 45°C and a flow rate of 300 µl/min with a 4 µl injected sample volume.
- 35 The mobile phase A is 15 mM ammonium acetate in 100% water, and solvent B was 100% acetonitrile.
- 36 The linear gradient for energy metabolites is set as follows: 0-18 min: 0% B to 40% B, 18-18.1 min: 40% B to 90% B, 18.1-23 min: 90% B.

- 37 ESI-MS/MS conditions are set as follows: sheath gas temperature: 350°C; ion source gas 1: 45 psi; ion source gas 2: 45 psi; curtain gas: 30 psi; ionSapary Voltage Floating: 4500 V; nebulizer pressure: 40 psi.
- 38 The mass spectrometer is operated with a dwell time of 200 ms.
- 39 To construct the metabolite MRMLibrary, each metabolite standard (100 mg/ml) is first analyzed in ESI negative mode via flow injection to get the optimal MRM transition parameters. 
- 40 The retention time of each metabolite is determined by measuring the corresponding MRM(Q1/Q3) transition individually on the column. 
- 41 For **absolute quantification** of metabolite, a standard mixture sample containing all metabolites are measured together with biological samples in each experiment, and is used as the reference sample for peak alignment and metabolite identification. 
- 42 LC-MS quantification of energy metabolites are achieved with six-point standard curves using SUCCINIC ACID (D6, 98%, DLM-831-5, CIL) diluted in a relevant matrix matched to the analytical sample. 
- 43 Data acquisition and processing are accomplished using Multiquant software. 

GC-MS for isotype-labeled metabolomics (amino acids)

- 44 Cells cultured with isotype-labeled amino acids are grinded in liquid nitrogen and re-suspended in cold 50% aqueous methanol and inserted in dry ice for 30 m
- 45 Chloroform is added and vortexed for 30 second before centrifugation for 15 m at 14,000 rpm (4 °C), supernatant is transferred to new tubes to evaporate.
- 46 O-Isobutylhydroxylamine hydrochloride was added to the dried pellet and incubated for 20 m at 85 °C.
- 47 After cooling, *N-tert*-butyldimethylsilyl-*N*-methyltrifluoroacetamide (MTBSTFA) is added and samples are re-incubated for 60 m at 85 °C before centrifugation for 15 m at 14,000 rpm (4 °C).
- 48 The supernatant is used for GC/MS analysis. A Shimadzu QP-2020 GC-MS is programmed with an injection temperature of 250°C injection. 



- 49 The GC flow rate with helium carrier gas is 0.92 ml/min.
- 50 The GC column used is a 30 m x 0.25 mm x 0.25 mm DB-5ms.
- 51 GC-MS interface temperature is 300°C and (electron impact) ion source temperature is set at 200 °C, with 70 V ionization voltage.
- 52 The mass spectrometer is set to scan m/z range 50-700, with 1 kV detector.
- 53 To determine ^{13}C labeling, the mass distribution for known fragments of metabolites was extracted from the appropriate chromatographic peak.
- 54 For each fragment, the retrieved data comprise mass intensities for the lightest isotopomer (without any heavy isotopes, M0), and isotopomers with increasing unit mass relative to M0.
- 55 The ratio of universal labeled [^{13}C] amino acids relative to total level of AA is analyzed and Z-score normalized.

