


Aug 11, 2025

🌐 Masson's Trichrome staining for histology

DOI

<https://dx.doi.org/10.17504/protocols.io.e6nvwbwwwvmk/v1>



Samson Dowland¹

¹Charles Perkins Centre Histology Facility, The University of Sydney, NSW, 2006, Australia

The University of Sydney



Samson Dowland

The University of Sydney

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.e6nvwbwwwvmk/v1>

Protocol Citation: Samson Dowland 2025. Masson's Trichrome staining for histology. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.e6nvwbwwwvmk/v1>



License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: March 05, 2025

Last Modified: August 11, 2025

Protocol Integer ID: 123814

Keywords: histology, trichrome, histochemistry, staining, microscopy, pathology, histopathology, histological staining method, haematoxylin as the nuclear stain, nuclear stain, plasma stain, staining intensity, histology masson, contrast between tissue component, original red dyes for the plasma stain, fibre stain, other tissue component, tissue component, collagen yellow, collagen from muscle, collagen, embedded tissue, iron haematoxylin, green fibre stain, similar to the original masson, masson, contrast

Abstract

Masson's trichrome is a histological staining method commonly used to distinguish collagen from muscle, cytoplasm and other tissue components.

The method presented here remains close to the original used by Masson in 1929.

While there are many variants of the original method, there is commonality in that they all utilise three components - a nuclear stain (usually an iron haematoxylin), a fibre stain (to specifically dye the collagen) and a plasma stain (which acts as a counterstain). They also consistently use a polyacid to differentiate the plasma staining.

In this method it is noted that, if using routine formalin-fixed paraffin-embedded tissues, including the optional mordanting step will increase staining intensity and enhance contrast between tissue components. The same effect can also be achieved by using Bouin's fixative.

Nuclear stain: This method uses Weigert's haematoxylin as the nuclear stain by staining the nuclei dark purple-black. Some variants use Verhoeff's haematoxylin to additionally highlight elastin by staining it black.

Plasma stain: This method uses the original red dyes for the plasma stain, specifically a mixture of Acid fuchsin and Ponceau de xylidene. The latter has a number of alternate names, including Xylidine Ponceau, Xylidine Ponceau 2R, Ponceau 2R, Ponceau G, Ponceau Red, Acid Red 26. All have the same C.I. 16150.

Common variants use a mixture of Acid fuchsin and Biebrich scarlet, which provides a similar result, with a slightly darker red.

Heteropolyacid: Like the original, this method uses a 1% phosphomolybdic acid (PMA) solution for the heteropolyacid. Most variants use the same, but some recommend combining with phosphotungstic acid (PTA) in a 2.5% PMA + 2.5% PTA solution.

Fibre stain: The fibre stain used here is Methyl blue which stains collagen a vibrant, royal blue. This is similar to the original Masson's method, which was thought to use Aniline blue water soluble (a mixture of Water blue and Methyl blue). Some other variants use green fibre stains, such as Light green or Fast Green FCF. The Fast Green FCF resists fading better than Light green, which fades over time, causing sections to lose their contrast. One uncommon variant uses Metanil yellow to stain the collagen yellow.



Protocol materials

- ⊗ Hydrochloric Acid **Banksia Scientific Company (AJAX FineChem) Catalog #AJA256**
- ⊗ Ethanol **POCD Scientific Catalog #ETHABS70WV5**
- ⊗ Ethanol **Hurst Scientific Catalog #EIMS95-5L**
- ⊗ Xylene **POCD Scientific Catalog #XYL5P**
- ⊗ DPX Mountant for histology **Merck MilliporeSigma (Sigma-Aldrich) Catalog #06522**
- ⊗ Ethanol **POCD Scientific Catalog #ETHABS20M**
- ⊗ Haematoxylin **Fronine Laboratory Supplies Catalog #JL555S**
- ⊗ Iron (III) chloride hexahydrate **Merck MilliporeSigma (Sigma-Aldrich) Catalog #F2877**
- ⊗ Acid Fuchsin **Fronine Laboratory Supplies Catalog #JL301**
- ⊗ Ponceau de Xylidene **E. Gurr Ltd**
- ⊗ Glacial acetic acid **Ajax Finechem (ThermoFisher Scientific) Catalog #AJA1-500ML**
- ⊗ Phosphomolybdic acid **Merck MilliporeSigma (Sigma-Aldrich) Catalog #221856**
- ⊗ Methyl Blue **Merck MilliporeSigma (Sigma-Aldrich) Catalog #M5528**
- ⊗ Picric acid **May & Baker LTD**



Deparaffinisation and rehydration

32m

1 Dewax and rehydrate paraffin sections

32m

Mordanting (optional)

1h 15m

2 Incubate in saturated picric acid solution at 57 °C

1h



Picric acid **May & Baker LTD** Saturated picric acid solution:

[M] 1.3 Mass Percent in distilled water

Safety information

Picric acid is an unstable explosive when dry, where it becomes highly sensitive to shock, heat and friction. Ensure it remains hydrated at all times.

Picric acid reacts with metals to form unstable picrate salts which are more explosive. Do not use metal containers, lids or racks.

Picric acid is toxic if swallowed, inhaled or absorbed through the skin.

Consult the MSDS and follow all recommended safety precautions.

3 Cool to Room temperature

10m

4 Wash until no longer yellow

5m



Nuclear staining

13m

5 Stain in Weigert's Haematoxylin (combine solution A & B immediately before use)

10m

Solution A: 1 g Haematoxylin in 100 ml of 95% Ethanol

Solution B: 1.6 g Ferric chloride hexahydrate (Iron III chloride; $\text{FeCl}_3 \cdot 6\text{H}_2\text{O}$) + 95 ml distilled water + 1 ml concentrated HCl

Haematoxylin **Fronine Laboratory Supplies Catalog #JL555S** CI: 75290

Ethanol **Hurst Scientific Catalog #EIMS95-5L**



Iron (III) chloride hexahydrate **Merck MilliporeSigma (Sigma-Aldrich) Catalog #F2877**



Hydrochloric Acid **Banksia Scientific Company (AJAX FineChem) Catalog #AJA256**

- 6 Wash with running tap water to remove excess stain

2m



- 7 Differentiate with acid alcohol (1-2 dips) to clear the slides

[1% Hydrochloric Acid (HCl) in 70% Ethanol]



Hydrochloric Acid **Banksia Scientific Company (AJAX FineChem) Catalog #AJA256**



Ethanol **POCD Scientific Catalog #ETHABS70WV5**

- 8 Wash in tap water

1m



- 9 Examine slides under microscope to assess nuclear staining. Nuclei should be dark purple-black and cytoplasm or collagenous tissue components should be grey. If nuclei are under-stained, return to step 5. If sections are over-stained, return to step 7.



Plasma stain

11m

- 10 Stain in Acid Fuchsin + Ponceau de xylidine solution (make fresh, lasts 2-3 days)

5m

Solution: 1 g Acid Fuchsin + 1 g Ponceau de xylidine + 99 ml distilled water + 1 ml Glacial acetic acid



Acid Fuchsin **Fronine Laboratory Supplies Catalog #JL301** CI: 42685



Ponceau de Xylidene **E. Gurr Ltd** CI: 16150



Glacial acetic acid **Ajax Finechem (ThermoFisher Scientific) Catalog #AJA1-500ML**

- 11 Wash in distilled water



Heteropolyacid (differentiation of plasma stain)

11m

- 12 Differentiate and mordant in 1% PMA solution 2 – 6 mins (check progress of differentiation under the microscope every 2 mins, until collagen is colourless)

6m



PMA Solution: 1 g Phosphomolybdic acid + 100 ml distilled water

Phosphomolybdic acid **Merck MilliporeSigma (Sigma-Aldrich) Catalog #221856**

13 Drain – do not rinse



Fibre stain

2m 30s

14 Stain in Methyl blue solution

2m

Solution: 2 g Methyl blue + 2.5 ml Glacial acetic acid + 100 ml distilled water (solution can be re-used, lasts many years)

Methyl Blue **Merck MilliporeSigma (Sigma-Aldrich) Catalog #M5528** CI: 42780

Glacial acetic acid **Ajax Finechem (ThermoFisher Scientific) Catalog #AJA1-500ML**

15 Wash in distilled water



16 Differentiate with 1% Glacial acetic acid

30s

Glacial acetic acid **Ajax Finechem (ThermoFisher Scientific) Catalog #AJA1-500ML**

Dehydration, clearing and mounting

11m 30s

17 95% ethanol

30s

Ethanol **Hurst Scientific Catalog #EIMS95-5L**

18 100% ethanol

30s

Ethanol **POCD Scientific Catalog #ETHABS20M**

Note: ensure whole slide and slide rack is fully immersed in ethanol to remove all traces of water, preventing water contamination of the xylene



19 100% ethanol

30s

 Ethanol **POCD Scientific Catalog #ETHABS20M**

Note: ensure whole slide and slide rack is fully immersed in ethanol to remove all traces of water, preventing water contamination of the xylene

20 100% Xylene

5m

 Xylene **POCD Scientific Catalog #XYL5P**

21 100% Xylene

5m

 Xylene **POCD Scientific Catalog #XYL5P**

22 Coverslip with a resinous mounting medium, such as DPX

 DPX Mountant for histology **Merck MilliporeSigma (Sigma-Aldrich) Catalog #06522**

Results

23

Expected result

Nuclei: Purple-black
Cytoplasm: Red
Keratin: Red
Erythrocytes: Red
Muscle: Red
Collagen: Blue

