

Aug 23, 2018

Mate-pair large libraries preparation for assembly of the *Lateolabrax maculatus* genome

 [GigaScience](#)

 In 1 collection

DOI

<https://dx.doi.org/10.17504/protocols.io.ss2eege>

GigaScience Press

BGI



Chang Li

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.ss2eege>

External link: <https://doi.org/10.1093/gigascience/giy114>

Protocol Citation: Chang Li: Mate-pair large libraries preparation for assembly of the *Lateolabrax maculatus* genome. protocols.io <https://dx.doi.org/10.17504/protocols.io.ss2eege>

Manuscript citation:

hangwei Shao, Chang Li, Na Wang, Yating Qin, Wenteng Xu, Qun Liu, Qian Zhou, Yong Zhao, Xihong Li, Shanshan Liu, Xiaowu Chen, Shahid Mahboob, Xin Liu, Songlin Chen, Chromosome-level genome assembly of the spotted sea bass, *Lateolabrax maculatus*, *GigaScience*, Volume 7, Issue 11, November 2018, giy114, <https://doi.org/10.1093/gigascience/giy114>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: August 22, 2018

Last Modified: August 24, 2018

Protocol Integer ID: 14906

Keywords: lateolabrax maculatus genome this protocol, lateolabrax maculatus genome, lateolabrax maculatus, pair large libraries preparation for assembly, genome, process of the mate, large libraries preparation, pair

Abstract

This protocol is used to clarify the process of the mate-pair large libraries preparation for the L. maculatus.

Materials

STEP MATERIALS

 T4 DNA polymerase **Enzymatics**

 Biotin dNTP Mix **Invitrogen - Thermo Fisher**

 T3 DNA ligase **Enzymatics**

 M280 streptavidin beads **Invitrogen - Thermo Fisher**

 T4 DNA polymerase **Enzymatics**

 Biotin dNTP Mix **Invitrogen - Thermo Fisher**

 T3 DNA ligase **Enzymatics**

 M280 streptavidin beads **Invitrogen - Thermo Fisher**

Protocol materials

 M280 streptavidin beads **Invitrogen - Thermo Fisher**

 T4 DNA polymerase **Enzymatics**

 Biotin dNTP Mix **Invitrogen - Thermo Fisher**

 T3 DNA ligase **Enzymatics**

Genomic DNA interruption

- 1 The genomic DNA was fragmented using a Covaris E220 ultrasonicator (Covaris, Brighton, UK) to obtain ~2kb (for 2 kb library) and a Hydroshear (GeneMachines, CA, USA) to obtain ~5 kb, ~10 kb and ~20kb fragments (for 5 kb, 10 kb and 20 kb libraries respectively).

Equipment

new equipment	NAME
Covaris	BRAND
1	SKU
E220 ultrasonicator	SPECIFICATIONS

Equipment

new equipment	NAME
GeneMachines	BRAND
1	SKU
Hydroshear	SPECIFICATIONS

End-repair

- 2 Repaired using T4 DNA polymerase, (ENZYMATICS, Beverly, the U.S.) 30 min at 20 °C.

 T4 DNA polymerase **Enzymatics**

 00:30:00

 20 °C



Biotin Label

- 3 Add Biotin Label by Biotin dNTP Mix (5mM) 30 min at 20 °C.

 Biotin dNTP Mix **Invitrogen - Thermo Fisher**

 00:30:00

 20 °C

Fragment selection

- 4 These fragments were further selected into size ranges of 2–2.4 kb, 5–5.5 kb, 10–11 kb or 20–23 kb by agarose gel electrophoresis.

Fragment cyclizing

- 5 The T4 and T3 DNA ligase were used to connect the ring. And then, Covaris LE220 was used to cyclizing DNA fragments.

 T3 DNA ligase **Enzymatics**

End-repair

- 6 Fragmented DNA labeled with biotin was captured on M280 streptavidin beads (Invitrogen, CA, USA), followed by end repair (30 min. at 20°C, 1000 rotation per minute, rpm, vibrate for 15 sec. per 2 min.), A-tailing (30 min. at 37°C, 1000 rpm vibrate for 15 sec. per 2 min.).

 M280 streptavidin beads **Invitrogen - Thermo Fisher**

 00:30:00

 20 °C

 00:00:15

 00:30:00

 37 °C

 00:00:15

Add adapter

- 7 Adaptor ligation (1h at 20°C, 1000rpm vibrate for 15 sec per 2 min.).

 01:00:00

 20 °C

 00:00:15

PCR amplification

- 8 PCR amplifications on beads 95°C 3 min., (98 °C 20 sec., 60 °C 15 sec., 72 °C 45 sec.) for N cycles, 72 °C 10 min., 4°C hold (For 2 kb library, N=16; For 5 kb library, 10 kb library and 20 kb library, N=18)] using Enzymatics (MA, USA) and NEB (MA, USA) reagent.

 00:03:00 95 °C