



# MBP-Clu-tail purification from Escherichia coli cells

DOI

[dx.doi.org/10.17504/protocols.io.6qpvr35xzvmk/v1](https://dx.doi.org/10.17504/protocols.io.6qpvr35xzvmk/v1)

Andreas Bracher<sup>1</sup>, F Ulrich Hartl<sup>1</sup>

<sup>1</sup>Department of Cellular Biochemistry, Max Planck Institute of Biochemistry



Patricia Yuste-Checa

MPI Biochemistry

## Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account



DOI: <https://dx.doi.org/10.17504/protocols.io.6qpvr35xzvmk/v1>

**Protocol Citation:** Andreas Bracher, F Ulrich Hartl . MBP-Clu-tail purification from Escherichia coli cells. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.6qpvr35xzvmk/v1>

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Created:** January 26, 2024

**Last Modified:** June 01, 2024

**Protocol Integer ID:** 94482

**Keywords:** ASAPCRN, tail purification from escherichia coli cell, fusion protein maltose binding protein, maltose binding protein, tail purification, fusion protein, purification, escherichia coli cell, protein, mbp, clu

## Funders Acknowledgements:

Aligning Science Across Parkinson's

Grant ID: ASAP-000282

## Abstract

This protocol details how to efficiently purify the fusion protein Maltose binding protein (MBP)-Clu-tail (204-238) from *Escherichia coli*.

## Attachments



MBP-Clu(204-238)

pur...

1.3MB

## Materials

### Buffers

#### ▪ Binding buffer:













	A	B
	Tris/HCl pH 7.4	20 mM
	NaCl	200 mM
	EDTA	1 mM

- Elution buffer: Binding buffer +  10 millimolar (mM) maltose (final concentration)






## Troubleshooting




## His<sub>6</sub>-Ubiquitin-GFP-Clu-tail expression and cell lysis



- 1 Express MBP-Clu-tail in E. coli BI21 (DE3) codon+RIL cells cultured in  1 L LB Medium containing  2  $\mu$ L glucose at  37 °C with  1 millimolar (mM) IPTG during  02:00:00 . 2h
- 2 Centrifuge culture and keep pellet. 
- 3 Re-suspend pellet in  25 mL volume of ice-chilled binding buffer, add Complete protease inhibitor cocktail (Roche) and  1 millimolar (mM) phenylmethylsulfonyl fluoride (PMSF).
- 4 Lyse cells by ultrasonication in ice bath (10 cycles of  00:00:30 ultrasonication with  00:01:30 intermittent cooling). 2m
- 5 Clear lysate by centrifugation at  22000 rpm in a JA25.50 rotor at  4 °C .

## Amylose affinity chromatography



- 6 Load supernatant onto a  20 mL Amylose Resin (New England Biolabs) column previously equilibrated with binding buffer by gravity flow at  4 °C .
- 7 Wash the column with 12 CV of ice-chilled binding buffer. 
- 8 Elute MBP-Clu-tail protein with 12× 3 mL of ice-chilled Elution buffer. Collect fractions of  3 mL volume. Store fractions  On ice .
- 9 Analyze eluted fractions by SDS-PAGE and Coomassie blue staining.



10 Pool fractions containing MBP-Clu-tail.

11 Concentrate pool to less than  10 mL volume by ultrafiltration using 10 kDa cut-off spin concentrator at  4 °C .

## Size exclusion chromatography

12 Apply concentrate on a HiLoad 26/600 Superdex-200 (Cytiva 28-9893-36) column equilibrated with PBS. Develop the column at  4 °C and collect  10 mL fractions.


13 Analyze eluted fractions by SDS-PAGE and Coomassie blue staining.



- 14 Merge fractions with MBP-Clu-tail peak. Concentrate to 1.5 mL volume by ultrafiltration using 10 kDa cut-off spin concentrator at  $4^{\circ}\text{C}$  , aliquot and flash-freeze purified MBP-Clu-tail in liquid nitrogen for storage at  $-70^{\circ}\text{C}$  .

#### Note

MBP-Clu-tail appears as a double band. In contrast to the lower band, the upper band is sensitive to cleavage by furin, suggesting that a protease from *E. coli* partially cleaves close to the furin site in MBP-Clu-tail.



#### Note

Concentrations were determined by absorbance at 280 nm using absorbance coefficients of  $66,350 \text{ M}^{-1} \text{ cm}^{-1}$  or  $1.645 \text{ L g}^{-1} \text{ cm}^{-1}$  for MBP-Clu-tail.

**Approximate yield:** From 1 L of culture around 15 mg of pure MBP-Clu-tail were obtained.