

Nov 26, 2024

Measurement of *Auxenochlorella protothecoides* (UTEX 250) Chlorophyll Fluorescence Kinetics (OJIP) and Photosynthetic Efficiency (Fv/Fm) using the Aqua Pen (AP110/C).

DOI

<https://dx.doi.org/10.17504/protocols.io.yxmvme595g3p/v1>

Dimitrios Camacho¹, Sabeeha S. Merchant¹

¹Department of Molecular and Cell Biology, University of California, Berkeley, California 94720, USA

Merchant Lab UC Berkel...



Dimitrios J. Camacho

Department of Molecular and Cell Biology, University of Cali...

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.yxmvme595g3p/v1>

Protocol Citation: Dimitrios Camacho, Sabeeha S. Merchant 2024. Measurement of *Auxenochlorella protothecoides* (UTEX 250) Chlorophyll Fluorescence Kinetics (OJIP) and Photosynthetic Efficiency (Fv/Fm) using the Aqua Pen (AP110/C)..

protocols.io <https://dx.doi.org/10.17504/protocols.io.yxmvme595g3p/v1>



License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: July 22, 2024

Last Modified: November 26, 2024

Protocol Integer ID: 103884

Keywords: chlorophyll fluorescence kinetic, chlorophyll fluorescence induction kinetic, measurement of auxenochlorella protothecoide, photosynthetic efficiency, auxenochlorella protothecoide, optimization of the super pulse intensity, super pulse intensity

Funders Acknowledgements:

National Institutes of Health (NIH): Molecular Basis of Cell Function T32 Training Grant

Grant ID: 5T32GM007232-44

US Department of Energy (DOE), Office of Biological and Environmental Research (BER): Systems Engineering of *Auxenochlorella protothecoides*: from Photosynthesis to Biofuels and Bioproducts

Grant ID: DE-SC0023027

University of California, Berkeley, Chancellor's Fellowship

Grant ID: N/A

Disclaimer

This protocol was adapted from the Photon Systems Instruments AquaPen-C AP110-C manual linked below:
https://handheld.psi.cz/documents/AquaPen_Manual-verze_02_2021.pdf

Abstract

This protocol describes a method for examining the chlorophyll fluorescence induction kinetics, "Kautsky effect," of *Auxenochlorella protothecoides* (UTEX 250) cells. Photosynthetic efficiency (Fv/Fm) measurements are also acquired during the OJIP protocol. This protocol details the optimization of the super pulse intensity setting for a control sample and the application of those settings to probe experimental samples.

Guidelines

Please read the entirety of the protocol before attempting. A summary of the step is highlighted in bold at the beginning of each step, where applicable. Detailed descriptions are provided in the rest of the step.

Normalize the sample input by cell density, chlorophyll content, or optical density (O.D. 750) to ensure relatively equal and complete light penetration of the samples. We recommend using 10^6 cells / mL in 3 mL.



Materials

Aqua pen Model: AP110/C, Photon Systems Instruments (firmware: 1.2.1.6)

Aqua pen cuvettes and caps, cat. no. PSI16107, Photon Systems Instruments

Box for dark adaptation

Green safelight

Computer

Installed FluorPen Software: <https://psi.cz/support/downloads/hd001/ap001/>

RAININ P1000 and tips

RAININ P20 and tips

Hemocytometer

Microscope

Kimwipes

Safety warnings

❗ The cuvettes are not disposable and should not be used for any other purpose.

Do not use acetone or alcohol to clean cuvettes.

Do not wash cuvettes in a dishwasher.

Do not scratch cuvettes.

Wash cuvettes with soapy water and air dry.

Take extreme caution when connecting the charging/ data transfer USB cable to the aqua pen. The connector is prone to irreversible damage. For this reason, we recommend keeping the cable attached and coiled.

If possible, keep the AquaPen in a protective case.

Super pulse intensity optimization of control sample

- 1 Optimize the super pulse intensity based on the wild type, nutrient replete, or other untreated control samples of your experiment.

Wild-type *Auxenochlorella* (UTEX 250) cells cultured in replete TAP medium typically have an optimal super pulse setting of 70% during the exponential phase of growth with $100 \mu\text{mol photons} \cdot \text{m}^{-2} \cdot \text{s}^{-1}$ light. If you are using cells that have been cultured in other conditions, follow the steps in this section to achieve the optimal super pulse intensity for your sample type.

Normalize the sample input by cell density, chlorophyll content, or optical density (O.D. 750) to ensure relatively equal and complete light penetration of the samples. We recommend using 10^6 cells / mL in 3 mL.

- 2 Count the cell density of your control sample using a hemocytometer. See (Camacho & Merchant, 2024) protocol for detailed instructions on how to count *Auxenochlorella protothecoides* cells using a hemocytometer.
- 3 Dilute the control sample to 10^6 cells / mL using the medium in which the cells were cultured.
- 4 Pipette 3 mL of the diluted cell suspensions to 10 cuvettes. You will test different super pulse light intensities on these samples to determine the optimum super pulse light intensity, also known as "saturating pulse".

Note

Hold the cuvettes at the top 1 cm. Do not touch the cuvette in the light path area.

- 5 Place the cuvettes in a dark cuvettes' holder/ box.

Move the box of cuvettes to a dark room with a green safelight installed.

Turn the green safelight off and incubate the control samples in the dark for 15 min.

- 6 **Set the pulse color to blue (455 nm).**


Press and hold the "SET" button for 3 sec to turn the aqua pen on.



Make sure that the date and time are correct and that the battery is sufficiently charged.
Press the **"MENU"** button to toggle through the options until the cursor points to **"Settings"**



Press the **"SET"** button to enter into the settings menu
Press **"MENU"** until the **"Measur. color"** option is selected




Press **"SET"** to enter the **"Measur. color"** setting.
Press **"SET"** until **455 nm** (Blue) is selected.



7 Change the super pulse intensity.

Press the **"MENU"** button to exit the **"Measur. color"** setting and back to the settings page.

Press **"MENU"** until "Super pulse" is selected. Press the **"SET"** button to enter into the **"Super pulse"** settings.



Set the super pulse intensity to 10%

In the later steps of this protocol, you will increase the super pulse intensity by increments of 10%.

The percentage indicated is the percentage of the maximum super pulse intensity of $3000 \mu\text{mol photons} \cdot \text{m}^{-2} \cdot \text{s}^{-1}$ blue light.


	A	B
	%	$\mu\text{mol photons} / \text{m}^2 / \text{s}$
	10	300
	20	600
	30	900
	40	1200

	A	B
	50	1500
	60	1800
	70	2100
	80	2400
	90	2700
	100	3000

8 Obtain an OJIP measurement of one control sample at a super pulse setting of 10%

Make sure that you are in a dark room. Turn the green safelight on and minimize the exposure of the sample to the green light.


Press MENU until you reach the **<<Return>>** option to exit the settings page.



Press "**MENU**" until you reach the "**Measurement**" option.

Press "**SET**" to access the measurement options.

Press "**MENU**" until you reach the "**OJIP**" option.



Mix the cells by pipetting up and down 5 times or by inverting the capped cuvette 5 times. Cells will settle at the bottom of the cuvette and will yield inaccurate results if not resuspended properly right before the measurement.

Wipe the outside of the cuvette with a Kimwipe.

Make sure that the cuvette is dry and clean before inserting it into the Aqua Pen.

Press the "**SET**" button to acquire the OJIP measurement. The results will not be displayed. You will see a percentage indicating the progress of the measurement. The measurement will take a few seconds. Record the counter number and the time that the measurement was taken.

9 Obtain an OJIP measurement of a new control sample at a super pulse setting of 20%. Repeat and increase by increments of 10% until all 10 samples have been measured.

After the measurement of the first sample, remove the cuvette from the AquaPen and set it aside.

Note

Make sure to remove the sample from the Aqua Pen before changing the settings. Holding the AquaPen at a weird angle, while a cuvette is still inside, may cause the sample to leak into the AquaPen's internal electrical components.

Press "MENU" to until the **<<Return>>** option is selected.

Repeat step 7 to change the super pulse intensity to 20%.

Use a new dark adapted control sample for the OJIP measurement.
Do not reuse samples that have just been measured.

Repeat step 8 to obtain the new OJIP measurement.

Repeat steps 7-9, increasing the superpose intensity by 10% until all 10 cuvettes from step 4 have been measured.

10 Export the data using the FluorPen Software.

Connect the USB cable to the AquaPen and your computer.

Make sure that the FluorPen software is installed. See materials for a link to download the software.

You may need to download and install a driver (also included in the link).

The FluoroPen software will only work on windows computers.

Open the FluorPen program and connect to the device.

If it is your first time using this software on your computer, you will need to register the AquaPen's serial number into the program.

Click "**Help**" in the toolbar then select "**Register**".

Enter the serial number into the dialogue box.

Each time you open up the software, the aqua pen will not be recognized and you will have to set up the connection.

Click on **Setup** in the tool bar and click on "Device ID" (see image below).




Image of FluorPen software and Setup tab opened.

The AquaPen will connect and the download icon should become visible.

Click on the Download icon.

A table with all data will automatically display. See image below. Each column is an independent measurement.

Index	18	19	20	21	22	23	24	25	26	27
Time	14:08:29 25.11.2024	14:09:24 25.11.2024	14:10:26 25.11.2024	14:11:10 25.11.2024	14:11:58 25.11.2024	14:12:39 25.11.2024	14:13:26 25.11.2024	14:14:19 25.11.2024	14:15:04 25.11.2024	14:15:47 25.11.2024
	OJIP	OJIP	OJIP	OJIP	OJIP	OJIP	OJIP	OJIP	OJIP	OJIP
Bdg	181	209	209	260	260	292	260	292	292	292
Fo	1380	3330	4776	6762	8778	10956	14174	16613	19312	23831
Fj	2369	5523	8534	11866	15085	18856	23082	27309	31081	36412
Fj	2212	6770	10159	13862	17523	21717	26366	31048	35372	41452
Fm	3754	7527	10957	14857	18531	22758	27374	32121	36542	42655
Fv	2174	4197	6181	8095	9753	11802	13200	15508	17720	18824
Vj	0.363	0.523	0.608	0.631	0.647	0.669	0.675	0.695	0.683	0.668
Vi	0.751	0.820	0.871	0.880	0.897	0.912	0.924	0.931	0.932	0.936
Fm/Fo	2.376	2.260	2.294	2.157	2.111	2.077	1.931	1.933	1.892	1.790
Fv/Fo	1.376	1.260	1.294	1.197	1.111	1.077	0.931	0.933	0.892	0.790
Fv/Fm	0.579	0.558	0.564	0.545	0.526	0.519	0.482	0.483	0.472	0.441
Mo	0.318	0.675	1.122	1.414	1.680	1.961	2.118	2.340	2.377	2.418
Area	3167640	5953893	4925601	5963410	7861858	7632374	10318201	5710770	6142038	5483154
Fix Area	3611768	7260667	10803220	14661337	18355594	22510302	27095552	31724098	35996776	41926088
H&Ch Area	2031832	4940000	6027410	7899608	9577946	11554741	12936119	15111762	16685549	18906492
Sm	1457.148	1418.607	796.894	736.678	806.096	646.702	781.682	368.247	356.473	291.285
N	1278.006	1831.963	1470.797	1651.753	2094.036	1894.918	2453.505	1249.202	1240.766	1053.823
Phi_Po	0.579	0.538	0.564	0.545	0.526	0.519	0.482	0.483	0.472	0.441
Phi_Po	0.637	0.477	0.392	0.369	0.353	0.331	0.325	0.310	0.317	0.332
Phi_Eo	0.369	0.266	0.221	0.201	0.196	0.171	0.157	0.150	0.149	0.146
Phi_Do	0.421	0.442	0.436	0.455	0.474	0.481	0.518	0.517	0.528	0.559
Phi_Pav	957.166	974.599	973.827	975.138	982.933	981.705	987.719	978.903	979.377	976.642
Phi_Abs	1.595	0.497	0.255	0.170	0.123	0.094	0.069	0.060	0.056	0.048
ABS/RC	1.514	2.316	3.272	4.115	4.936	5.650	6.509	7.026	7.382	8.198
TRo/RC	0.877	1.291	1.846	2.242	2.598	2.930	3.139	3.392	3.481	3.618
Description										

Table with data is automatically displayed after data is downloaded.

The index corresponds to the counter value from step 8. Verify that the index and time of the measurement matches.



FluorPen

File Device Setup Help

Untitled - 1

Index	18
Time	14:08:29 25.11.2024
	OJIP
	Bckg 181
	Fo 1580
	Fj 2369
	Fi 3212
	Fm 3754
	Fv 2174
	Vj 0.363
	Vi 0.751
	Fm/Fo 2.376
	Fv/Fo 1.376
	Fv/Fm 0.579
	Mo 0.318
	Area 3167840
Value	Fix Area 3611768
	HACH Area 2031832
	Sm 1457.148
	Ss 1.140
	N 1278.006
	Phi_Po 0.579
	Psi_o 0.637
	Phi_Eo 0.369
	Phi_Do 0.421
	Phi_Pav 957.166
	Pi_Abs 1.595
	ABS/RC 1.514
	TRo/RC 0.877
	...
Description	

Data from a single measurement. Fv/Fm value is highlighted in a red box.

11 Determine the optimal super pulse intensity.

Record the Fv/Fm values for each cuvette and their corresponding super pulse light intensities.



Plot the Fv/Fm values (y axis) against the super pulse intensity (x axis).

The super pulse intensity with the highest Fv/Fm value is your optimal super pulse intensity.

If the highest Fv/Fm value is measured at 100% super pulse intensity, it is possible that the sample is too dense and will need to be diluted.


Measurement of OJIP and Fv/Fm from experimental samples.

- 12 Change the super pulse intensity to the optimal intensity (See step 7).
- 13 Determine the cell density of your samples (see steps 1 and 2).

Dilute each sample to 10^6 cells / mL in 3 mL of medium and transfer to a cuvette (same as step 3).
- 14 Dark adapt the samples for 15 min (see step 5).
- 15 Obtain the OJIP measurement using the optimized super pulse intensity (similar to Step 8).
- 16 Export the data (see step 10).

Clean up

- 17 Turn the AquaPen off by pressing and holding the **"MENU"** button.
You may also turn the instrument off by navigating back to the home page and selecting **"Turn off"**



- 18 Remove the final sample from the AquaPen.

Empty the cuvettes into an autoclavable liquid waste container and autoclave or kill the cells with 10% bleach.

- 19 Clean the cuvettes and caps by soaking in a soapy water solution.
Do not use bleach, acetone, or any alcohol solution to clean the cuvettes.
Do not wash cuvettes in a dish washer.

Air dry the Cuvettes and caps and store them safely in a box.

Protocol references

Dimitrios Camacho, Sabeeha Merchant 2024. Determination of *Auxenochlorella protothecoides* Cell Density With a Hemocytometer. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.rm7vzjk12lx1/v1>