

Sep 30, 2019 Version 2

Measurement of Left Ventricular Hemodynamic Parameters in Intact Mice V.2

DOI

dx.doi.org/10.17504/protocols.io.7tdhni6



E. Dale Abel¹

¹University of Iowa

Diabetic Complications Consortium Tech. support email: rmcindoe@augusta.edu



Lili Liang

OPEN ACCESS



DOI: dx.doi.org/10.17504/protocols.io.7tdhni6

External link: https://www.diacomp.org/shared/document.aspx?id=43&docType=Protocol

Protocol Citation: E. Dale Abel 2019. Measurement of Left Ventricular Hemodynamic Parameters in Intact Mice. **protocols.io** https://dx.doi.org/10.17504/protocols.io.7tdhni6

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: September 30, 2019

Last Modified: September 30, 2019

Protocol Integer ID: 28229

Keywords: Measurement of Left Ventricular Hemodynamic Parameters, cardiovascular



Abstract

Summary:

This protocol describes the procedure used by the DiaComp for measuring left ventricular hemodynamic parameters in intact mice.

Diabetic Complication:



Cardiovascular



1 Measurement of in vivo cardiac contractile performance using a micro-manometer tipped catheter:

Indices of left ventricular contractile performance are determined in mice under isoflurane anesthesia by percutaneously introducing a 1.4 Fr. Millar catheter into the carotid artery and gently advancing this into the left ventricular cavity. Pressure waves are recorded captured and analyzed with pClamp 8.2 software. This technique provides the following indices of left ventricular function.

- (1) Left ventricular systolic pressure (LVSP) Units: mmHg
- (2) Left ventricular end diastolic pressure (LVEDP) Units: mmHg
- (3) Heart Rate Units: beats per minute
- (4) Positive dP/dt First derivative of LV pressure. This index is a measure of the rate of pressure development in the left ventricle. Units: mmHg/msec.
- (5) Negative dP/dt First derivative of LV pressure. This index is a measure of the rate of ventricular relaxation, or the rate of pressure decay in the LV Units: mmHg/msec.