


Oct 18, 2019

## Measurement of XylE (Catechol 2,3-Dioxygenase) enzyme activity by microplate reader

DOI

[dx.doi.org/10.17504/protocols.io.7suhnew](https://dx.doi.org/10.17504/protocols.io.7suhnew)



Ji Gao<sup>1</sup>

<sup>1</sup>Tsinghua University



Ji Gao

Tsinghua University

### Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.7suhnew>

**Protocol Citation:** Ji Gao 2019. Measurement of XylE (Catechol 2,3-Dioxygenase) enzyme activity by microplate reader. protocols.io <https://dx.doi.org/10.17504/protocols.io.7suhnew>

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** September 30, 2019


**Last Modified:** October 18, 2019

Protocol Integer ID: 28212

**Keywords:** microplate reader, enzyme activity by microplate reader, dioxygenase, simple measurement of xyle, measurement of xyle, enzyme activity, dihydroxybenzene, product of the enzyme, enzyme, catechol, photometric mode of microplate reader, xyle

## Abstract

Simple measurement of XylE (Catechol 2,3-Dioxygenase) enzyme activity by microplate reader. Catechol 2,3-Dioxygenase can catalyze catechol (1,2-Dihydroxybenzene) to 2-HMS, which has a high absorbance at 377 nm, so we may use photometric mode of microplate reader to measure the concentration of 2-HMS, product of the enzyme.



## Guidelines

This is part of a tryout for [our iGEM project](#).

The microplate we used is listed here.

Equipment	
Varioskan LUX Multimode Microplate Reader	NAME
microplate reader	TYPE
Thermo Scientific	BRAND
VL0000D2	SKU
<a href="https://www.thermofisher.com/order/catalog/product/VL0000D2?SID=srch-srp-VL0000D2">https://www.thermofisher.com/order/catalog/product/VL0000D2?SID=srch-srp-VL0000D2</a> <sup>LINK</sup>	

Perform the measurement at 🌡️ Room temperature or 🌡️ 25 °C .

## Materials

### MATERIALS

⊗ double distilled water (ddH2O)

⊗ 12-Dihydroxybenzene **Adamas-beta Catalog #17253**

⊗ Corning™ 96-Well Clear Bottom Black or White Polystyrene Microplates **Fisher Scientific Catalog #07-200-565**

## Troubleshooting

## Safety warnings

- ⚠️ Catechol can be hazardous if not properly operated.  
Please refer to <https://pubchem.ncbi.nlm.nih.gov/compound/catechol#section=Safety-and-Hazards>



## Before start


Use 1M 100 millimolar (mM) catechol water solution as stock, and use 1M 10 millimolar (mM) as working solution.

We use LB cell culture (BL21 strain in our project) to perform measurement.

Add cell culture of a good state (logarithm phase) to 96-well plate.



## Manual steps

- 1 Add  100  $\mu\text{L}$  cell culture per well.

### Note

It is recommended to use replicates and controls to avoid mistakes or deviation.



- 2 Turn on the microplate reader and computer. Run the software to set the protocol and plate layout.

## Protocol for the instrument


- 3 Shake the plate for  00:00:05 at  600 rpm .

- 4 Measure the absorbance of cell culture at 600 nm as the estimation of cell amounts.

- 5 Measure the absorbance of cell culture at 377 nm as the baseline of A377 before the reaction.



- 6 Plate out and add  2.5  $\mu\text{L}$  catechol (  10 millimolar (mM) ) to every well, then plate in immediately.

### Note

The working concentration of substrate is  0.25 millimolar (mM) here. Other working concentration could also work.

- 7 Shake the plate for  00:00:05 at  600 rpm .



- 8 Measure the absorbance of cell culture at 377 nm, continuously (like  00:00:10 per reading) or after a certain period of time (like  00:01:30 ).

## Manual steps

- 9 Export the data from the software.

### Safety information


Properly dispose of the contaminated cell culture and microplates!

- 10 Normalize the A<sub>377</sub> value by optical density at 600 nm.

Make a Time-A<sub>377</sub> plot like this:

### Expected result

#### Normalized measurement of Xyle-FUS activity



Or simply subtract the absorbance before the reaction from the value after a specific period of time:

## Expected result

