

Feb 04, 2022

Measuring nitrate/nitrite (NO_x) concentrations in water samples

DOI

dx.doi.org/10.17504/protocols.io.b4nnqvde

Jacob Waldbauer¹, Amy Amy Zimmerman²

¹University of Chicago; ²Pacific Northwest National Lab



Jacob Waldbauer

OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.b4nnqvde

Protocol Citation: Jacob Waldbauer, Amy Amy Zimmerman 2022. Measuring nitrate/nitrite (NO_x) concentrations in water samples. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.b4nnqvde>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: February 04, 2022

Last Modified: February 04, 2022

Protocol Integer ID: 57774

Abstract

Colorimetric (i.e., UV-visible spectroscopy) assay for measurements of nitrate/nitrite down to 0.5 $\mu\text{mol}\cdot\text{L}^{-1}$ concentration. Vanadium(III) reduces nitrate to nitrite prior to measurement with traditional Griess reaction. Nitrate and nitrite concentrations can be measured for the same sample. NOTE: This protocol is written for measurement in 96-well plates.



Materials

Spectrophotometer (or plate reader)

200 and 1000 μ L filter tips

200 and 1000 μ L pipettes

microcentrifuge tubes

Tube racks

15 mL conical tubes (if necessary)

Vortexer

96-well microplate with lid, clear

saturated VCl_3

2% SULF

0.1% NEDD

- 1 Making standards (sodium nitrate and sodium nitrite)
 - 1.1 Prepare 200 μM stock solution: Dilute 1:500 from 0.1 M stock solution \Rightarrow 20 μL + 9.980 mL nanopure water.
 - 1.2 Dilute the stock solution to the following concentrations in nanopure water: 0, 10, 20, 30, 40, 60, 100, 150, 200 μM .
- 2 Making Reagents and Solutions.
 - 2.1 Saturated vanadium(III) chloride (VCl_3) solution: Dissolve 100 mg in 12.5 mL 1 M HCl. Remove excess precipitates with nylon-66 syringe filter. DO NOT HEAT. Store solution at 4°C protected from light. Use within 2 weeks. NOTE: Appears viable after 2 months. Development of a lighter blue color indicates oxidation— authors suggest discarding solution, but light blue color observed in our lab shortly after preparation, and worked fine.
 - 2.2 2% (w/v) Sulfanilamide (SULF) solution: Add 1 g to 50 mL HCl. Heat to dissolve (in water bath set to $\sim 55^\circ\text{C}$). Filter to remove trace particulates. Store solution at 4°C protected from light. Stable for several months. Discard if colored.
 - 2.3 0.1% (w/v) N-(1-Naphthyl)ethylenediamine dihydrochloride (NEDD) solution: Add 50 mg to 50 mL nanopure water. Heat to dissolve (in water bath set to $\sim 55^\circ\text{C}$). Filter to remove trace particulates. Store solution at 4°C protected from light. Stable for several months. Discard if colored.
- 3 Assay set-up.
 - 3.1 Calculate total volume of SULF/NEDD mixture needed to run triplicate wells for each standard and sample: ($\#$ standards + $\#$ samples + 1 extra) \times 3 \times 80 μL = total vol (μL). a. If quantifying both total NO_x and NO_2^- , double the $\#$ samples. b. Account for both sample (no added SULF/NEDD mixture) and reagent blanks (nanopure plus reagents, i.e., 0. μM NO_x).
 - 3.2 Prepare the total volume of SULF/NEDD mixture needed by mixing equal parts of the reagents (i.e., if need 3 mL total volume, add 1.5 mL SULF to 1.5 mL NEDD in 15 mL tube). Vortex to mix.
 - 3.3 Aliquot 80 μL of each sample/standard to the appropriate wells of a clear 96-well plate.

- 3.4 Add 80 μL of VCl_3 to each well for total NO_x concentrations and 80 μL of nanopure water to each well for nitrite concentrations. Mix briefly by pipetting (use multichannel pipet).
- 3.5 Quickly add 80 μL of SULF/NEDD mixture to each well. Mix briefly by pipetting.
- 3.6 Incubate at room temperature for 30 minutes (45 minutes maximum).
- 3.7 Measure absorbance on plate reader.
- 4 Reading plates.
 - 4.1 Turn on Tecan Infinite 200 PRO plate reader 20–30 minutes prior to use.
 - 4.2 Once warmed up, open the iControl software on MLCLab-PC.
 - 4.3 Open file “Miranda_NOx_96well” (or “Miranda_NOx_24well” if appropriate).
 - 4.4 Load the plate—check whether the “plate with cover” box is checked (since using clear plates for this, can be read with lid on).
 - 4.5 Read absorbance at 520 ± 9 nm (25 flashes). Program automatically opens an Excel file that documents read parameters and data.
- 5 Analyzing data.
 - 5.1 Subtract the absorbance values of the samples mixed with diluting solution/medium instead of the SULF/NEDD mixture (sample blanks) from the corresponding reacted sample absorbances (= corrected sample absorbance).
 - 5.2 Subtract the average absorbance of the nanopure water tubes (i.e., 0 μM NO_x) mixed with diluting solution/medium instead of the SULF/NEDD mixture (standard blanks) from



the absorbances of all the standards.

- 5.3 Plot corrected absorbance (y) vs. concentration (x) for all standards to establish a standard curve with linear regression.
- 5.4 Use the equation of the standard curve to calculate sample concentration from fluorescence.