



Feb 24, 2023

Measuring photophysiological traits of diatoms from Rapid Light Curves using a Water-PAM

DOI

dx.doi.org/10.17504/protocols.io.36wgg4zmovk5/v1



Phoebe Argyle^{1,2}, Jana Hinnert³, Nathan G. Walworth⁴, Sinead Collins⁵, Naomi M. Levine⁴, Martina A. Doblin^{1,6}

¹Climate Change Cluster, University of Technology Sydney, Ultimo NSW Australia;

²Ministry of Marine Resources, Cook Islands;

³Institute of Coastal Ocean Dynamics, Helmholtz-Zentrum Hereon, 21502, Geesthacht, Germany;

⁴Department of Biological Sciences, University of Southern California, Los Angeles, CA 90089-0371;

⁵Institute of Evolutionary Biology, University of Edinburgh, Edinburgh, UK;

⁶Sydney Institute of Marine Science, Mosman, NSW Sydney



Phoebe Argyle

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.36wgg4zmovk5/v1>

Protocol Citation: Phoebe Argyle, Jana Hinnert, Nathan G. Walworth, Sinead Collins, Naomi M. Levine, Martina A. Doblin 2023. Measuring photophysiological traits of diatoms from Rapid Light Curves using a Water-PAM . **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.36wgq4zmovk5/v1>

Manuscript citation:

Argyle, P.A., Hinnert, J., Walworth, N.G., Collins, S., Levine, N.M., Doblin, M.A., 2021. A High-Throughput Assay for Quantifying Phenotypic Traits of Microalgae. *Frontiers in Microbiology* 12(2910).

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: August 24, 2021

Last Modified: February 24, 2023

Protocol Integer ID: 52690

Keywords: diatom, microalgae, photophysiology, fluorometry, traits, measuring photophysiological trait, photophysiological traits of diatom, throughput assay for quantifying phenotypic trait, photosynthetic activity, ik for microalgae culture, throughput assay, quantifying phenotypic trait, diatom plasticity, phenotypic traits of microalgae, rapid light curve protocol, rapid light curve, phenotypic trait, multivariate trait analysis, biological axis, microalgae culture, set of biological axis

Funders Acknowledgements:

Moore Foundation Marine Microbes Initiative

Grant ID: MMI 7397

Abstract

This is a brief protocol for how to measure ETR_{max}, alpha (α) and I_k for microalgae cultures using the Rapid Light Curve protocol included in the WinControl software used with a Walz Water-PAM instrument.

Based on the protocol developed in:

Ralph, P.J., Gademann, R., 2005. Rapid light curves: a powerful tool to assess photosynthetic activity. *Aquat. Bot.* 82(3), 222-237.

This method was used in:

Argyle, P. A., Walworth, N. G., Hinnert, J., Collins, S., Levine, N. M., & Doblin, M. A. (2021). Multivariate trait analysis reveals diatom plasticity constrained to a reduced set of biological axes. *ISME Communications*, 1(1), 59.

Argyle, P. A., Hinnert, J., Walworth, N. G., Collins, S., Levine, N. M., & Doblin, M. A. (2021). A high-throughput assay for quantifying phenotypic traits of microalgae. *Frontiers in microbiology*, 12, 706235.

Image Attribution


Michaela Larsson, UTS Sydney



Materials

Equipment	
Water-PAM	NAME
PAM fluorometer	TYPE
WALZ	BRAND
WATER-PAM-II	SKU
https://www.walz.com/products/chl_p700/water-pam/introduction.html	LINK

Troubleshooting

Prepare culture

- 1 Take a  2 mL aliquot of microalgae culture into a glass cuvette compatible with the WATER-PAM instrument.

When the overall culture volume is limited (i.e. <2 mL of culture available for the measurement), take  0.5-1 mL and dilute with sterile seawater to a final volume of  2 mL .

Notes:

The optimal volume for the WATER-PAM is  2-3 mL .

Make measurements

- 2 Set up the WATER-PAM according to the manufacturer's instructions.

Equipment

Water-PAM

NAME

PAM fluorometer

TYPE

WALZ



BRAND

WATER-PAM-II

SKU

https://www.walz.com/products/chl_p700/water-pam/introduction.html^{LINK}

- 3 Initiate a Rapid Light Curve by clicking Program > Light Curve in WinControl.




To achieve a successful RLC, there should be at least 3 points on the rise, multiple points to estimate the ETRmax, and at least one point on the decline at the highest irradiance (see instrument manual).

If the curve does not achieve these criteria, this may require optimization of the instrument settings such as the gain and irradiance steps.

Modify the Light Curve Intensity (LI) to change the first irradiance step of the curve, keeping in mind there will always be 8 steps to the curve.

NB: Do not adjust the gain settings in between samples as this renders the data incomparable.

The curve width is usually set to 10 seconds.



NB: If comparing different cultures, it is important to conduct all measurements within a similar time window based on the photoperiod. I.e. within a 1-2 hour window 2 hours after the onset of light.