

Jan 13, 2020

# Measuring relative reactivity of mouse TCRs against a mouse cancer cell line

 In 1 collection

DOI

[dx.doi.org/10.17504/protocols.io.ba8gihtw](https://dx.doi.org/10.17504/protocols.io.ba8gihtw)

Bulent Arman Aksoy<sup>1</sup>, Pinar Aksoy<sup>1</sup>, Elinor Gottschalk<sup>1</sup>, Jeff Hammerbacher<sup>1</sup>

<sup>1</sup>Medical University of South Carolina

Hammer Lab

Tech. support phone: +18437924527 email: [arman@hammerlab.org](mailto:arman@hammerlab.org)



**Bulent Arman Aksoy**

Medical University of South Carolina



OPEN  ACCESS



**DOI:** [dx.doi.org/10.17504/protocols.io.ba8gihtw](https://dx.doi.org/10.17504/protocols.io.ba8gihtw)

**Protocol Citation:** Bulent Arman Aksoy, Pinar Aksoy, Elinor Gottschalk, Jeff Hammerbacher 2020. Measuring relative reactivity of mouse TCRs against a mouse cancer cell line. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.ba8gihtw>

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** January 10, 2020

**Last Modified:** January 13, 2020

**Protocol Integer ID:** 31720

**Keywords:** TCR, reactivity, bioactivity, T cell, cytotoxicity, OT-I, mouse, MC38B



## Abstract

This protocol repurposes Promega's T Cell Activation Bioassay workflow to be able to test relative mouse TCR reactivity against a cell line. This specific protocol uses MC38 as the target as it doesn't normally present SIINFEKL and have good H2Kb and H2Db expression levels. The reactivity will be in relative to the positive control (OT-I reactivity against SIINFEKL-pulsed cells) and the negative control (OT-I reactivity against unpulsed cells).

## Materials

### MATERIALS

- ✕ EcoRI-HF - 10,000 units **New England Biolabs Catalog #R3101S**
- ✕ NotI-HF - 2,500 units **New England Biolabs Catalog #R3189L**
- ✕ HiScribe T7 ARCA mRNA Kit (with Tailing) - 20 rxns **New England Biolabs Catalog #E2060S**
- ✕ UltraPure™ Phenol:Chloroform:Isoamyl Alcohol (25:24:1, v/v) **Thermo Fisher Scientific Catalog #15593031**
- ✕ Nuclease-Free Water
- ✕ Ethanol (100%, Molecular Biology Grade) **Fisher Scientific Catalog #BP2818500**
- ✕ T Cell Activation Bioassay (NFAT) **Promega Catalog #J1621**
- ✕ SpectraMax i3 Multi-Mode Microplate Detection Platform **Molecular Devices Catalog #i3x**
- ✕ Corning™ RPMI 1640 Medium (Mod.) 1X with L-Glutamine **Fisher Scientific Catalog #MT10041CV**
- ✕ Fetal Plus® **Atlas Biologicals Catalog # FP-0500-A**
- ✕ Penicillin-Streptomycin (10,000 U/mL) **Thermo Fisher Scientific Catalog #15140122**
- ✕ Neon™ Transfection System **Thermo Fisher Scientific Catalog #MPK5000**
- ✕ Chloroform **Fisher Scientific Catalog #C298-4**
- ✕ Neon™ Transfection System 100 µL Kit **Thermo Fisher Catalog #MPK10096**
- ✕ pcDNA3.1( )-OTI-TCRA **addgene Catalog #131035**
- ✕ pcDNA3.1( )-OTI-TCRB **addgene Catalog #131036**
- ✕ Cd8a (NM\_001081110) Mouse Tagged ORF Clone **OriGene Catalog #MR227539**
- ✕ Cd8b1 (NM\_009858) Mouse Tagged ORF Clone **OriGene Catalog #MR225204**
- ✕ CELL CULTURE MICROPLATE 96 WELL PS F-BOTTOM (CHIMNEY WELL) WHITE CELLSTAR® TC LID WITH CONDENS **greiner bio-one Catalog #655083**

### STEP MATERIALS

- ✕ EcoRI-HF - 10,000 units **New England Biolabs Catalog #R3101S**
- ✕ NotI-HF - 2,500 units **New England Biolabs Catalog #R3189L**
- ✕ pcDNA3.1( )-OTI-TCRA **addgene Catalog #131035**
- ✕ pcDNA3.1( )-OTI-TCRB **addgene Catalog #131036**
- ✕ Cd8a (NM\_001081110) Mouse Tagged ORF Clone **OriGene Catalog #MR227539**
- ✕ Cd8b1 (NM\_009858) Mouse Tagged ORF Clone **OriGene Catalog #MR225204**
- ✕ T Cell Activation Bioassay (NFAT) **Promega Catalog #J1621**
- ✕ SpectraMax i3 Multi-Mode Microplate Detection Platform **Molecular Devices Catalog #i3x**

## Protocol materials

- ✕ Neon™ Transfection System **Thermo Fisher Scientific Catalog #MPK5000**
- ✕ Cd8b1 (NM\_009858) Mouse Tagged ORF Clone **OriGene Catalog #MR225204**
- ✕ UltraPure™ Phenol:Chloroform:Isoamyl Alcohol (25:24:1, v/v) **Thermo Fisher Scientific Catalog #15593031**
- ✕ Neon™ Transfection System 100 µL Kit **Thermo Fisher Catalog #MPK10096**
- ✕ Fetal Plus® **Atlas Biologicals Catalog # FP-0500-A**
- ✕ SpectraMax i3 Multi-Mode Microplate Detection Platform **Molecular Devices Catalog #i3x**
- ✕ Ethanol (100%, Molecular Biology Grade) **Fisher Scientific Catalog #BP2818500**
- ✕ EcoRI-HF - 10,000 units **New England Biolabs Catalog #R3101S**
- ✕ Chloroform **Fisher Scientific Catalog #C298-4**
- ✕ Cd8a (NM\_001081110) Mouse Tagged ORF Clone **OriGene Catalog #MR227539**
- ✕ Penicillin-Streptomycin (10,000 U/mL) **Thermo Fisher Scientific Catalog #15140122**
- ✕ Cd8b1 (NM\_009858) Mouse Tagged ORF Clone **OriGene Catalog #MR225204**
- ✕ NotI-HF - 2,500 units **New England Biolabs Catalog #R3189L**
- ✕ pcDNA3.1( )-OTI-TCRA **addgene Catalog #131035**
- ✕ T Cell Activation Bioassay (NFAT) **Promega Catalog #J1621**
- ✕ Corning™ RPMI 1640 Medium (Mod.) 1X with L-Glutamine **Fisher Scientific Catalog #MT10041CV**
- ✕ Nuclease-Free Water
- ✕ SpectraMax i3 Multi-Mode Microplate Detection Platform **Molecular Devices Catalog #i3x**
- ✕ NotI-HF - 2,500 units **New England Biolabs Catalog #R3189L**
- ✕ pcDNA3.1( )-OTI-TCRA **addgene Catalog #131035**
- ✕ pcDNA3.1( )-OTI-TCRB **addgene Catalog #131036**
- ✕ EcoRI-HF - 10,000 units **New England Biolabs Catalog #R3101S**
- ✕ pcDNA3.1( )-OTI-TCRB **addgene Catalog #131036**
- ✕ T Cell Activation Bioassay (NFAT) **Promega Catalog #J1621**
- ✕ CELL CULTURE MICROPLATE 96 WELL PS F-BOTTOM (CHIMNEY WELL) WHITE CELLSTAR® TC LID WITH CONDENS **greiner bio-one Catalog #655083**
- ✕ HiScribe T7 ARCA mRNA Kit (with Tailing) - 20 rxns **New England Biolabs Catalog #E2060S**
- ✕ Cd8a (NM\_001081110) Mouse Tagged ORF Clone **OriGene Catalog #MR227539**
- ✕ pcDNA3.1( )-OTI-TCRA **addgene Catalog #131035**
- ✕ pcDNA3.1( )-OTI-TCRB **addgene Catalog #131036**



✕ Cd8a (NM\_001081110) Mouse Tagged ORF Clone **OriGene Catalog #MR227539**

✕ Cd8b1 (NM\_009858) Mouse Tagged ORF Clone **OriGene Catalog #MR225204**

✕ NotI-HF - 2,500 units **New England Biolabs Catalog #R3189L**

✕ EcoRI-HF - 10,000 units **New England Biolabs Catalog #R3101S**

✕ T Cell Activation Bioassay (NFAT) **Promega Catalog #J1621**

✕ SpectraMax i3 Multi-Mode Microplate Detection Platform **Molecular Devices Catalog #i3x**

## Before start


Make sure you have some familiarity with


- Plasmid propagation, midi-prepping, restriction, and purification
- In vitro transcription and RNA handling
- Electroporation
- Basic cell culture maintenance


## Preparation of electroporation material

- 1 Order, clone, and midi-prep all the TCR and mouse CD8 plasmids:

 pcDNA3.1( )-OTI-TCRA **addgene Catalog #131035**

 pcDNA3.1( )-OTI-TCRB **addgene Catalog #131036**

 Cd8a (NM\_001081110) Mouse Tagged ORF Clone **OriGene Catalog #MR227539**

 Cd8b1 (NM\_009858) Mouse Tagged ORF Clone **OriGene Catalog #MR225204**

and make sure they are of good quality for further applications.

- 2 Linearize plasmids using the corresponding enzymes right at the end of their inserts





Preferred enzyme for the OT-I plasmids is **EcoRI**:



 EcoRI-HF - 10,000 units **New England Biolabs Catalog #R3101S**

and the preferred enzyme for the mouse CD8s is **NotI**:

 NotI-HF - 2,500 units **New England Biolabs Catalog #R3189L**

We recommend the following restriction reaction:














-  50 µg of plasmid DNA
-  25 µL of the corresponding restriction enzyme
-  25 µL of the CutSmart Buffer (10X)
- Top the reaction with nuclease-free water to  250 µL

Incubate the reaction at  37 °C for at least  01:00:00 .



### Note

It is very important to fully linearize the plasmid to prevent potential off-running mRNAs. Based on the incubation time, the amount of enzyme can be reduced but the linearization should always be quality-checked via running the product on agarose gel when in doubt.

- 3 Extract the linearized DNA via the standard phenol:chloroform extraction protocol:


1. Add  250  $\mu\text{L}$  of nuclease-free water so that the final volume for the restriction reaction is  500  $\mu\text{L}$
  2. Add  500  $\mu\text{L}$  of phenol:chloroform and vortex well
  3. Spin at  14000 rpm, 4°C, 00:05:00
  4. Transfer the (top) aqueous layer to a new tube, add  500  $\mu\text{L}$  chloroform, and vortex well
  5. Spin at  14000 rpm, 4°C, 00:05:00
  6. Transfer the (top) aqueous layer to a new tube, add  1000  $\mu\text{L}$  absolute EtOH, and mix well by inverting the tube a few times
  7. Keep the sample at  -20 °C for at least  00:30:00
  8. Spin at  14000 rpm, 4°C, 00:30:00
  9. Discard the supernatant without losing the pellet and add  500  $\mu\text{L}$  70% EtOH to wash the pellet
  10. Spin at  14000 rpm, 4°C, 00:10:00
  11. Discard the supernatant, remove all the residual alcohol, and resuspend the pellet in  50  $\mu\text{L}$  nuclease-free water
  12. Quality check the final DNA solution and estimate the concentration via Nanodrop.
- 4 In vitro transcribe mRNA using the linearized templates using NEB's **mRNA synthesis** and LiCl isolation protocols.

#### Note


We recommend starting with  10  $\mu\text{g}$  of linearized template for each product and scaling the NEB's recommended reaction by 10X. Our preferred final elution volume is  250  $\mu\text{L}$ , which should yield **1.5 - 2  $\mu\text{g}/\mu\text{L}$**  mRNA.

Store the IVT'ed mRNA at  -80 °C for future use.

## Culturing and expanding effector and target cells (Day -3)




- 5 Thaw the Jurkat-NFAT cells that come with the T cell bioactivity kit:
-  T Cell Activation Bioassay (NFAT) **Promega Catalog #J1621**



and culture them by seeding 5 million cells in  50 mL of Jurkat media within a T75 flask for at least 3 days or until they reach a density of **1.5 million cells per mL**.

#### Note

Jurkat media:

-  500 mL of **Corning™ RPMI 1640 Medium (Mod.) 1X with L-Glutamine**
-  50 mL of **Fetal Plus®**
-  5 mL of **Penicillin-Streptomycin (10,000 U/mL)**

- 6 Thaw and start culturing MC38 cells. Seeding ~2 million cells in a T75 flask and culturing them for at least three days should yield enough cells for the co-culture.


#### Note

We have been using MC38 cells as our effectors but the choice of cell line is up to the experimenter. Ideally, the cell line doesn't present the SIINFEKL peptide on its own (without pulsing) so that we can use it as a negative control when co-cultured with the OT-I TCR. This protocol assumes, the cell line is of adherent nature so any suspension cell line could require some customization.





## Co-culture setup (Day 0)

- 7 Electroporate Jurkats with mouse CD8 and OT-I subunits

### 7.1

Fill 4 wells of a 6-well culture plate with  6 mL of warm Jurkat media. We will be using this plate as our recovery plate after the electroporation.

### 7.2

1. Collect 20 million Jurkats
2. Spin them down at  350 x g, 4°C, 00:05:00
3. Re-suspend them in  25 mL of PBS (first wash)
4. Spin them down at  350 x g, 4°C, 00:05:00
5. Re-suspend them in  25 mL of PBS (second wash)





6. Re-suspend them in  of R buffer

### 7.3 Add

1.  of mouse CD8A mRNA (~  )
2.  of mouse CD8B mRNA (~  )
3.  of mouse OT-I alpha mRNA (~  )
4.  of mouse OT-I beta mRNA (~  )

onto the  Jurkat cell suspension in R buffer and make sure you mix them well

- 7.4 Using Neon's the 100 uL tips, electroporate cells at **1350 V 10 ms 3 pulse** setting. Use one tip for one reaction. Change the E2 solution every 6 electroporation reactions. Recover at most two reactions within each 6-well-plate well. Electroporate at least 8 reactions, which should yield 16 million Jurkat cells.

#### Equipment

new equipment

NAME

Thermo Fisher Scientific Neon™ Transfection System

BRAND

MPK5000S

SKU

- 7.5 Let the electroporated Jurkat cells recover for at least 12 hours.

- 8 Seed MC38s into two solid-bottom white plates through 2-fold dilutions:

	1	2	3	4	5	6	7	8	9	10	11	12
A	-	-	-	-	-	-	-	-	-	-	-	-
B	-	-	R1	R1	R1	R1	R1	R1	R1	R1	-	-
C	-	-	R2	R2	R2	R2	R2	R2	R2	R2	-	-
D	-	-	R3	R3	R3	R3	R3	R3	R3	R3	-	-
E	-	-	R4	R4	R4	R4	R4	R4	R4	R4	-	-
F	-	-	R5	R5	R5	R5	R5	R5	R5	R5	-	-
G	-	-	R6	R6	R6	R6	R6	R6	R6	R6	-	-
H	-	-	-	-	-	-	-	-	-	-	-	-
	MC38B	100K	50K	25K	12.5K	6.3K	3.2K	1.6K	0.8K			

Plate set up for the target cells. Have at most 6 replicates (row-wise) and prepare 2-fold serial dilutions for each replicate. We will assume that the cells will replicate once overnight so the numbers will double on the day of the co-culture.

### Note

Although the final volume doesn't matter that much since we will be aspirating the media before setting the co-culture, when in doubt you can go with 60  $\mu$ L of media per well. To have 100K cells in 60  $\mu$ L of media, the serial dilution should start roughly at 3.3 million cells per mL concentration and we will be needing roughly 1.5 million cells per plate.

8.1 Let the targets cells attach and settle down for at least 12 hours.

## Assay luciferase activity (Day 1)

9 Replenish Jurkats with fresh media:

1. Collect and combine all electroporated cells into a single 50-mL falcon tube
2. Spin 350 x g, 4°C, 00:05:00
3. Discard the supernatant
4. Re-suspend in 10 mL of fresh Jurkat media (~ 1.3 million cells per mL)
5. Split the cell suspension into two (5 mL each) 15-mL falcon tubes
6. Label and pulse one of the tubes with the **SIINFELK** peptide at 10 micromolar ( $\mu$ M) (9.63 ug/mL)

10 Aspirate the media from the MC38-seeded plates

11 Add 75  $\mu$ L (~100K) of the electroporated Jurkats onto each well. Label the plates as **pulsed** or **unpulsed** accordingly.







## Note

To reduce the chances of cross-contamination, always start by setting up the co-culture for the unpulsed condition.

This should give us the following plate setup:

	1	2	3	4	5	6	7	8	9	10	11	12
A	-	-	-	-	-	-	-	-	-	-	-	-
B	-	-	R1	R1	R1	R1	R1	R1	R1	R1	-	-
C	-	-	R2	R2	R2	R2	R2	R2	R2	R2	-	-
D	-	-	R3	R3	R3	R3	R3	R3	R3	R3	-	-
E	-	-	R4	R4	R4	R4	R4	R4	R4	R4	-	-
F	-	-	R5	R5	R5	R5	R5	R5	R5	R5	-	-
G	-	-	R6	R6	R6	R6	R6	R6	R6	R6	-	-
H	-	-	-	-	-	-	-	-	-	-	-	-
		MC38B	200K	100K	50K	25K	12.5K	6.3K	3.2K	1.6K		
		Jurkat	100K	100K	100K	100K	100K	100K	100K	100K		

Plate setup for the coculture condition. Jurkat concentration is kept fixed but the target cells are titrated down from 2:1 target:effector ratio using 2-fold serial dilution.

- 12 Co-culture for at least  06:00:00 .
- 13 Take the plates and the luciferase substrates out and let them equilibrate at  Room temperature for  00:10:00 .
- 14 Add  75  $\mu$ L of the luciferase reagent onto each well and let the reactions run for at  00:10:00 .
- 15 Measure the luciferase activity using a standard plate reader with luminescence reading capability  
 SpectraMax i3 Multi-Mode Microplate Detection Platform **Molecular Devices Catalog #i3x**

## Dataset

OT-I TCR reactivity against SIINFEKL-pulsed or -unpulsed MC38s


NAME

<https://docs.google.com/spreadsheets/d/1exKy0eE89bW9RymnPyf7ww-TIs1-V3EMyZ224h-xs58/edit?usp=sharing>

LINK

## Expected result

OT-I TCR reactivity against MC38B (proxied via Jurkat-NFAT luminescence assay)



Median luminescence at different target:effector ratios across pulsed (positive) and unpulsed (negative) control samples. Each TCR that will be tested will produce a reactivity metric relative to these two controls.