



Jun 04, 2018

Version 2

# Measuring spectral reflectance and transmittance (350-2500 nm) of small and/or narrow leaves using the Spectra Vista Corporation (SVC) DC-R/T Integrating Sphere V.2

 Version 1 is forked from [Measuring spectral reflectance and transmittance \(350-2500 nm\) of large leaves using an integrating sphere](#)



DOI

[dx.doi.org/10.17504/protocols.io.qp7dvrn](https://dx.doi.org/10.17504/protocols.io.qp7dvrn)

Etienne Laliberté<sup>1</sup>, Raymond Soffer<sup>2</sup>

<sup>1</sup>Université de Montréal; <sup>2</sup>National Research Council

Canadian Airborne Biodiversity Observatory  
Tech. support email: [jocelyne.ayotte@umontreal.ca](mailto:jocelyne.ayotte@umontreal.ca)



Etienne Laliberté

Université de Montréal

## Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.qp7dvrn>



**External link:** [www.caboscience.org](http://www.caboscience.org)

**Protocol Citation:** Etienne Laliberté, Raymond Soffer 2018. Measuring spectral reflectance and transmittance (350-2500 nm) of small and/or narrow leaves using the Spectra Vista Corporation (SVC) DC-R/T Integrating Sphere. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.qp7dvrn>

**Manuscript citation:**

Noda, H. M., T. Motohka, K. Murakami, H. Muraoka, and K. N. Nasahara. 2013. Accurate measurement of optical properties of narrow leaves and conifer needles with a typical integrating sphere and spectroradiometer. Plant, Cell & Environment 36:1903–1909. <https://doi.org/10.1111/pce.12100>

**License:** This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** In development

**We are still developing and optimizing this protocol.**

**Created:** June 04, 2018

**Last Modified:** June 04, 2018

**Protocol Integer ID:** 12767

**Keywords:** leaf spectral reflectance, leaf spectroscopy protocol, canadian airborne biodiversity observatory, measuring spectral reflectance, reflectance measurement, spectral reflectance, measurements of adaxial reflectance, individual leaf, reflectance, carnegie airborne observatory, range field spectroradiometer, leaf, canopy plant, adaxial reflectance, narrow leaf, spectra vista corporation, sphere user manual, integrating sphere user manual, transmittance

## Abstract

Here we describe the standardised protocol used by the **[Canadian Airborne Biodiversity Observatory](#)** (CABO) to measure leaf spectral reflectance and transmittance, using an integrating sphere fitted to a portable full-range field spectroradiometer, for the special case where an individual **leaf is too small and/or too narrow** to entirely cover the reflectance or transmission port of the integrating sphere. Briefly, three arrays of mature, healthy and sunlit leaves from a canopy plant are arranged on a custom sample mount, and are then used for measurements of adaxial reflectance and transmittance. Reflectance measurements are referenced to a calibrated Spectralon® disk and corrected for stray light. Our leaf spectroscopy protocol is an adaptation of that of **[Noda et al. \(2013\)](#)** to the SVC DC-R/T sphere, and also builds from protocols by the **[Carnegie Airborne Observatory](#)** and the **[SVC](#)** integrating sphere user manual.



## Guidelines

### Handling Spectralon®

- **Do not touch Spectralon®** (e.g. sphere interior, reference disks, plugs) with your fingers.
- **Do not use canned air** to remove dust on the Spectralon® disk; canned air contains chemicals that can alter Spectralon®'s optical properties.
- **Do not attempt to clean Spectralon®** in the field, other than **blowing surface dust only on the Spectralon® reference disk or sphere plugs** using the Canless Air Duster System; cleaning Spectralon® requires a special procedure that should only be done in the lab.

### Equipment

- Spectra Vista Corporation **HR-1024i** full-range (350-2500 nm) field spectroradiometer
- Spectra Vista Corporation 3-inch Spectralon® DC-R/T **Sphere**
- Semi-rugged laptop or PDA running the SVC Scan software
- **Canless Air Duster System O<sub>2</sub> Hurricane** (**never use canned air**) to remove dust from the surface of the Spectralon® reference disk
- Plastic containers with lids to temporarily store leaf arrays during measurements (optional)

### Consumables

- Nitrile gloves for handling leaves
- Whatman No. 2 filter paper (110 mm diameter)
- Acetate sheets (to make thin plastic sample mounts)
- Manila file folders (to make thin cardboard sample platforms)
- Scotch™ Magic Tape

## Materials

### MATERIALS

- ✂ Manila File Folders **Staples Catalog #116723**
- ✂ Write-On Transparency Film **Staples Catalog #954144**
- ✂ Whatman™ Qualitative Filter Paper: Grade 2 Circles (110 mm diameter) **Fisher Scientific Catalog #09-810E**
- ✂ Scotch™ Magic Tape **Staples Catalog #14172**

## Troubleshooting

## Safety warnings

- ! The lamp of the integrating sphere can get **very hot** and should be handled from its slotted base to avoid burns.



## Before start

Consult the user manual of the spectroradiometer and the integrating sphere to set up the instrument.

The instrument should be set up in the shade, sheltered as much as possible from the elements.

All canopy plants selected for measurements should have already been tagged, identified, and georeferenced before spectroscopy measurements start.

The spectroscopist should be positionned as close as possible to the sampled plants to minimise time from collection to measurement.

The spectroscopist should be in a comfortable position and have enough room around the instrument to spread leaf samples around without the risk of mixing up individual leaves during handling.

**Mature, fully-developped, healthy-looking leaves from the sunlit (>3 h per day of direct sunlight) portion of the canopy are selected** for spectral measurements from the bulk leaf sample (often one of a few branches). Leaves should be collected from the uppermost surface of the branch (i.e. receiving the most direct sunlight).



## Instrument Set-Up

- 1 Install the integrating sphere onto the spectroradiometer.

### Note

Follow the SVC integrating sphere manual p. 6-16.

- 2 Position the lamp over the sphere **primary light entrance port**.


### Note

Make sure lamp is secured in locked position.


### Safety information

The lamp can get very hot. Grab it by the slotted heat shield.

- 3 Power the spectroradiometer on and **warm up for >15 min.**

 00:15:00 Spectroradiometer warm-up period

- 4 Power the integrating sphere lamp and **warm up for >5 min.**

 00:05:00 Integrating sphere warm-up period

- 5 Screw the tethered light trap on the **transmission port** sample holder.

#### Note

The light trap can stay on the sample holder for the entire measurement session.

- 6 Check lamp alignment.

#### Note

Use a thin piece of paper at the exit of the reflectance sample port (empty port) to ensure the light beam under-fills and is centered in the reflectance port. **If it is not, then proceed to lamp alignment** as described in the SVC integrating sphere user manual, p. 29-32.

## Software Set-Up

- 7 Ensure that the spectroradiometer is connected to the computer via Bluetooth.

#### Note

The default Bluetooth password for SVC HR-1024i spectroradiometers is **hr1024i**.

- 8 Launch SVC Scan.



- 9 Select **auto-save mode** using *File > Data Options > Auto Save*.


- 10 Select *File > New* to create a new working folder in your shared '**spectra**' Google Drive CABO folder.

#### Note

If you do not yet have a shared Google Drive folder for your project, contact the CABO data manager to create one ([etienne.laliberte@umontreal.ca](mailto:etienne.laliberte@umontreal.ca) or [jeremy.goimard@umontreal.ca](mailto:jeremy.goimard@umontreal.ca)).

Your shared CABO Google Drive folder should be named *YYYYstarted-YourLastName-Degree-University*, e.g. 2018-Girard-MSc-UdeM.


- 11 Create a new folder named *YYYY-MM-DD-SiteID-SpectroradiometerSerialNumber* (without spaces) within that '**spectra**' folder, and keep the default file base name (= *grMMDDYY*).



**Note**

The base file name is the prefix for all files that will be saved within that folder. In auto-save mode, each file will be added an incremental step, e.g., *gr053018\_0000.sig*, *gr053018\_0001.sig*, etc.


- 12 Go to *Control > Setup Instrument...* to connect the spectroradiometer.

**Note**

Different instruments may be configured under different virtual COM ports.

- 13
1. Choose the *RAW DN* optic
  2. Select 5 sec scan time
  3. Select *Auto Integration*
  4. Ensure the date and time are correct
  5. Check *Integration Scales Raw DN Data*




Setup
✕

**Computer Port Control**

☒ List Only Configured Ports
 

COM4

☐ List All Possible Ports
 

COM4

Disconnect

Connect

**HR-1024i Settings**

Optic:

RAW DN

**Scan Timing**

☒ Specify Total Scan Time (Sec):
 

5

☐ Specify Coadds:
 

Si: 4

Swir1: 32

Swir2: 32

**Integration Time (mSec)**

☒ Auto Integration
 

Si: 50

Swir1: 20

Swir2: 6

Date: 06/04/18

Time: 11:18:30

Stored Scans: 0


Target Photo Acquisition: None

☒ Integration Scales RAW DN Data

Cancel

OK

- 14 In *Control* > *Overlap / Matching Settings*:
1. Select *Preserve Overlapped Detector Data*
  2. Select *Matching Type* > *None*



**Overlap / Matching Settings**

**Detector Overlap Controls**

☒ Preserve Overlapped Detector Data  
☐ Remove Overlapped Detector Data

Si->Swir1 Transition Wvl: 990  
Swir1->Swir2 Transition Wvl: 1900

**Detector Matching Controls**

**Matching Type**  
☒ None  
☐ Radiance  
☐ Reflectance

**Matching Region**  
Beginning Wvl: 980  
Ending Wvl: 1000

**Matching Factor Limits**  
Min: 0.80 Max: 1.20

**NIR-SWIR Overlap Algorithm**  
☒ Use The NIR-SWIR Overlap Algorithm Within The Matching Region

Set Default Values Cancel OK


Review Protocol Summary Diagram

- 15 Review the **document** summarising the different sphere configurations (A–I), and the scans that need to be recorded in each configuration.



## Prepare Sample Mounts and Holders

- 16 Build a few sample mounts from thin plastic (e.g. acetate sheets).




### Note

The sample mount is made of a thin plastic sheet with a square window that is larger than the edge of the sample port lip (e.g. 2.5 cm × 2.5 cm).

### Note

Fig. 3a from Noda et al. (2013). <https://doi.org/10.1111/pce.12100>

- 17 Build a few sample platforms from thin cardboard (e.g. file folders).



#### Note


The sample platform is made of cardboard and will eventually attach to the integrating sphere to create a flat surface.

#### Note

Fig. 3c from Noda et al. (2013). <https://doi.org/10.1111/pce.12100>

## Prepare Leaf Arrays

- 18 Fix leaves onto the sample mount with tape to make leaf array #1.



#### Note

Fig. 3b from Noda et al. (2013). <https://doi.org/10.1111/pce.12100>

#### Note

Try to position as much leaf material as possible in the center of the port where the intensity of the light beam will be highest.

Leaves should be separated by about 1 mm to avoid multiple scattering among them.

If leaves are shorter than the sample mount hole (i.e. <2.5 cm long), fix leaves by their petiole with tape in one (or two) rows.

- 19 Fix leaves onto the sample mount with tape to make leaf array #2.

**Note**

Try to position as much leaf material as possible in the center of the port where the intensity of the light beam will be highest.

Leaves should be separated by about 1 mm to avoid multiple scattering among them.

If leaves are shorter than the sample mount hole (i.e. <2.5 cm long), fix leaves by their petiole with tape in one (or two) rows.

- 20 Fix leaves onto the sample mount with tape to make leaf array #3.

**Note**

Try to position as much leaf material as possible in the center of the port where the intensity of the light beam will be highest.

Leaves should be separated by about 1 mm to avoid multiple scattering among them.

If leaves are shorter than the sample mount hole (i.e. <2.5 cm long), fix leaves by their petiole with tape in one (or two) rows.

- 21 Fix sample mount holding leaf array #1 onto the sample platform with tape.

**Note**

The sample mount holding the leaf array should be sandwiched between the two sides of the cardboard sample platform.

- 22 Fix sample mount holding leaf array #2 onto the sample platform with tape.

**Note**

The sample mount holding the leaf array should be sandwiched between the two sides of the cardboard sample platform.

- 23 Fix sample mount holding leaf array #3 onto the sample platform with tape.

**Note**

The sample mount holding the leaf array should be sandwiched between the two sides of the cardboard sample platform.

## Configuration A: Reflectance Mode, Leaf Array, Reference

- 24 Screw the tethered light trap on the **reflectance port** sample holder.

### Note

The light trap can stay on the sample holder for the entire measurement session.


- 25 Place the tethered calibrated Spectralon® reflectance standard over the **reflectance port**.

### Note

Place the standard over the reflectance port so that the light beam shines directly on its reflective surface (= facing inside of the sphere).

- 26 Position leaf array #1 over the **transmission port** so that the **adaxial** (upper) surfaces of the leaves face into the sphere.

- 27 Collect a '**Reference Scan**' in this configuration.



#### Note

This corresponds to the **reference radiance** in reflectance mode for the **leaf array** ( $R_{\text{ref},a}$ ). The reference data will be automatically saved in all successive target scan files until a new 'Reference Scan' is made.

- 28 Collect a '**Target Scan**' in the same configuration and **save the file**.


### Configuration B: Reflectance Mode, Leaf Array + Paper, Reference

- 29 Place a filter paper directly behind the sample platform holding leaf array #1.

#### Note

The leaf array #1 should remain in the same position.

- 30 Collect a '**Target Scan**' in this configuration and **save the file**.



#### Note

This corresponds to the **reference radiance** for the **leaf array + paper** in reflectance mode ( $R_{\text{ref,ap}}$ ).

### Configuration C: Reflectance Mode, Paper, Reference


- 31 Carefully remove the sample platform holding leaf array #1 from the transmission port.

#### Note

The filter paper should remain in the same position.

- 32 Collect a '**Target Scan**' in this configuration and **save the file**.



**Note**

This corresponds to the **reference radiance** for the **filter paper** in reflectance mode ( $R_{\text{ref,p}}$ ).


**Configuration D: Reflectance Mode, Stray light**

- 33 Remove the filter paper from the transmission port sample holder.
- 34 Remove the tethered calibrated Spectralon® reflectance standard from the reflectance port.

35 Place the tethered calibrated Spectralon® reflectance standard over the **transmission port**.

36 Collect a '**Target Scan**' in this configuration and **save the file**.

**Reflectance mode**  
Stray light radiance ( $R_{str}$ )



**Note**

This corresponds to the **stray light radiance** in reflectance mode ( $R_{\text{str}}$ ).


**Configuration E: Reflectance Mode, Paper, Target**

- 37 Place the filter paper over the **reflectance port**.

**Note**

Focus the measurements on the same area of the paper used in previous measurements.

- 38 Collect a '**Target Scan**' in this configuration and **save the file**.




**Note**

This corresponds to the **target radiance** of the **filter paper** in reflectance mode ( $R_{\text{tar,p}}$ ).

**Configuration F: Reflectance Mode, Leaf Array #1 + Paper, Target**

- 39 Position the sample platform containing leaf array #1 in front the filter paper over the **reflectance port** with the **adaxial** (upper) surfaces of the leaves face into the sphere.
- 40 Collect a '**Target Scan**' in this configuration and **save the file**.

**Note**

This corresponds to the **target radiance** of **leaf array #1 + paper** in reflectance mode ( $R_{\text{tar,ap,1}}$ ).


## Configuration G: Reflectance Mode, Leaf Array #1, Target

- 41 Remove the filter paper from the reflectance port.

### Note

The sample platform holding leaf array #1 should remain in the same position.

- 42 Collect a '**Target Scan**' in this configuration and **save the file**.




### Note

This corresponds to the **target radiance** of **leaf array #1** in reflectance mode ( $R_{tar,a,1}$ ).

## Configuration F: Reflectance Mode, Leaf Array #2 + Paper, Target

- 43 Carefully replace leaf array #1 by leaf array #2.
- 44 Place the filter paper directly behind leaf array #2 over the **reflectance port**.
- 45 Collect a '**Target Scan**' in this configuration and **save the file**.



#### Note

This corresponds to the **target radiance of leaf array #2 + paper** in reflectance mode ( $R_{\text{tar,ap},2}$ ).


### Configuration G: Reflectance Mode, Leaf Array #2, Target

- 46 Remove the filter paper from the reflectance port.

Note

The sample platform holding leaf array #2 should remain in the same position.

- 47 Collect a '**Target Scan**' in this configuration and **save the file**.




Note

This corresponds to the **target radiance** of **leaf array #2** in reflectance mode ( $R_{tar,a,2}$ ).

## Configuration F: Reflectance Mode, Leaf Array #3 + Paper, Target

- 48 Carefully replace leaf array #2 by leaf array #3.
- 49 Place the filter paper directly behind leaf array #3 over the **reflectance port**.
- 50 Collect a '**Target Scan**' in this configuration and **save the file**.



#### Note

This corresponds to the **target radiance** of **leaf array #3 + paper** in reflectance mode ( $R_{\text{tar,ap},3}$ ).

### Configuration G: Reflectance Mode, Leaf Array #3, Target


- 51 Remove the filter paper from the reflectance port.



**Note**

The sample platform holding leaf array #3 should remain in the same position.

- 52 Collect a '**Target Scan**' in this configuration and **save the file**.

**Note**

This corresponds to the **target radiance** of **leaf array #3** in reflectance mode ( $R_{tar,a,3}$ ).

**Configuration A: Reflectance Mode, Leaf Array, Reference**

- 53 Remove leaf array #3 from the reflectance port.


- 54 Place the tethered calibrated Spectralon® reflectance standard over the **reflectance port**.

**Note**

Place the standard over the reflectance port so that the light beam shines directly on its reflective surface (= facing inside of the sphere).

- 55 Position leaf array #1 over the **transmission port** so that the **adaxial** (upper) surfaces of the leaves face into the sphere.

- 56 Collect a '**Target Scan**' in this configuration and **save the file**.



#### Note

This second reference radiance measurement is only used to assess system stability in reflectance mode.

### Configuration H: Transmittance Mode, Reference

- 57 Remove leaf array #1 from the **transmission port**.
- 58 Position the sample platform containing leaf array #1 over the **reflectance port** with the **abaxial** (lower) surface of the leaves facing into the sphere.
- 59 Remove the tethered calibrated Spectralon® reflectance standard from the sphere **transmission port**.



60 Remove the light trap from the transmission port sample holder.

61 Position the lamp over the sphere **transmission port**.

Note

Make sure lamp is secured in locked position.

Safety information

The lamp can get very hot. Grab it by the slotted heat shield.

62 Install the Spectralon® plug over the **primary light port**.

Note


Ensure that the curved plug is placed the correct way to match the curvature of the sphere.

63 Collect a '**Reference Scan**' in this configuration.

## Transmittance mode

### Leaf array

### Reference radiance ( $T_{ref}$ )



#### Note

This corresponds to the **reference radiance** in transmittance mode for the **leaf array** ( $T_{ref}$ ). The reference data will be automatically saved in all successive target scan files until a new 'Reference Scan' is made.

64 Collect a '**Target Scan**' in this configuration and **save the file**.

## Configuration I: Transmittance Mode, Target

65 Carefully remove leaf array #1 from the **reflectance port**.



66 Gently pull lamp and transmission port sample holder away from the sphere.


#### Safety information

The lamp can get very hot. Grab it by the slotted heat shield.

67 Place the sample platform holding leaf array #1 over the **transmission port** with the abaxial (lower) surfaces of the leaves facing into the sphere.

68 Release the transmission sample holder and move lamp back to its locked position.

69 Collect a '**Target Scan**' under this configuration and **save the file**.

**Note**

This corresponds to the **target radiance** in **transmittance** mode for leaf array #1 ( $T_{\text{tar},1}$ ).

- 70 Carefully replace leaf array #1 by leaf array #2.
- 71 Collect a '**Target Scan**' under this configuration and **save the file**.

**Note**

This corresponds to the **target radiance** in **transmittance** mode for leaf array #2 ( $T_{\text{tar},2}$ ).

72 Carefully replace leaf array #2 by leaf array #3.

73 Collect a '**Target Scan**' for leaf array #3 in this configuration and **save the file**.

**Note**

This corresponds to the **target radiance** in **transmittance** mode for leaf array #3 ( $T_{\text{tar},3}$ ).

## Configuration H: Transmittance Mode, Reference

74 Remove leaf array #3 from the transmission sample port holder.

75 Position the sample platform containing leaf array #1 over the **reflectance port** with the **abaxial** (lower) surface of the leaves facing into the sphere.


76 Collect a '**Target Scan**' in this configuration and **save the file**.



## Transmittance mode

### Leaf array

### Reference radiance ( $T_{\text{ref}}$ )



### Note

This second reference radiance is only used to assess system stability in transmittance mode.

## Calculating Absolute Adaxial Reflectance of Leaf Array

77 The equation (Noda et al. 2013; eqn. 9) for **adaxial reflectance** of leaf array  $i$ ,  $\rho_{a,i}$  is

$$\rho_{a,i} = [(R_{\text{tar},a,i} - R_{\text{str}}) \div (R_{\text{ref},a} - R_{\text{str}})] \times \rho_{\text{ref}} \times [1 \div (1 - G_{r,i})]$$

where

$R_{tar,a,i}$  is the target radiance of leaf array  $i$  (adaxial side) in reflectance mode,  
 $R_{ref,a}$  is the reference radiance used for all leaf arrays in reflectance mode,  
 $R_{str}$  is the stray light radiance in reflectance mode,  
 $\rho_{ref}$  is the absolute reflectance of the calibrated Spectralon® reflectance standard, and  
 $G_{r,i}$  is the gap fraction in reflectance mode for leaf array  $i$ , which is calculated **at 400 nm** (Noda et al. 2013; eqn. 13) by

$$G_{r,i} = [ ( (R_{tar,ap,i} - R_{str}) \div (R_{ref,ap} - R_{str}) ) - ( (R_{tar,a,i} - R_{str}) \div (R_{ref,a} - R_{str}) ) ] \times ( \rho_{ref} \div \rho_p )$$

where

$R_{tar,ap,i}$  is the target radiance of leaf array  $i$  + filter paper in reflectance mode,  
 $R_{ref,ap}$  is the reference radiance used for all leaf arrays + filter paper in reflectance mode, and  
 $\rho_p$  is the absolute reflectance of the filter paper, which is calculated (Noda et al. 2013; eqn. 3) by

$$\rho_p = [(R_{tar,p} - R_{str}) \div (R_{ref,p} - R_{str})] \times \rho_{ref}$$

where

$R_{tar,p}$  is the target radiance of the filter paper in reflectance mode, and  
 $R_{ref,p}$  is the reference radiance of the filter paper in reflectance mode.

## Calculating Absolute Adaxial Transmittance of Leaf Array

78 The equation for **adaxial transmittance** of leaf array  $i$ ,  $\tau_{a,i}$  is

(modified from eqn. 15 of Noda et al. 2013 for the SVC sphere)

$$\tau_{a,i} = [ (T_{tar,i} \div T_{ref}) - G_{t,i} ] \times [ 1 \div (1 - G_{t,i}) ]$$

where

$T_{tar,i}$  is the target radiance of leaf array  $i$  in transmittance mode,  
 $T_{ref}$  is the reference radiance used for all leaf arrays in transmittance mode,



$G_{t,i}$  is the gap fraction in transmittance mode for leaf array  $i$ , which is calculated **at 400 nm** (since  $\tau_{a,i} = 0$  at that wavelength) by

$$G_{t,i} = T_{\text{tar},i} \div T_{\text{ref}}$$