

Jun 25, 2019

Media and reagents for *Seminavis robusta* cultivation and experiments

 In 2 collections

DOI

<https://dx.doi.org/10.17504/protocols.io.33xgqpn>

Lev Tsypin¹, Aaron Turkewitz²

¹California Institute of Technology; ²University of Chicago

Protist Research to Opti...



Lev Tsypin

Stanford

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.33xgqpn>

Protocol Citation: Lev Tsypin, Aaron Turkewitz 2019. Media and reagents for *Seminavis robusta* cultivation and experiments. protocols.io <https://dx.doi.org/10.17504/protocols.io.33xgqpn>

Manuscript citation:



License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: June 12, 2019

Last Modified: June 25, 2019

Protocol Integer ID: 24407

Keywords: Diatom, marine, protist, cultivation, reagents for seminavis robusta cultivation, seminavis robusta cultivation, seminavis robusta, culturing condition, reagent

Abstract

This protocol describes the culturing conditions used to grow *S. robusta* in the lab for all experiments.

Guidelines

It matters whether you use NanoPure/MilliQ water or glass distilled water. In glass distilled water-based media, *S. robusta* tends to disadhere from the bottom of the culture flask and float more readily than in MilliQ-based media.

Materials

MATERIALS

- ⊗ Bioactif artificial sea salt **Tropic Marin Catalog #10304**
- ⊗ Guillard's (F/2) Marine Water Enrichment Solution **Merck MilliporeSigma (Sigma-Aldrich) Catalog #G9903-500ML**
- ⊗ Sodium bicarbonate **Merck MilliporeSigma (Sigma-Aldrich) Catalog #S6014**
- ⊗ Penicillin G sodium salt **Merck MilliporeSigma (Sigma-Aldrich) Catalog #P3032-25MU**
- ⊗ Gentamycin sulfate **Merck MilliporeSigma (Sigma-Aldrich) Catalog #345814-1GM-M**
- ⊗ Ampicillin sodium salt **Fisher Scientific Catalog #BP1760-25**
- ⊗ Streptomycin sulfate **MP Biomedicals Catalog #100556**
- ⊗ Agar **Fisher Scientific Catalog #BP1423-500**
- ⊗ 175 cm² flasks with sealing caps **Thermo Scientific Catalog #156502**
- ⊗ 25 cm² flasks with sealing caps **USA Scientific Catalog #CC7682-4325**
- ⊗ Stericup **Merck Millipore (EMD Millipore) Catalog #SCGPT05RE**
- ⊗ Cell Scraper **Fisher Scientific Catalog #08-100-241**
- ⊗ D-sorbitol **Merck MilliporeSigma (Sigma-Aldrich) Catalog #M4125**
- ⊗ Mannitol **Merck MilliporeSigma (Sigma-Aldrich) Catalog #S1876**



Preparation of F/2 medium

- 1 Dissolve  34.5 g Tropic-Marin Bioactif (TMB) in  1 L MilliQ water
- 2 Add  80 mg sodium bicarbonate
- 3 Autoclave
- 4 After cooling, add  20 mL 50x Guillard's supplement
- 5 Filter using Stericup into pre-sterilized bottles
- 6 Store at 4C in the dark

Maintenance antibiotics (200x stock)

- 7 Dissolve  1 g penicillin G sodium salt ,  1 g ampicillin sodium salt ,
 0.2 g streptomycin sulfate , and  0.1 g gentamycin sulfate in a final volume of
 20 mL F/2 medium .
- 8 Aliquot and store at -20 C.

Agar plates

- 9 Dissolve  17.25 g TMB in  1 L MilliQ



- 10 Add  80 mg sodium bicarbonate
- 11 Add  0 μ L  15 g agar , dissolve while stirring and heating (may require boiling)
- 12 Autoclave
- 13 After cooling below 60 C, add  20 mL 50x Guillard's supplement , as well as selective antibiotics, if desired.
- 14 Pour approx.  19 mL per plate
- 15 Store in the dark at 4 C

Electroporation buffers

- 16  1 Molarity (M) sorbitol in MilliQ,  0.77 Molarity (M) mannitol in MilliQ, or  0.77 Molarity (M) mannitol with  8 % volume F/2 in MilliQ.