



Aug 26, 2023

🌐 Medium fractionation and EV preparation

DOI

<https://dx.doi.org/10.17504/protocols.io.81wgbx63ylpk/v1>

wusj¹, Nancy C. Hernandez Villegas², schekman¹

¹Department of Molecular and Cell Biology, Howard Hughes Medical Institute, University of California, Berkeley, Berkeley, United States;

²Helen Wills Neuroscience Institute, University of California, Berkeley, Berkeley, United States



Nancy C. Hernandez Villegas

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.81wgbx63ylpk/v1>

Protocol Citation: wusj , Nancy C. Hernandez Villegas, schekman 2023. Medium fractionation and EV preparation. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.81wgbx63ylpk/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: August 01, 2023

Last Modified: August 26, 2023

Protocol Integer ID: 85757

Keywords: synuclein, characterization of the extracellular alpha, extracellular alpha, medium fractionation, ev preparation, ev preparation this protocol

Funders Acknowledgements:

Nancy C Hernandez Villegas

Grant ID: NIH training program T32GM139780

Randy Schekman

Grant ID: Howard Hughes Medical Institute

Abstract

This protocol describes the characterization of the extracellular alpha-synuclein and DNAJC5

Materials

For specific information about reagents and materials to perform an EV preparation, refer to a previous paper published by our lab www.bio-protocol.org/e3706



Medium fractionation

- 1 Centrifuge the conditioned medium at 1,500 x g for  00:20:00 at  4 °C in a Sorvall RC 6+ centrifuge with fixed angle rotor F14S-6X250y FiberLite 20m
- 2 Pour the supernatant into a new container without disturbing the pellet.
- 3 Centrifuge the supernatant at 10,000 x g for  00:30:00 at  4 °C in a Sorvall RC 6+ centrifuge with fixed angle rotor F14S-6X250y FiberLite 30m
- 4 The supernatant fractions at each step were collected and treated with methanol/chloroform to precipitate proteins which were then collected by centrifugation
- 5 Pellet fractions were resuspended in sample buffer to achieve a 20-fold concentration.
- 6 The sedimented fractions at each step were also collected and resuspended in sample buffer.
- 7 All the fractions were analyzed by immunoblot.

Extracellular vesicle (EV) preparation

- 8 Centrifuge the conditioned medium at 1,500 x g for  00:20:00 at  4 °C in a Sorvall RC 6+ centrifuge with fixed angle rotor F14S-6X250y FiberLite 20m
- 9 Pour the supernatant into a new container without disturbing the pellet.
- 10 Centrifuge the supernatant at 10,000 x g for  00:30:00 at  4 °C in a Sorvall RC 6+ centrifuge with fixed angle rotor F14S-6X250y FiberLite 30m
- 11 Pour the supernatant into a new container without disturbing the pellet.



- 12 Add 35 ml of the collected supernatant into a single 38.5 ml ultra-clear tube. Repeat this step 5 times to fill a total of six 38.5 ml ultra-clear tubes.
- 13 Centrifuge at $\sim 100,000 \times g$ (29,500 RPM) for  01:30:00 using a SW32 Ti rotor at  4 °C at maximum acceleration and brake. 1h 30m
- 14 Aspirate the supernatant, put the tubes upside down on paper towel to dry any residue of medium. Aspirate the residue by vacuum system if needed.
- 15 Resuspend the pellet by adding 500 uL of PBS (pH 7.4) buffer to the bottom of each tube, put the tubes on a compatible rack, and rigorously shake the tubes on orbital shaker in cold room for  00:30:00 . 30m
- 16 Combine the resuspend pellet ($\sim 500 \text{ uL} \times 6 \text{ tubes} = 3\text{mL}$; +1mL PBS) and centrifuge again at $120,000 \times g$ (36,500 RPM) x  01:00:00 at  4 °C in SW-55 rotor 1h
- 17 Resuspend the pellet again in 200uL PBS Buffer and add 1mL of 60% (1.8M) sucrose buffer (20mM Tris 7.4, 150 mM NaCl) and vortex to mix the sample evenly. The final sucrose concentration should be above 50% measured by refractometer.
- 18 Aliquots of 40% (5 ml) and 10% (2 ml) sucrose buffer were sequentially overlaid above the sample. The tubes were then centrifuged at $150,000 \times g$ for  16:00:00 at  4 °C in an SW41 Ti swinging-bucket rotor (Beckman Coulter) 16h



- 19 After centrifugation, 0.5 ml fractions were collected from top to bottom and samples were analyzed by SDS-PAGE and immunoblot.