



Aug 27, 2025

Medium throughput protein expression and purification for protein design using heat lysis

DOI

<https://dx.doi.org/10.17504/protocols.io.dm6gpq6b5lzp/v1>

Ida K. Grene^{1,2,3}, Nadja Joachim⁴, Francisca Pinheiro^{1,2,3}, Vili Lampinen^{1,2,3}, Kaare Teilum⁴, Magnus Kjaergaard^{1,2,3}

¹Molecular Biology and Genetics, Aarhus University, Denmark;

²The Danish National Research Foundation Center for Proteins in Memory (PROMEMO), Aarhus University;

³The Danish Research Institute for Translational Neuroscience (DANDRITE), Nordic EMBL Partnership for Molecular Medicine, Aarhus University;

⁴University of Copenhagen



Ida Kjærsgaard Grene

Aarhus University

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.dm6gpq6b5lzp/v1>

Protocol Citation: Ida K. Grene, Nadja Joachim, Francisca Pinheiro, Vili Lampinen, Kaare Teilum, Magnus Kjaergaard 2025.

Medium throughput protein expression and purification for protein design using heat lysis. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.dm6gpg6b5lzp/v1>

Manuscript citation:

Pinheiro F, Nowak JS, Zueva E, Pheasant EC, Grene IK, Lampinen V, et al. Screening de novo designed protein binders in unpurified lysate using flow induced dispersion analysis. *Protein Science*. 2025; 34(10):e70286. <https://doi.org/10.1002/pro.70286>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: May 07, 2025

Last Modified: August 27, 2025

Protocol Integer ID: 217837

Keywords: de novo protein design, protein expression, heat lysis, Medium throughput, screening, bacterial expression, protein purification, small scale, multiwell plates, purification for protein design, medium throughput protein expression, pure lysate sample of the protein, protein design, heat lysis this protocol, using heat lysi, protein, purification step, purification, preparation of heat, using heat, pure lysate sample, stable de novo, lysi, heat

Funders Acknowledgements:

Danish National Research Foundation

Grant ID: DNRF133

Independent Research Fund Denmark

Grant ID: 026-00069B

Lundbeck Foundation

Grant ID: R449-2023-1396

Carlsberg Foundation

Grant ID: CF20-0610

Abstract

This protocols describes medium throughput protein expression in *E. coli* using heat-lysis as a built in purification step. It is intended for preparation of heat-stable *de novo* designed proteins for screening *in vitro*.¹⁻³ It provides a small but relatively pure lysate sample of the protein within 2 working days.

Guidelines

Remember to always keep track of the samples when transferring between plate types.

Gene of interest should be in a pET system or another IPTG inducible expression vector.⁴

Materials

Equipment:

Suggested products will be in parantheses

- Flowbench, appropriate for working with bacteria
- Water bath
- Shaking incubator, preferably with sticky surface for securing plates
- Centrifuge with adaptor for multiwell plates
- Multichannel pipettes taking 1-1000 μL volumes
- Heatblock that can evenly heat 96 wells to 95°C (Eppendorf® ThermoMixer® with deepwell plate smartblock)
- Platerreader for absorption at 600nm
- SDS-page setup

Materials:

Suggested products will be in parantheses

- Plasmid DNA at $\approx 100 \text{ ng}/\mu\text{L}$. (We use synthetic genes cloned into the NdeI-XhoI sites of pET28a+ and stored in 96-well plates)
- Chemically competent *E. coli* BL21 (DE3) strain
- Lennox Broth medium (Pr. liter: 10.0 g Tryptone, 5.0 g Sodium Chloride, 5.0 g Yeast Extract)
- Appropriate antibiotic (e.g. kanamycin for pET28a+)
- ZYM-5052 autoinduction medium⁵
- Sterile reservoirs for multichannel pipettes.
- Sterile 96-well plates (UNIPLATE Collection and Analysis Microplate, 96-well, 2 ml, natural polypropylene, round well bottom)
- Sterile 24-well plates (UNIPLATE Collection and Analysis Microplate, 24-well, 10 ml, natural polypropylene, round well bottom)
- Breathable film (AeraSeal™ film)
- 96-well plate able to withstand 95°C (Eppendorf® Deepwell Plate 96/1000 μL)
- 96-well filter plate (AcroPrep™ Advance filter plates for Lysate Clearance - 1 mL, 3.0 μm Glass Fiber/0.2 μm Supor membrane)
- 96-well clear plate for the plate reader
- Lysis buffer (pH 7.4 phosphate buffer + 300 mM NaCl or similar)
- Polyacrylamide gels, loading dye, appropriate size marker and running buffer
- BCA Protein Assay Kit

If going for chemical lysis:

- Bacterial protein extraction reagent (B-PERT™ Complete Bacterial Protein Extraction Reagent)

If going for purification:

- 96-well IMAC plate (Cytiva His MultiTrap FF)
- Imidazole stock for buffers

Before start

If you have not bought your 96- and 24-well plates, medium and reservoirs sterile, make sure to autoclave them or sterilize in another way before starting.

If keeping the medium in the fridge, bring it to room temperature before adding to the cells.

Prepare the lysis buffer before starting day 2.



Transformation (day 1):

1h 53m

1 Before starting:

15m

- Thaw competent *E. coli* BL21 On ice
- Thaw plasmid DNA Room temperature
- Cool a sterile 96-well plate On ice
- Preheat the water bath 42 °C

2 Add 10 µL *E. coli* BL21 per well in the 96-well plate.

10m

3 Add 1 µL plasmid DNA solution to the bacteria. Stir with pipette tip to connect droplets, do not pipette up and down.

10m

Note

Alternatively, spin the plate after adding the DNA to make sure that cells and DNA come into contact.

4 Incubate the plate 00:15:00 On ice .

15m



5 Heat-shock bacteria in water bath by holding the plate half submerged in the water 42 °C 00:00:30 . Put directly back on ice afterwards.

1m

6 Add 50 µL sterile LB medium to each well.

2m

7 Incubate transformation mixtures 400 rpm, 37°C, 01:00:00 .

1h



Expression (day 1)

12h 22m

8 Prepare autoinduction medium with appropriate antibiotic in a sterile reservoir.

2m



9 Add 2 mL autoinduction medium to each well in 4 sterile 24-well plates.

5m

10 Add 50 μ L of the transformation mixture to each well in the 24-well plates. Resuspend any pelleted bacteria before transferring.

15m




Note

To keep the layout compatible with a multichannel pipette, this template can be used (assuming samples 1-96)

96 well plate:

	1	2	3	4	5
A	1	2	3	4	5
B	13	14	15	16	17
C	25	26	27	28	29
D	37	38	39	40	41
E	49	50	51	52	53
F	61	62	63	64	65
G	73	74	75	76	77
H	85	86	87	88	89
	6	7	8	9	10
A	6	7	8	9	10
					



D

47

48

E

59

60

F

71

72

G

83

84

H

95

96

24-well plate 1:

	1	2	3	4	5
A	1	2	3	4	5
B	13	14	15	16	17
C	25	26	27	28	29
D	37	38	39	40	41

	6
A	6
B	18
C	30
D	42

24-well plate 2:

	1	2	3	4	5
A	7	8	9	10	11
B	19	20	21	22	23
C	31	32	33	34	35
D	43	44	45	46	47


	6
A	12
B	24
C	36
D	48

24-well plate 3:

	1	2	3	4	5
A	49	50	51	52	53
B	61	62	63	64	65
C	73	74	75	76	77
D	85	86	87	88	89

	6
A	54
B	66
C	78
D	90

24-well plate 4:



11 Cover plates with breathable film.



- 12 Incubate cultures  400 rpm, 37°C  Overnight .

12h



Harvest (day 2)

47m

- 13 Add 90 μ L LB per well in a 96-well plate.

2m

Note

Should be plate reader compatible.



- 14 Add 10 μ L of culture per well.

5m

Note

Make sure to resuspend if pelleted.



- 15 Measure OD600 in the plate reader. Remember to include a couple of wells with just LB as blank controls.

10m

Note

OD of the culture will then be 10x this number. This step is qualitative and meant as a quick check of the growth. The values vary between experiments.



- 16 Take a 15 μ L sample of each production from diluted 96-well plate samples and save for SDS-PAGE analysis.

10m

- 17 Pellet bacteria in the 24-well plate  2500 x g, Room temperature, 00:15:00

15m

- 18 Remove the supernatant. Freeze 24-well plate or proceed to lysis.

5m





- 19 Decide on your lysis protocol. This will depend on the stability of your proteins of interest.

STEP CASE

If you have heat stable proteins

16 steps

This protocol will remove many of the endogenous *E. coli* protein in your samples⁶.

- 20 Preheat Thermomixer  95 °C 15m
- 21 Resuspend the pellets in 500 µL lysis buffer and move to a heat-stable 96-well plate. 10m
- #### Note
- If the lysate is very viscous you can increase the volume, but this will decrease the final protein concentration.
- 22 Heat-lyse the bacteria  95 °C  00:15:00 . 15m
- 23 Optional to ease filtration: Centrifuge  3000 x g, Room temperature, 00:10:00 and take only the supernatant to the next step. 10m
- 24 Transfer the lysate to a 96-well filter plate stacked on a collection plate (96 deep-well plate). 5m
- 25 Centrifuge  1500 x g, Room temperature, 00:05:00 . If the liquid in some wells has not gone fully through, resuspend with pipette and centrifuge again. 5m
- 26 Take 15 µl samples from cleared lysates for SDS-PAGE. 15m
- 27 Estimate protein concentration in the lysates using the Bicinchoninic Acid (BCA) Method. *



Affinity chromatography purification (day 2, optional)

1h 10m

- 28 Prepare 200 mL binding buffer (20 mM imidazole) and 60 mL elution buffer (500 mM imidazole). Open the plate according to instructions and place it on top of a collection plate (96 deep-well). 15m
- 29 Centrifuge the plate to remove the storage solution from the resin 5m
 500 x g, 4°C, 00:02:00
- 30 Add 500 µL deionized water/well  500 x g, 4°C, 00:02:00 . 5m
- 31 Add 500 µL binding buffer/well and mix briefly to equilibrate the resin, centrifuge  500 rpm, 4°C, 00:02:00 . Repeat once. 10m
- 32 Apply filtered lysate to the wells, mix briefly, and incubate  00:03:00 5m
- 33 Remove the flow-through by centrifuging 5m
 100 x g, 4°C, 00:04:00 , or until the wells are empty
- 34 Add 500 µL binding buffer/well and mix briefly to wash out unbound sample, centrifuge  500 x g, 4°C, 00:02:00 10m
Repeat once.
- 35 Add 200 µL of elution buffer/well and mix for  00:01:00 . Change collection plate and centrifuge the plates  500 x g, 4°C, 00:02:00 and collect the fractions. 15m
Repeat twice.

Note

If required, change collection plate between each elution (to prevent unnecessary dilution of the target protein).

Protocol references

1. Watson JL, Juergens D, Bennett NR, Trippe BL, Yim J, Eisenach HE, Ahern W, Borst AJ, Ragotte RJ, Milles LF, Wicky BIM, Hanikel N, Pellock SJ, Courbet A, Sheffler W, Wang J, Venkatesh P, Sappington I, Torres SV, Lauko A, De Bortoli V, Mathieu E, Ovchinnikov S, Barzilay R, Jaakkola TS, DiMaio F, Baek M, Baker D. De novo design of protein structure and function with RFdiffusion. <https://doi.org/10.1038/s41586-023-06415-8>
2. Dauparas J, Anishchenko I, Bennett N, Bai H, Ragotte RJ, Milles LF, Wicky BIM, Courbet A, de Haas RJ, Bethel N, Leung PJY, Huddy TF, Pellock S, Tischer D, Chan F, Koepnick B, Nguyen H, Kang A, Sankaran B, Bera AK, King NP, Baker D. Robust deep learning-based protein sequence design using ProteinMPNN. <https://doi.org/10.1126/science.add2187>
3. Pacesa M, Nickel L, Schellhaas C, Schmidt J, Pyatova E, Kissling L, Barendse P, Choudhury J, Kapoor S, Alcaraz-Serna A, Cho Y, Ghamary KH, Vinué L, Yachnin BJ, Wollacott AM, Buckley S, Westphal AH, Lindhoud S, Georgeon S, Goverde CA, Hatzopoulos GN, Gönczy P, Muller YD, Schwank G, Swarts DC, Vecchio AJ, Schneider BL, Ovchinnikov S, Correia BE. BindCraft: one-shot design of functional protein binders. <https://doi.org/pii:2024.09.30.615802.10.1101/2024.09.30.615802>
4. Rosenberg, A. H. et al.. Vectors for selective expression of cloned DNAs by T7 RNA polymerase. [https://doi.org/10.1016/0378-1119\(87\)90165-X](https://doi.org/10.1016/0378-1119(87)90165-X)
5. Studier, F William. Protein production by auto-induction in high density shaking cultures. <https://doi.org/10.1016/j.jpep.2005.01.016>
6. Christoph Kalthoff. A novel strategy for the purification of recombinantly expressed unstructured protein domains. [https://doi.org/10.1016/S1570-0232\(02\)00908-X](https://doi.org/10.1016/S1570-0232(02)00908-X)

Citations

Francisca Pinheiro, Jan S. Nowak, Elena Zueva, Emily C. Pheasant, Ida Kjærsgaard Grene, Vili Lampinen, Magnus Kjaergaard. Screening de novo designed protein binders in unpurified lysate using flow induced dispersion analysis
<https://doi.org/10.1101/2025.06.17.6601270>

Acknowledgements

We thank Nadja Joachim and Kaare Teilum for sharing their protein expression protocol.