



May 22, 2024

## Melatonin ELISA

DOI

[dx.doi.org/10.17504/protocols.io.81wgbx35ylpk/v1](https://dx.doi.org/10.17504/protocols.io.81wgbx35ylpk/v1)

daniel.dautan daniel<sup>1,2</sup>, Per Svenningsson<sup>1,2</sup>

<sup>1</sup>Department of Clinical Neuroscience, Karolinska Institutet, 171 76 Stockholm, Sweden;

<sup>2</sup>Aligning Science Across Parkinson's (ASAP) Collaborative Research Network, Chevy Chase, MD 20815, USA



Jacquelyn Haytayan

Weill Cornell Medicine

### Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.81wgbx35ylpk/v1>

**Protocol Citation:** daniel.dautan daniel, Per Svenningsson 2024. Melatonin ELISA. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.81wgbx35ylpk/v1>

**License:** This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

We use this protocol and it's working



**Created:** February 15, 2024

**Last Modified:** September 23, 2024

**Protocol Integer ID:** 95309

**Keywords:** ASAPCRN, melatonin, sleep, melatonin elisa measurement of mouse plasma melatonin, melatonin elisa measurement, mouse plasma melatonin, using elisa kit, elisa kit, mouse plasma

**Funders Acknowledgements:**

**Aligning Science Across Parkinson's**

Grant ID: 020608

## Abstract

**Measurement of mouse plasma melatonin** using ELISA Kit (Enzo Life Sciences, ENZ-KIT150-0001, NY, US) according to manufacturer instructions.








## Materials

Enzo Life Sciences, ENZ-KIT150-0001

## Troubleshooting



## Melatonin Extraction










- 1 Mix  200  $\mu\text{L}$  of plasma with an equal volume of cold Ethyl Acetate. Vortex gently.
- 2 Allow layers to separate on ice for  00:03:00 . Vortex again and incubated  On ice for  00:02:00 .
- 3 Spin samples at 1000g for  00:10:00 at  4  $^{\circ}\text{C}$  . Transfer the organic layer to a new tube.
- 4 Dry samples and resuspend in  220  $\mu\text{L}$  of 1X stabilizer.

5m

10m

## ELISA





2h

- 5 Add  100  $\mu\text{L}$  of standards' working solutions and samples to provided 96-well plate in duplicates.
- 6 Immediately after add  50  $\mu\text{L}$  of melatonin tracer to each well (except blanks) followed by  50  $\mu\text{L}$  of 1X melatonin antibody (except blanks).
- 7 Cover plates with the provided plate sealer. Incubate for  01:00:00 at  37  $^{\circ}\text{C}$  with 500 rpm shaking.
- 8 Decant the solution from each well. Add  400  $\mu\text{L}$  of wash solution to each well.
- 9 Decant the solution from each well and pat dry against clean absorbent paper.
- 10 Repeat wash step 3 times.
- 11 Add  200  $\mu\text{L}$  of melatonin conjugate solution to each well (except blanks). Cover plate with the sealer. Incubate for  00:30:00 at  Room temperature .

1h

30m



- 12 Perform the wash step as described above. Add  200  $\mu\text{L}$  of substrate reagent to each well. Incubate for  00:30:00 at  37 °C protected from light. 30m
- 13 Add  50  $\mu\text{L}$  of stop solution to each well. Measure the optical density using a microplate reader with absorbance set to 450 nm.