

Apr 28, 2024

Version 4

MERS-CoV Mpro fluorescence dose response V.4

DOI

<https://dx.doi.org/10.17504/protocols.io.eq2ly7r1rlx9/v4>

Haim Barr^{1,2}, Noa Lahav^{1,2}

¹The Weizmann Institute of Science; ²ASAP Discovery

Haim Barr: General acknowledgement: The Wohl Drug Discovery institute, The Nancy and Stephen Grand Israel National Center for Personalized Medicine.;

ASAP Discovery



Noa Lahav

Weizmann Institute

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.eq2ly7r1rlx9/v4>

Protocol Citation: Haim Barr, Noa Lahav 2024. MERS-CoV Mpro fluorescence dose response. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.eq2ly7r1rlx9/v4> Version created by [Mary-Ann Xavier](#)

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: April 12, 2024

Last Modified: October 07, 2024

Protocol Integer ID: 98133

Keywords: mpro fluorescence dose response, biochemical assay, direct enzyme activity measurement method, treatments for viral infectious disease, viral infectious disease, fluorescence intensity, final experiment concentration, cov, like protease, mer, final assay

Funders Acknowledgements:

National Institutes of Health/National Institute Of Allergy and Infectious Diseases (NIH/NIAID)

Grant ID: U19AI171399

Disclaimer

The content is solely the responsibility of the authors and does not necessarily represent the official views of the National Institutes of Health.

Abstract

This is a functional, biochemical assay used to identify treatments for viral infectious disease in MERS-CoV 3C-like protease.

Utilizing a direct enzyme activity measurement method, the experiment was performed in a 384-well plate reading the fluorescence intensity. This assay tested the mode of action of inhibition.

Final Experiment Concentrations

	A	B	C
	Reagent	Concentration	Units
	MERS Mpro	50	nM
	MERS Substrate peptide	550	nM
	HEPES pH=7.3	20	mM
	NaCl	50	mM
	BSA	0.1	mg/ml
	Triton X-100	0.01	% (v/v)
	TCEP	1	mM

For more information, please see the Materials section.

Guidelines

Plate Information:

Total Assay Volume: 20 μ L

Compounds Top Assay Concentration: 100 μ M

Dilution Factor: 2

Dose Response Points: 12

Number of Replicates: 2

Backfill with DMSO: Yes

Materials

Assay Buffer Reagents (Concentration listed are the stock concentrations)

- [M] 40 millimolar (mM)  HEPES Buffer (pH 7.3) Fisher Scientific Catalog #BP299-1 (or similar)
- [M] 100 millimolar (mM)  Sodium Chloride Fisher Scientific Catalog #S271 (or similar)
- [M] 10 mg/mL  BSA-Molecular Biology Grade - 12 mg New England Biolabs Catalog #B9000S (or similar)
- [M] 10 % volume  Triton X-100 Merck MilliporeSigma (Sigma-Aldrich) Catalog #T8787-50ML (or similar)
- [M] 1000 millimolar (mM)  TCEP HCl P212121 Catalog #SV-TCEP (or similar)

(all components are added fresh to the assay buffer before each experiment)

Additional Reagents

- [M] 507000 nanomolar (nM) **MERS Mpro Enzyme***

***Note:** The original MERS Mpro stock enzyme had a concentration of [M] 507000 nanomolar (nM) but when a new stock solution was made/delivered on 2023-02-21 the new stock had a concentration of

[M] 478000 nanomolar (nM) Both stock solutions were diluted with fresh assay buffer to create a

[M] 100 nanomolar (nM) solution before each experiment

- [M] 750000 nanomolar (nM) **MERS Substrate***

***Note:** MERS Substrate (5-FAM)-GVLQSGLV-K(Dabcyl)-K-NH₂ Stock was purchased from Peptide 2.0 and dissolved in DMSO with an original concentration of [M] 750000 nanomolar (nM) however it was diluted with the same assay buffer used in the experiment before conducting each experiment to yield a concentration of

[M] 1100 nanomolar (nM)

Safety warnings

 Please be sure to wear proper Personal Protective Equipment (PPE) while performing this experiment.

Before start

Note: Inhibitor compounds stock concentration is **20 mM**. Compounds are pre-dispensed into 384 plates and stored at -20 C until use.

Prepare Reagents

1 **PREPARE** all of the reagents/buffers required for this experiment.

A	B	C	D	E
Reagent	Stock Concentration	Concentration Loaded into Combi	Final Concentration	Units
MERS-Cov Mpro Enzyme (original stock)	507000	100	50	nM
20230221 MERS Mpro Enzyme	478000	100	50	nM
MERS Substrate	750000	1100	550	nM
Assay Buffer				
HEPES (pH 7.3)	40	20	20	mM
Sodium Chloride	100	50	50	mM
BSA	10	0.1	0.1	mg/mL
Triton X-100	10	0.01	0.01	% by volume
TCEP	1000	1	1	mM

For more information, please see the Materials section

Prepare 384 Well Plate

15m

2 **PRIME** with **Assay Buffer** by Multi-Drop Combi Tube Dispensing Cassette by selecting the **PRIME** button on the Combi Dispenser until the tubes are filled completely.

2.1 **DISPENSE**  10 µL Assay Buffer to Columns **1** and **23** of assay plate

- **Note:** These will represent the ***inhibitor control columns*** (Contain: Substrate, Assay Buffer, DMSO, **no experimental compounds**)



- 2.2 **EMPTY** Multi-Drop Combi Tube Dispensing Cassette (by selecting the **EMPTY** button on the Combi Dispenser until the tubes of the cassette are emptied).
- Discard Assay Buffer discharged from the cassette.
- 3 **PRIME** with [M] 100 nanomolar (nM) MERS MPro Enzyme by Multi-Drop Combi Tube Dispensing Cassette by selecting the **PRIME** button on the Combi Dispenser until the tubes were filled completely.
- **Note:** Be sure to cycle dispensing several times on a clean plate lid (This confirms there are no bubbles in the Dispensing Cassette).
- 3.1 **DISPENSE** 10 µL [M] 100 nanomolar (nM) MERS MPro Enzyme to Columns **2** through **22** and Column **24**
- Note:**
- [M] 100 nanomolar (nM) MERS MPro Enzyme is two times the final concentration for the assay. It is diluted to be a final concentration of [M] 50 nanomolar (nM) MERS MPro Enzyme .
 - Column 2 and Column 24 are **neutral control columns** (Contain: Enzyme, Substrate, DMSO, no compounds)
- 3.2 **EMPTY** Multi-Drop Combi Tube Dispensing Cassette (by selecting the **EMPTY** button on the Combi Dispenser until the tubes of the cassette are emptied.)
- Discard the [M] 100 nanomolar (nM) MERS MPro enzyme discharged from the cassette.
- 4 **CENTRIFUGE** 15000 rpm, Room temperature, 00:01:00 plate to remove bubbles 1m
- 5 **INCUBATE** plate for 00:15:00 at Room temperature 15m
- ▲ Make sure the plate is protected from light!
- 6 **PRIME** with **Assay Buffer** by Multi-Drop Combi Tube Dispensing Cassette by selecting the **PRIME** button on the Combi Dispenser until the tubes were filled completely. Then, **EMPTY** the Multi-Drop Combi Tube Dispensing Cassette.
- 7 **PRIME** with [M] 1100 nanomolar (nM) MERS Substrate by Multi-Drop Combi Tube Dispensing Cassette by selecting the **PRIME** button on the Combi Dispenser until the tubes were filled completely.



- **Note:** Be sure to cycle dispensing several times on a clean plate lid (This confirms there are no bubbles in the Dispensing Cassette).

8 **DISPENSE**  10 µL 1100nM MERS Substrate into Columns **1 through 24** (the full plate)

Note:

- [M] 1100 nanomolar (nM) MERS Substrate is two times the final concentration for the assay. It is diluted to be a final concentration of [M] 50 nanomolar (nM) MERS Substrate

9 **CENTRIFUGE**  15000 rpm, Room temperature, 00:01:00 plate in plate centrifuge to remove bubbles

1m

10 **INCUBATE** plate for  01:00:00 at  Room temperature

1h

▲ Make sure the plate is protected from light!

Recommended: Clean the Multi-Drop Combi Reagent Dispenser during this incubation step

Read Plate Fluorescence

11 **READ** and **RECORD** the plate Relative fluorescence units (RFU) via the "**MERS Protocol**" on the **PHERASTAR FS Control Software**.

- Software is a standard Fluorescence Assay set for Optimal excitation wavelength 485 nm, emission wavelength 528 nm, and a Gain of 300.



Equipment

PHERASTAR FS

Microplate reader

BMG LABTECH

0471B0001A

https://www.bmglabtech.com/en/pherastar-fsx/?utm_term=pherastar%20plate%20reader&utm_campaign=usa.roi.products&utm_source=FwOoFR_5EUHUaAIkREALw_wcB

Expected result

Gain 300 should yield ~10,000 RFU in full reaction and ~6,000 RFU in Buffer Control