

Mar 25, 2020

## Measuring JF dye kinetics in the brain of C57/Bl mice

 In 2 collections

DOI

[dx.doi.org/10.17504/protocols.io.5vmg646](https://dx.doi.org/10.17504/protocols.io.5vmg646)

Boaz Mohar<sup>1</sup>

<sup>1</sup>HHMI Janelia Research Campus



Boaz Mohar

HHMI Janelia Research Campus

### Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.5vmg646>

**Protocol Citation:** Boaz Mohar 2020. Measuring JF dye kinetics in the brain of C57/Bl mice. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.5vmg646>

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

We use this protocol and it's working



**Created:** July 26, 2019

**Last Modified:** March 25, 2020

**Protocol Integer ID:** 26253

**Keywords:** jf dye kinetics in the brain, measuring jf dye kinetic, systemic dye injection, jf dye, clearance kinetics of jf dye, retro orbital injection, clearance of the dye, dye, bl mice the purpose, orbital injection, bl mice, mice, mouse, double exponential decay function, systemic delivery

## Abstract

The purpose of this protocol is to compare the intake and clearance kinetics of JF dye systemic delivery using retro orbital injections.

We image a mouse with a cranial window before and after systemic dye injection.

We want to fit the clearance of the dye with a double exponential decay function:

$$y = a * e^{-\frac{1}{b*x}} + c * e^{-\frac{1}{d*x}} + e$$

## Troubleshooting



## Baseline

- 1 Use a mouse with a window (doesn't have to be in ALM)
- 2 Move the mouse to an induction chamber with 3% ISO and flow of ~2L.min. Wait until the mouse takes ~ 1 breaths per second
- 3 Move the mouse to the microscope and clean the window with Q-tips and 70% ethnol. 3m
- 4 Select if you want to use a 4x or 20x objective.  
If using a 20x add water above the window.  
Lower the objective using room light or white light illumination and focus on the top of the cortex.
- 5 Pick the appropriate filter cube and illumination source for the dye
- 6 Setup the imaging conditions in a way that baseline image takes ~ 10-20% of the dynamic range but is above read noise level. We use 50ms exposure with 40 frames per timepoint totaling ~2s duration per timepoint. Or 200ms with 10 frames if signal is lower.
- 7 Record baseline image at a few illumination settings to make sure that the first timepoints after injection are not saturated.  
For example with a LED use 100mA, 500mA 1A and 2A driving current for the same acquisition settings.  
When imaging the first image after injection try using 500mA but change if it is too dark / saturated.
- 8 Record a white light reference image so you don't have to bleach your signal to return to the same field of view. If possible record the objective's location.

## Dye injection and clearance

- 9 Take the animal out of the imaging rig and back to the induction box if it is breathing faster than 1/s
- 10 Prepare a fresh dye aliquot

- 11 Inject into the **retro-orbital** sinus
- 12 Return the animal to the rig and return to the same field of view using the white light image as reference
- 13 Start imaging using the fluorescence light source at fast interval at first (~1-3min) to catch the rise and fast decay phases and depending on the dynamics start reducing the intervals (5-20min) after ~ 1h you should get a sense of the amount of change you see and start separating time points by more than 30min.
- 14 After ~4h you can usually take the animal out of the rig and perform recordings at >2h intervals.
- 15 After 24h usually no more dynamics is observed.

## Analysis

- 16 Load images and mean each timepoint, save the mean image
- 17 Load the mean images in **Fiji** and using the stack sorter move the first (baseline) image to be the last.
- 18 Use the Linear Stack Alignment with SIFT plugin to align the stack of images.
- 19 Save the aligned stack and load in Matlab
- 20 Define a mask (ellipse or otherwise) of the imaged region of interest (imellipse). Median the pixels there, subtract the baseline timepoint
- 21 Fit a double exponent model to the rest of the median pixels (excluding baseline) using the file creation date as x axis.  
$$y = \mathbf{a} * e^{-\frac{1}{\mathbf{b} * x}} + \mathbf{c} * e^{-\frac{1}{\mathbf{d} * x}}$$



Add a lower bound of 0 to all parameters.

If cell are expressing a HaloTag, use:

$$y = \mathbf{a} * e^{-\frac{1}{\mathbf{b} * x}} + \mathbf{c} * e^{-\frac{1}{\mathbf{d} * x}} + \mathbf{e}$$

First point used will be the peak (Most times the first as the rise is very fast)