

Jun 01, 2025

Version 2

🌐 Metabarcoding for Fish and Crustaceans in Diet Samples Using 2-PCR protocol with Unique Dual Indexing V.2

🔗 Version 1 is forked from [Metabarcoding Fecal Swabs or Stomach Contents for Fish and Crustaceans using 2-PCR protocol and Illumina MiSeq](#)



DOI

<https://dx.doi.org/10.17504/protocols.io.261ge5dnwg47/v2>

Eldridge Wisely¹

¹University of Arizona



Eldridge Wisely

University of California, San Diego

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.261ge5dnwg47/v2>

Protocol Citation: Eldridge Wisely 2025. Metabarcoding for Fish and Crustaceans in Diet Samples Using 2-PCR protocol with Unique Dual Indexing. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.261ge5dnwg47/v2> Version created by **Eldridge Wisely**

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: May 23, 2024

Last Modified: June 01, 2025

Protocol Integer ID: 100414

Keywords: rna gene of crustacean, pcr protocol with unique dual indexing, crustaceans in diet sample, including unique molecular index, mitochondrial 12s rna gene, mitochondrial16, unique molecular index, rna gene, indexing primer, 170bp region of the mitochondrial16, second pcr step, pcr protocol, pcr, resulting pcr product, fish, unique dual indexing, pcr product, unique dual index, metabarcoding, removal of pcr, illumina sequencer, indexing, diet sample, addition to the unique dual index

Funders Acknowledgements:

National Science Foundation

Grant ID: OCE #2128592

Abstract

This protocol describes a method to metabarcode a 170bp region of the mitochondrial16S rRNA gene of crustaceans and a 163-185bp region of the mitochondrial 12S rRNA gene of fishes. These regions are subjected to PCR separately in multiple replicates and the resulting PCR products are pooled by sample and then indexed with Nextera-style unique dual indexes which are compatible with all Illumina sequencers including NovaSeq. This protocol differs from the protocol it was forked from only in the indexing primers which save time in the lab, reduce the potential for human-error in indexing, and facilitate the removal of PCR-errors introduced in the second PCR step by including Unique Molecular Indexes (UMIs) in addition to the Unique Dual Indexes to eliminate tag-jumps.

Image Attribution

Haley Capone

Guidelines

The PCR conditions described here are different from the PCR conditions described by Miya et al., and Berry et al. in their respective publications introducing the primers used here. This difference is due to the use of the Takara High Fidelity PCR EcoDry Premix in this protocol.

Materials

96-well PCR plates
Adhesive foil PCR plate covers

1.5mL tubes

Illumina® DNA/RNA UD Indexes Set A, Tagmentation (96 Indexes, 96 Samples)
20091654

Optional additional Indexes for more samples to be sequenced in the same sequencing run:

Illumina® DNA/RNA UD Indexes Set B, Tagmentation (96 Indexes, 96 Samples)
20091656

Illumina® DNA/RNA UD Indexes Set C, Tagmentation (96 Indexes, 96 Samples)
20091658

Illumina® DNA/RNA UD Indexes Set D, Tagmentation (96 Indexes, 96 Samples)
20091660

PCR machine

Equipment to run gels
optionally: equipment for fluorometric quantification

Equipment

96-well Magnetic Rack Separator

Magnetic Rack Separator

Sergi Lab Supplies

B08134P9RT

https://www.amazon.com/Magnetic-Separator-Protein-Purification-Format/dp/B08134P9RT/ref=asc_df_B08134P9RT/?tag=&linkCode=df0&hvadid=416872221972&hvpos=&hvnetw=g&hvrnd=12953200023550024012&hvpon e=&hvptwo=&hvqmt=&hvdev=c&hvdvcm dl=&hvlocint=&hvlocphy=903024

NAME

TYPE

BRAND

SKU

L
I
N
K

Equipment

Magnetic Rack for for 1.5 mL Tubes

NAME

Magnetic Rack for DNA, RNA Purification; for 1.5 mL centrifuge Tubes

TYPE

Sergi Lab Supplies

BRAND

B0BZWXZM22

SKU

https://www.amazon.com/Magnetic-Rack-Purification-centrifuge-Tubes/dp/B0BZWXZM22/ref=asc_df_B0BZWXZM22/?tag=hyprod-20&linkCode=df0&hvadid=652498086131&hvpos=&hvnetw=g&hvrand=6716034042841103246&hvpone=&hvpstwo=&hvmqmt=&hvddev=c&hvdvcmdl=&hvlocint=&hvlocphy=9

L
I
N
K

Protocol materials

⊗ MiFish-F Integrated DNA Technologies, Inc. (IDT) Catalog #custom

⊗ MiFish-R Integrated DNA Technologies, Inc. (IDT) Catalog #custom

⊗ Crustacean16S-F Integrated DNA Technologies, Inc. (IDT) Catalog #custom

⊗ Crustacean16S-R Integrated DNA Technologies, Inc. (IDT) Catalog #custom

⊗ Agencourt AMPure XP Beckman Coulter Catalog #A63880

⊗ Buffer EB Qiagen Catalog #19086

⊗ 2x Kapa HiFi Hotstart Readymix Kapa Biosystems Catalog #KK2602

⊗ Nuclease-free water Integrated DNA Technologies, Inc. (IDT) Catalog #11-05-01-14

Troubleshooting

Before start

Work in a pre-PCR lab, as separated as possible from post-PCR products.

Clean work area with 10% bleach solution before beginning work for the day, then change gloves so that no bleach carryover to your samples or reactions occurs.



Prepare Primers

- 1 Order metabarcoding primers with diversity spacers and Illumina overhang sequences (Illumina, 2013):

 MiFish-F Integrated DNA Technologies, Inc. (IDT) Catalog #custom (Miya et al., 2015):

TCGTCGGCAGCGTCAGATGTGTATAAGAGACAGNNNGTCGGTAAACTCGTGCCAGC

 MiFish-R Integrated DNA Technologies, Inc. (IDT) Catalog #custom (Miya et al., 2015):

GTCTCGTGGGCTCGGAGATGTGTATAAGAGACAGNNNNCATAGTGGGGTATCTAATCCCAG
TTTG

 Crustacean16S-F Integrated DNA Technologies, Inc. (IDT) Catalog #custom Berry et al., 2017):

TCGTCGGCAGCGTCAGATGTGTATAAGAGACAGNNNGGACGATAAGACCCTATA

 Crustacean16S-R Integrated DNA Technologies, Inc. (IDT) Catalog #custom

(Berry et al., 2017):

GTCTCGTGGGCTCGGAGATGTGTATAAGAGACAGNNNATTACGCTGTTATCCCTAAAG

We got ours from <https://www.idtdna.com/> as 4 nmole Ultramer DNA Oligos, with standard desalting.

- 2 Briefly centrifuge primer tubes, then reconstitute primers to [M] 100 micromolar (μM) stock solutions by adding 40uL of

 Nuclease-free water Integrated DNA Technologies, Inc. (IDT) Catalog #11-05-01-14

to each tube of 4nmole primers.

Incubate at  Room temperature for a minimum of  00:05:00 , then vortex gently and centrifuge briefly.

5m




- 3 Make [M] 5 micromolar (μM) working solutions of each primer by adding  95 μL of
-  Nuclease-free water **Integrated DNA Technologies, Inc.**
(IDT) **Catalog #11-05-01-14**
- and  5 μL of primer stock solution for each  100 μL of primer that you intend to use within the next week or so.

Create Plate Map

- 4 Determine which sample will go into each well. This should be the same for each primer set and each replicate. Include at least one extraction control (you can combine aliquots of the extraction controls from each round of DNA extraction into one tube, and use that as your single extraction control), and include a PCR negative control for each plate of PCR. See example below of 21 samples, a field negative sample, a combined extraction control, and a PCR negative.

	1	2	3	4	5
A	Sample8_rep1	Sample16_rep1	Sample24_rep1	PCR neg_rep1	Sample8_rep2
B	Sample7_rep1	Sample15_rep1	Sample23_rep1	Ext Con_rep1	Sample7_rep2
C	Sample6_rep1	Sample14_rep1	Sample22_rep1	Field Neg_rep1	Sample6_rep2
D	Sample5_rep1	Sample13_rep1	Sample21_rep1	Sample29_rep1	Sample5_rep2
E	Sample4_rep1	Sample12_rep1	Sample20_rep1	Sample28_rep1	Sample4_rep2
F	Sample3_rep1	Sample11_rep1	Sample19_rep1	Sample27_rep1	Sample3_rep2
G	Sample2_rep1	Sample10_rep1	Sample18_rep1	Sample26_rep1	Sample2_rep2
H	Sample1_rep1	Sample9_rep1	Sample17_rep1	Sample25_rep1	Sample1_rep2
	6	7	8	9	10
A			PCR neg_rep2	Sample8_rep3	



- 4.1 Do not mix sample types between invasively sampled methods (fecal swabs, or stomach contents) and non-invasively sampled methods (eDNA from water or sediment) in the same PCR procedure.



MiFish Takara PCR Conditions

- 5 Make your MiFish Mastermix:

For each **PCR replicate of each sample** you intend to process (+10% overage), mix:

 0.7 μ L  5 micromolar (μ M)

 MiFish-F Integrated DNA Technologies, Inc. (IDT) Catalog #custom

 0.7 μ L  5 micromolar (μ M)

 MiFish-R Integrated DNA Technologies, Inc. (IDT) Catalog #custom

 22.6 μ L

 Nuclease-free water Integrated DNA Technologies, Inc. (IDT) Catalog #11-05-01-14

For a full plate of 96 reactions, multiply 105.6*the per-sample volumes in the recipe to make the mastermix.

- 5.1 Add  24 μ L of your MiFish metabarcoding mastermix to each well of

 Takara High Fidelity PCR EcoDry Premix Takara Bio Inc. Catalog #639280

- 5.2 Add  1 μ L DNA extracted from diet samples.

- 5.3 Mix and stir together with pipette tip, swirling to make sure the liquid is in the bottom, and bringing any bubbles to the surface of each reaction.

- 5.4 Cap each row of reaction tightly before beginning any other PCR reaction in the same room.

MiFish Takara PCR Conditions

- 6  95 $^{\circ}$ C for  00:01:00

35 cycles of:

3m 30s

🌡️ 95 °C for ⌚ 00:00:30

🌡️ 66 °C for ⌚ 00:01:00

followed by:

🌡️ 68 °C for ⌚ 00:01:00

Hold at 🌡️ 4 °C

Crustacean_16S Takara PCR Recipe

7 Make your Crustacean_16S Mastermix:

For each **PCR replicate of each sample** you intend to process (+10% overage), mix:

🧪 2 µL [M] 5 micromolar (µM)

🧬 Crustacean16S-F Integrated DNA Technologies, Inc. (IDT) Catalog #custom

🧪 2 µL [M] 5 micromolar (µM)

🧬 Crustacean16S-R Integrated DNA Technologies, Inc. (IDT) Catalog #custom

🧪 20 µL

🧬 Nuclease-free water Integrated DNA Technologies, Inc. (IDT) Catalog #11-05-01-14

For a full plate of 96 reactions, multiply 105.6*the per-sample volumes in the recipe to make the mastermix.

8 Add 🧪 24 µL of your Crustacean_16S metabarcoding mastermix to each well of

🧬 Takara High Fidelity PCR EcoDry Premix Takara Bio Inc. Catalog #639280

8.1 Add 🧪 1 µL DNA extract

8.2 Mix and stir together with pipette tip, swirling to make sure the liquid is in the bottom, and bringing any bubbles to the surface of each reaction.

- 8.3 Cap each row of reaction tightly before beginning any other PCR reaction in the same room.

Crustacean_16S Takara PCR Conditions

4m

- 9  95 °C for  00:01:00
- 35 cycles of:
-  95 °C for  00:00:30
-  50 °C for  00:01:00
-  68 °C for  00:00:30
- followed by:
-  68 °C for  00:01:00
- then hold at  4 °C

4m

Combine PCR Products of Biological Samples by Sample and Negatives by Primer

- 10 Gently vortex to mix and lightly centrifuge PCR products to get any bubbles from the bottom and any droplets off of the lids.
8-strip tubes can be put into a plate holder and spun down as a plate if this option is available.
- 11 Get a new sterile 96-well plate out of its packaging and immediately cover with foil (a 15-minute treatment under a UV light is helpful to sterilize before covering with foil).
- 11.1 For biological samples (not PCR negative controls): combine  10 µL of each of the 6 PCR products **by sample** into the new sterilized plate.
- Peel back one row of foil at a time and using a multtip pipette, open only one row of PCR products at a time to combine. Cap the resulting combined row of the new plate with a strip-cap and close the PCR products of the previous row before opening a new row.
- This should result in 60uL of an equal volume of MiFish and Crustacean-16S PCR products, each from the same original sample.
- 11.2 For negative controls: combine  20 µL of each of the 3 negative PCR controls **by primer**. You should have one MiFish combined negative and one Crustacean-16S negative, each with a total of  60 µL

11



Visualize PCR Products

- 12 Make a 1.7% to 2% agarose gel and run a representative sample of reactions on it to make sure the PCRs worked, producing bands in the 250-300bp range. Use a ladder than allows you to distinguish small bands between 100-500bp. Check some PCR negatives to see that they don't have bands. Be very careful opening the PCR plate wells at this point to avoid cross-contamination.
- 12.1 Run the gel at 100V until the dye band has traveled $\frac{3}{4}$ of the length of the gel, then visualize.
- 12.2 Any DNA smaller than 250 bp will be excess primer and oligos and will need to be cleaned with SPRI beads. Check that the initial PCR worked before continuing on with the library-building procedure. II

Prepare EtOH for bead cleanup, and bring beads to room temperature

12m 30s

- 13 Get AmpureXP beads out of the refrigerator, and bring to room temp, swirl to mix occasionally, or use a rocking platform.
- 14 Make fresh 80% EtOH so that you will have at least  200 μL of EtOH per well of the combined plate.
- 15 Get 2 sterile DNAase/RNAse free 96-well PCR plates out of their packaging and immediately cover with adhesive foil or if possible, UV clean the plates for  00:15:00, then immediately cover with adhesive foil. 15m

One plate will be for the bead-cleanup steps, and the other will be for the final, cleaned reactions.

Perform a 1.5x bead cleanup with Ampure XP beads.

12m 30s

- 16 In the bead-cleanup plate, do the following steps for one 8-sample row of the plate at a time, pulling back the foil cover for each row after the previous one has been completed.
- 16.1 Add 1.5x the sample volume of Ampure XP beads. 5m

In this case you have 60uL of combined PCR products per well, so you will add

 90 μL of room temperature, well-mixed

 Agencourt AMPure XP **Beckman Coulter Catalog #A63880** to each well, and pipette mix very thoroughly, by stirring and pipetting up and down ten times.

Incubate  00:05:00 at room temperature.

17 After incubation, place 96-well plate on a

7m

Equipment

96-well Magnetic Rack Separator

Magnetic Rack Separator

Sergi Lab Supplies

B08134P9RT

<https://www.amazon.com/Magnetic-Separator-Protein-Purification-Format/dp/B08134P9RT?tag=&linkCode=df0&hvadid=416872221972&hvpos=&hvnetw=g&hvrnd=12953200023>

for  00:02:00 to  00:05:00 or until liquid is clear.

18 remove and discard liquid from the row, being careful not to touch the beads with the pipette or to let the beads dry for more than 30 seconds.

Tip: If you do get beads in your pipette tip, just put the liquid and beads back into the well and wait until the solution clears before trying again.



18.1 Add  100 μL of fresh 80%EtOH to each well of beads, without disturbing the beads or removing the plate from the magnet. Incubate at  Room temperature for

30s

 00:00:30

18.2 Remove the EtOH, then immediately add another  100 μL of 80% EtOH to the wells, incubate for  00:00:30  Room temperature .

30s

18.3 Remove ALL EtOH, and let the row of beads dry just enough to lose some shine but not enough to start cracking. This should be approximately  00:00:30 to  00:01:00 1m 30s

18.4 Remove the plate with new row of cleaned beads from the magnetic plate, and add  30 μ L of  Buffer EB **Qiagen Catalog #19086** 5m
to each well of beads, pipette mixing each well thoroughly. Incubate  00:05:00 at  Room temperature

18.5 Place back on the magnetic rack for  00:01:00 until liquid is clear again. 1m

18.6 Roll back the foil on the final cleaned reactions plate row by row. Remove  28 μ L clear eluate from the bead-cleanup plate, and place in the appropriate wells of the final cleaned reactions plate. Immediately cover the wells containing the cleaned PCR product with 8-strip caps.

18.7 Uncover the next row of samples for cleaning and  until all rows are cleaned. 

Note

Safe stopping point. Samples can be stored at 4C after this step.

Prepare Indexing PCR

19 Work in a pre-PCR area, preferably a cleaned and UV-sterilized hood to prepare your indexing reactions before going into the post-PCR area to add the PCR1 products.


20 Thaw IDT for Illumina UDI-UMI indexes on ice. Thaw only as many plates as you will need.

21 Create an indexing plate map. When using IDT for Illumina UDI-UMI indexes, fill by numbered column, not by lettered row. Make note of the indexing plate (is it plate A or B or C or D?), and the lot number. Mark on a 

	1	2	3	4	5
A	S8	S16	S24		
B	S7	S15	S23		
C	S6	S14	S22	MiFishPCRneg	
D	S5	S13	S21	Crust16S-PCR neg	
E	S4	S12	S20	S28	
F	S3	S11	S19	S27	
G	S2	S10	S18	S26	
H	S1	S9	S17	S25	
	6	7	8	9	10
A					

B					
C					
D					
E					
F					
G					
H					

	11	12
A		
B		
C		



Example plate map with 28 samples and 2 primers. Can be filled from A to H or H to A depending on your preference.

22 Indexing PCR Mastermix Recipe:

6 μL 2x Kapa HiFi Hotstart Readymix **Kapa Biosystems Catalog #KK2602**

1.1 μL

Nuclease-free water **Integrated DNA Technologies, Inc. (IDT) Catalog #11-05-01-14**

per sample.

Multiply by number of wells *10% as explained above, to create master mix.

23 In a new, clean 96-well plate (UV before use if possible and prepare in a pre-PCR space):

Add 7.1 μL Indexing Mastermix to each well that will be used. You can choose to create a library negative (good practice) at this point if you have enough empty wells for one additional one.

24 Gently vortex thawed plate of UDI-UMI indexes.

Using a 10uL filter tip set to 2.4uL, press plunger to release air bubble, then pierce the foil of the indexing plate well corresponding to the well of the plate map that you intend to index (using a multtip pipette is recommended). Collect the  2.4 μL of indexing primers and place in their corresponding well in the indexing reaction plate.

25 Take the prepared indexing plate to the post-PCR space to add the cleaned PCR products.

Add  2.5 μL of cleaned PCR product to the indexing reactions according to the indexing plate map. Cap with 8-strip tube caps, gently vortex to mix, then briefly centrifuge to get any liquid off of the caps or sides of the reaction wells.

If you are able to do a library negative, add  2.5 μL of the water used in your PCR lab to clean and elute the samples to add as a library negative control sample.

Indexing PCR Conditions

26  95 °C  00:03:00

4m 50s

8 cycles of:

 98 °C  00:00:20

 65 °C  00:00:15

 72 °C  00:00:15

final extension of:

 72 °C  00:01:00

then hold  4 °C

Gel to check Indexing PCR Products

27 Visualize PCR products in a 1.7-2% gel. Bands should be around 350-400bp.

Combine and Clean all indexed samples from each plate

- 28 Combine 10uL of up to 70 indexed samples (library) into a single 1.5mL (or 1.7mL) tube. If there are more than 70 samples, you will need another tube.
- 29 Multiply the volume of the pooled libraries in each tube by 0.9 to get the volume of Ampure XP beads needed to clean up the reactions.

For 70 uL of combined libraries you will need 63uL of beads for a total of 133uL of beads+library pool.

Perform a 0.9x bead cleanup with Ampure XP beads

28m

- 30 In the 1.5mL tube of pooled libraries, add 0.9x sample volume of Ampure XP beads and pipette mix well. incubate  Room temperature for  00:10:00

10m

- 31 Make enough fresh 80% EtOH to have 3x the total volume of the beads+library pools plus a bit extra.

- 32 Place 1.5mL tube into a magnetic rack

5m

Equipment

Magnetic Rack for for 1.5 mL Tubes

Magnetic Rack for DNA, RNA Purification; for 1.5 mL centrifuge Tubes

Sergi Lab Supplies

BOBZWXZM22

<https://www.amazon.com/Magnetic-Rack-Purification-centrifuge-Tubes/dp/BOBZWXZI20&linkCode=df0&hvadid=652498086131&hvpos=&hvnetw=g&hvrnd=671603404284>

and incubate  Room temperature for  00:05:00

- 33 Discard liquid and add an equal or greater volume of 80% EtOH. Incubate  Room temperature for  00:01:00

1m

34 Repeat the ethanol wash two more times  , then after the third 80% EtOH wash, remove all EtOH and dry the beads slightly (just until no longer wet-looking but not cracking either).

35 Resuspend beads with  100 μ L  Buffer EB Qiagen Catalog #19086 by pipette mixing thoroughly. Incubate  Room temperature  00:10:00

10m

36 Place 1.5 mL tube back on magnet rack and wait until liquid is clear, approximately  00:02:00

2m

37 remove 100uL of the clear eluate from the tube with beads while on the magnet and place in a new 1.5mL tube.

38 If you had more than one 1.5mL tube, combine equal volumes of the resulting cleaned pooled libraries (the clear eluates) into a new tube. For example, combine 20uL from each cleaned pooled library tube into a new tube.

39 Quantify with Qubit Broad range and visualize in a gel, then send an aliquot of the cleaned, pooled libraries for sequencing. Check with the sequencing core you're working with for their minimum concentration and volume and try to exceed it by a decent (10-50%) margin if possible.

Optional Quantification and Visualization

40 If you have a Qubit available, quantify with the Broad Range chemistry so that you know that the aliquot you send for sequencing is at least the sequencing core's minimum concentration. If it's below that concentration,  [go to step #30](#) starting with your cleaned pool, and elute in  30 μ L instead of 100uL in step 33.

Send for Illumina Sequencing

28m

41 Send 100uL of cleaned, pooled libraries for sequencing on any Illumina platform. (This unique dual indexing strategy compatible with any Illumina sequencer, but for more than 1.5 plates of diet samples, NovaSeq will be more likely to provide sufficient sequencing coverage than MiSeq).

Protocol references

Berry, Tina E., Sylvia K. Osterrieder, Dáithí C. Murray, Megan L. Coghlan, Anthony J. Richardson, Alicia K. Grealy, Michael Stat, Lars Bejder, and Michael Bunce. 2017. "DNA Metabarcoding for Diet Analysis and Biodiversity: A Case Study Using the Endangered Australian Sea Lion (*Neophoca cinerea*)." *Ecology and Evolution* 7 (14): 5435–53. <https://doi.org/10.1002/ece3.3123>.

Miya, M., Y. Sato, T. Fukunaga, T. Sado, J. Y. Poulsen, K. Sato, T. Minamoto, et al. 2015. "MiFish, a Set of Universal PCR Primers for Metabarcoding Environmental DNA from Fishes: Detection of More than 230 Subtropical Marine Species." *Royal Society Open Science* 2 (7): 150088. <https://doi.org/10.1098/rsos.150088>.

16S Metagenomic Sequencing Library Preparation." 2013. Illumina.

https://support.illumina.com/content/dam/illumina-support/documents/documentation/chemistry_documentation/16s/16s-metagenomic-library-prep-guide-15044223-b.pdf