

Aug 03, 2025

Metabolic stability assay in human, rat, dog or mouse hepatocytes

DOI

<https://dx.doi.org/10.17504/protocols.io.14egnykbqv5d/v1>

Nick Lynch^{1,2}

¹Curlew Research; ²ASAP Discovery Consortium

ASAP Discovery



Nick Lynch

Curlew Research, ASAP Discovery Consortium

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.14egnykbqv5d/v1>

Protocol Citation: Nick Lynch 2025. Metabolic stability assay in human, rat, dog or mouse hepatocytes . **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.14egnykbqv5d/v1>

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: May 25, 2025

Last Modified: August 03, 2025

Protocol Integer ID: 218899

Keywords: metabolic stability, hepatocyte, ADME, DMPK, drug discovery, understanding metabolic stability, metabolic stability, drug development, hepatocyte, liver cell, drug compound, critical aspect of drug discovery, stability, drug discovery, liver, drug, assay, identifying compound

Disclaimer

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to protocols.io is not peer reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgement, advice, diagnosis, or treatment. Any action you take or refrain from taking, using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with protocols.io, can be held

Abstract

Metabolic stability, measured in hepatocytes, refers to the rate at which a drug compound is metabolized or broken down by liver cells. This is a critical aspect of drug discovery and development, as it helps determine how long a drug will remain in the body and at what concentration it will exert its therapeutic effect. In drug development, understanding metabolic stability is crucial for identifying compounds that are stable enough to reach their target and for optimizing chemical structures to improve metabolic stability.

Safety warnings

- ⚠ Always wear appropriate PPE for this protocol
Refer to Material Safety Data Sheets for additional safety and handling information.

Summary

- 1 The metabolic stability in the hepatocytes is performed in 96-well plate format
This protocol can be used with the following hepatocytes:
 - Human
 - Rat
 - Mouse
 - Dog

Sample Preparation

- 2 The compounds are incubated in 1 μ M concentration in a shaking incubator at 37°C
- 3 The relevant hepatocytes are diluted to a concentration of 0.5×10^6 cells /mL and equilibrated at 37°C for 10 minutes

Assay

- 4 Biotransformation is initiated by addition of compound and mixing
- 5 Standard time-points are used as 0, 5, 10, 20, 40, 60 min, and the final solvent concentrations are 0.99 % methanol and 0.01 % DMSO
- 6 At each specified time-point, a sample aliquot (25 μ L) is removed from the test incubation mixture and combined into a cassette of up to 4 compounds, in 300 μ L ice cold methanol containing internal standards (Dextromethorphan, Diazepam, Midazolam, Phenacetin and Naloxone), and mixing to stop the reaction
- 6.1 The quenched samples are then centrifuged to precipitate the protein.

Analytical

- 7 The supernatants are analyzed by LC-MS/MS using generic analytical methods to measure the test parent compound remaining at each time-point.



Data Analysis

8 Clearance and half-life are calculated from the elimination rate constant derived from plotting response ratio versus time

9 For each test compound injection,

Response rate = Test peak area / Internal standard peak area

The natural log of the test compound response ratio is plotted versus time.

The -ve slope of the semilog plot gives the elimination rate constant k (min^{-1}), from which intrinsic clearance Cl_{int} ($\mu\text{L}/\text{min}/10^6$ cells) and half-life $t_{1/2}$ (min) are calculated.

E.g. $t_{1/2}$ (min) = $\text{Ln}2 / k(\text{min}^{-1})$

Protocol references

https://cdn.prod.website-files.com/659db4dde7f13f395f49157c/65f05da357b031fab6bc2e50_Concept%20Life%20Sciences_LE_HEPATOCYTE_STABILITY%20_2_.pdf

Acknowledgements

Grateful to Concept Life Sciences for supplying the original protocol summary