



Mar 12, 2023

## Metagenomic Library Prep from fecal sample lysate

DOI

[dx.doi.org/10.17504/protocols.io.rm7vzbob4vx1/v1](https://dx.doi.org/10.17504/protocols.io.rm7vzbob4vx1/v1)

Noah Noah Snyder-Mackler<sup>1</sup>

<sup>1</sup>Arizona State University



smacklabasu Snyder-Mackler

ASU

### Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.rm7vzbob4vx1/v1>

**Protocol Citation:** Noah Noah Snyder-Mackler 2023. Metagenomic Library Prep from fecal sample lysate . **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.rm7vzbob4vx1/v1>

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** In development

We are still developing and optimizing this protocol

**Created:** January 31, 2023

**Last Modified:** March 12, 2023

**Protocol Integer ID:** 76174

**Keywords:** Microbiome, Metagenome, Library, Lysate, metagenomic library prep from fecal sample lysate, metagenomic library prep from fecal sample, metagenomic library prep, fecal sample lysate, using illumina dna prep kit, illumina dna prep kit, fecal sample

## Disclaimer

dual index sequences here

**Nextera 10bp dual index sequences**

## Abstract

Metagenomic library prep from fecal sample lysate using Illumina DNA prep kit (1/2 reactions).

## Attachments



SMack\_Lab\_nxtra\_inde..

.

18KB

## Guidelines

### Abbreviations

BLT: Bead-Linked Transposomes

TB1: Tagmentation Buffer 1

MM: Master Mix

TWB: Tagmentation Wash Buffer

TSB: Tagmentation Stop Buffer



EPM: Enhanced PCR Mix

SPB: Sample Purification Beads

RSB: Resuspension Buffer



## Materials

- PCR plates and covers
- Nuclease-Free water
-  Illumina® DNA Prep (M) Tagmentation (24 Samples IPB) **Illumina, Inc. Catalog #20060060**
  - Bead-Linked Transposomes (BLT)
  - Tagmentation Buffer 1 (TB1)
  - Tagmentation Wash Buffer (TWB)
  - Tagmentation Stop Buffer (TSB)
  - Enhanced PCR Mix (EPM)
  - Sample Purification Beads (SPB)
  - Resuspension Buffer (RSB)
- Freshly prepared  80% ethanol **Fisher Scientific**
- Index Adaptors (96)
- Magnetic stand

## Troubleshooting





## Preparation


- 1 Add  18  $\mu\text{L}$  of nuclease-free water and  2  $\mu\text{L}$  of sample to a PCR plate

## Tagmentation Reaction



15m







- 2 Take out TB1 and keep on  On ice . Turn on thermocycler to  37  $^{\circ}\text{C}$
- 3 Multiply by number of samples, for a 96 well plate use x100 and pipet mix together to create MM. Vortex BLT vigorously for 10 seconds to resuspend

	Reagent	x1	x100
	BLT	5 ul	500ul
	TB1	5 ul	500ul

- 4 Vortex and add  10  $\mu\text{L}$  of reaction MM to each well

### Note

We found this easiest using a repeater with a  0.5 mL tip for 50 aliquots of  10  $\mu\text{L}$

- 5 Pipette mix slowly 5 times using multichannel (should have  30  $\mu\text{L}$  of volume in each well
- 6 Cover and seal plate. Then run "TAG" on thermocycler
  - Choose the preheat lid option and set to  100  $^{\circ}\text{C}$
  - Set the reaction volume to  30  $\mu\text{L}$
  -  55  $^{\circ}\text{C}$  for  00:15:00
  - Hold at  10  $^{\circ}\text{C}$

15m



7 During TAG incubation take out TSB and place at  $37^{\circ}\text{C}$  and vortex before addition to samples

8 Add  $10\ \mu\text{L}$  of TSB to each sample well using repeater with  $0.5\ \text{mL}$  tip to stop enzymatic reaction

#### Note

Be careful because buffer is foamy and be sure to pipette mix after addition

9 Cover with thick plastic and run "PTC" protocol on thermocycler

15m

- Choose the preheat lid option and set to  $100^{\circ}\text{C}$
- Set the reaction volume to  $40\ \mu\text{L}$
- $55^{\circ}\text{C}$  for 00:15:00
- Hold at  $10^{\circ}\text{C}$

10 Place on magnetic stand until beads pellet and discard supernatant

#### Note

Take out primer plates to thaw at  $4^{\circ}\text{C}$

11 1. Wash 2x with  $100\ \mu\text{L}$  of TWB. Place on magnetic stand until beads pellet + remove and discard supernatant after each wash. After both washes add in another  $100\ \mu\text{L}$  of TWB to plate while on magnetic stand to prevent beads from drying while you prep PCR mix (don't pipette mix)

## Make PCR Master Mix

1h 10m 15s



12 Take out EPR and thaw on ice

13 1. Make master mix, multiply by number of samples (use x100 for a 96 well plate) (use a 5 ml tube for a full plate of MM)



	Reagent	x1	x100
	EPR	11 ul	1,056 ul
	Nuclease free water	11 ul	1,056 ul


14 Place on magnetic stand until beads pellet + remove and discard supernatant (TWB)















15 Add  20  $\mu\text{L}$  of MM using repeater and  5 mL attachment

16 Add in  2.5  $\mu\text{L}$  of i5 and i7 primers


17 Pipette mix, cover with thick plastic and run "BLT" 12 cycle on thermocycler (

1h 10m 15s

 01:00:00 )

- Choose the preheat lid option and set to  100 °C
-  68 °C for  00:03:00
-  98 °C for  00:03:00
- 12 cycles of:
  -  98 °C for  00:00:45
  -  62 °C for  00:00:30
  -  68 °C for  00:02:00
-  68 °C for  00:01:00
- Hold at  10 °C











#### Note

1. Take out SPB and RSB while PCR is running to thaw
2. You may need to place the RSB on the thermocycler at  37 °C to get it to thaw



## Bead Clean Up

10m



- 18 Grab a new PCR plate and add  60  $\mu\text{L}$  of nuclease-free water using a multichannel pipette and  45  $\mu\text{L}$  of SPB using a repeater with a .5 ml tip to each sample well
- 19 After PCR is complete spin down plate → place on magnetic stand → transfer  25  $\mu\text{L}$  of supernatant to new PCR plate containing water and SPB. Discard old plate.
- 20 Pipette mix 10 times and incubate for  00:05:00 at  Room temperature  
During this incubation take out a new PCR plate and add  15  $\mu\text{L}$  of SPB to each well using a repeater
- 21 After the incubation place plate on magnetic stand until beads pellet
- 22 Transfer  125  $\mu\text{L}$  to the new plate containing the  15  $\mu\text{L}$  of SPB. Discard old plate.
- 23 Pipette mix 10 times and incubate for  00:05:00 at  Room temperature
- 24 Place on magnetic stand + remove and discard supernatant

5m

5m


## Ethanol Washes

3m 30s

- 25 Keep plate on magnetic stand and prepare fresh 80% ethanol
- 26 While keeping the plate on the magnetic stand add  200  $\mu\text{L}$  of 80% ethanol without mixing
- 27 Incubate  00:00:30 + remove and discard supernatant. Repeat ethanol wash and discard supernatant.

30s

**Note**

Use a  20  $\mu\text{L}$  pipette to remove as much ethanol as possible

28 Air dry while on magnetic stand (  00:03:00 )


3m



**Note**

be sure to not let beads crack

**Elution**

2m

29 Remove beads from stand and add  32  $\mu\text{L}$  of RSB to beads

30 Pipette mix to resuspend and incubate at  Room temperature for  00:02:00

2m

31 Place on magnetic stand and transfer supernatant to new and final sturdy PCR plate.  
Seal, label, and store!