



Jun 04, 2025

🌐 MG-132 (Invitrogen, #tlrl-mg132-2) and hydroxychloroquine (Sigma, #H0915-5MG) assay

DOI

<https://dx.doi.org/10.17504/protocols.io.5qpvokrpxl4o/v1>

Jailyn Izu¹

¹NIMH



Christopher Bartley

Translational Immunopsychiatry Unit, NIMH

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.5qpvokrpxl4o/v1>

Protocol Citation: Jailyn Izu 2025. MG-132 (Invitrogen, #tlrl-mg132-2) and hydroxychloroquine (Sigma, #H0915-5MG) assay. protocols.io <https://dx.doi.org/10.17504/protocols.io.5qpvokrpxl4o/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working



Created: May 15, 2024

Last Modified: March 26, 2026

Protocol Integer ID: 99888

Keywords: role for various protein degradation pathway, various protein degradation pathway, tbk1 n455, hydroxychloroquine, tbk1, mg132, 5mg

Funders Acknowledgements:

National Institute of Mental Health

Grant ID: ZIA MH002982

National Institute of Mental Health

Grant ID: R01MH122471

National Institute of Neurological Disorders and Stroke

Grant ID: U01NS120836

Abstract

Protocol for evaluating the role for various protein degradation pathways on the half-life of TBK1 N455S

Materials

- COS-1 cells
- Cycloheximide (Thermo Fisher Scientific, Catalog #AC357420010)
- Dimethyl Sulfoxide (DMSO), filtered ()
- RNase free H₂O, filter (Quality biological, #351-029-721)
- M-PER Mammalian Protein Extraction Reagent (Thermo Fisher Scientific, Catalog #78505)
- Halt Protease Inhibitor Cocktail, EDTA-free (100X) (Thermo Fisher Scientific, Catalog #78425)
- MG-132 (Invitrogen, #tlrl-mg132-2)
- Hydroxychloroquine (Sigma, #H0915-5MG)



Create stock solutions of MG-132 and HCQ

- 1 Obtain  5 mg MG-132 from Invitrogen (#tlrl-mg132-2)
- 2 Reconstitute MG-132 with  250 μL DMSO (filtered) to a final concentration of  42 millimolar (mM) .
- 3 Create 25  10 μL aliquots and store in  $-20\text{ }^{\circ}\text{C}$. These are your stock solutions which are stable for up to 3 months. You will dilute these aliquots for the experiment.
- 4 Obtain  5 mg hydroxychloroquine sulfate (HCQ) from Sigma (#H0915-5MG)
- 5 Reconstitute HCQ with  250 μL RNase free H_2O (filtered) to a final concentration of  46 millimolar (mM) .
- 6 Create 25  10 μL aliquots and store in  $-20\text{ }^{\circ}\text{C}$. These are your stock solutions. You will dilute these aliquots for the experiment.

Prepare MG-132 and HCQ dilutions for the experiment

- 7 For MG-132, create a 1:21 dilution of 42 mM stock solution in media to obtain a  2 millimolar (mM) dilution. You will add a volume of this to the cells at each time point.

*Perform dilutions needed for the experiment the **day of**, as dilutions of the stock solution should not be stored for more than 24 hours.*
- 8 For HCQ, create a 1:21 dilution of 46 mM stock solution in H_2O to obtain a  5 millimolar (mM) dilution. You will add a volume of this to the cells at each time point.

*Perform dilutions needed for the experiment the **day of**, as dilutions of the stock solution should not be stored for more than 24 hours.*

Perform cycloheximide assay



- 9 Follow the protocol for the cycloheximide assay, treating with CHX + MG-132 OR HCQ at time points (READ PROTOCOL-SPECIFIC PARAMETERS BELOW BEFORE BEGINNING).

Protocol

NAME

Cycloheximide assay (Thermo Fisher Scientific, Catalog #AC357420010)

CREATED BY

Christopher Bartley

Preview

Protocol-specific parameters:

There will be 3 conditions: (1) CHX alone, (2) CHX + MG132, and (3) CHX + HCQ). Therefore,

- For Step 4 you will,

(1) Add  20 μ L of [M] 10 mg/mL CHX (for a final concentration in the wells of 100 ug/ml CHX) +  20 μ L H₂O (filtered) to the 24-hour time point, CHX alone wells only.

(2) add  20 μ L of [M] 10 mg/mL CHX +  20 μ L of [M] 2 millimolar (mM) MG-132 to the 24-hour time point, CHX + MG-132 wells only.

(3) add  20 μ L of [M] 10 mg/mL CHX +  20 μ L of [M] 5 millimolar (mM) HCQ to the 24-hour time point, CHX + HCQ wells only.

- For Step 5 you will,

(1) Add  20 μ L of [M] 10 mg/mL CHX (for a final concentration of [M] 10 ug/ml CHX) +  20 μ L H₂O to the 16-hour time point, CHX alone wells only.

(2) add  20 μ L of [M] 10 mg/mL CHX +  20 μ L of [M] 2 millimolar (mM) MG-132 to the 16-hour time point, CHX + MG-132 wells only.



(3) add $20\ \mu\text{L}$ of $10\ \text{mg/mL}$ CHX + $20\ \mu\text{L}$ of $5\ \text{millimolar (mM)}$ HCQ to the 16-hour time point, CHX + HCQ wells only.