

Jun 20, 2023

## MGH Harvard SenNet Processing murine lung for paired single cell RNA-seq and mass spec

DOI

[dx.doi.org/10.17504/protocols.io.e6nvwd6e9lmk/v1](https://dx.doi.org/10.17504/protocols.io.e6nvwd6e9lmk/v1)

gshipkovenska Shipkovenska<sup>1,2</sup>

<sup>1</sup>Massachusetts General Hospital; <sup>2</sup>Harvard Medical School

Cellular Senescence Net...



gshipkovenska Shipkovenska

MGH

### Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.e6nvwd6e9lmk/v1>

**Protocol Citation:** gshipkovenska Shipkovenska 2023. MGH Harvard SenNet Processing murine lung for paired single cell RNA-seq and mass spec. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.e6nvwd6e9lmk/v1>

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



**Protocol status:** Working

**We use this protocol and it's working**

**Created:** June 20, 2023

**Last Modified:** June 20, 2023

**Protocol Integer ID:** 83738

**Keywords:** single cell suspension of murine lung, mgh harvard sennet processing murine lung, paired single cell rna, single cell rna, murine lung, obtaining single cell suspension, single cell suspension

**Funders Acknowledgements:**

NIH

Grant ID: GRANT13844270

## Abstract

Protocol for obtaining single cell suspension of murine lung.



## Materials

PBS - Phosphate-Buffered Saline (10X)\_pH 7.4, RNase-free (Thermo Fisher Scientific; cat. no: AM9625)

1x PBS Buffer for washes and cell suspension

5ml Phosphate-Buffered Saline (10X) pH 7.4, RNase-free (Thermo Fisher Scientific; cat. no: AM9625).

45ml Invitrogen water

0.5% BSA in 1x PBS solution

Dissolve 0.250mg BSA in 1xPBS solution above.

DNase I (50,000x) 50,000 U/ml

Liberase Stock Solution (2.5mg/ml, 13 Wunsch units/ml)

To one 5mg vial, add 2ml cold Invitrogen ultrapure water. Mix on ice by periodically stirring. Can rock on rocking platform at 4°C for maximum 30min. Aliquot in 50-100ul aliquots (~40) and freeze at -20°C. Do not freeze thaw.

Enzyme master-mix

		Origin concentration	Volume for 4ml (per lung)
	Liberase TM (Roche) (Cat# 05 401 119 001)	2,500ug/mL	160ul
	DnaseI	50,000U/ml	5ul
	RPMI 1640 (has Mg and Ca)		3.640 mL
	FBS		200ul

RBC Lysis buffer (commercial, Thermo Fisher Scientific, 00-4333-57)f

100% FBS

DMSO

## Troubleshooting

## Dissection

- 1 Euthanize mice by CO2 method.
- 2 Perfuse lungs via right ventricle.
- 3 Dissect out trachea and lungs and place in HBSS on ice.

## Single cell suspension

- 4
- 5 Pipet 1ml digestion solution on the lid of a 10cm dish. Place the lung in it. Using a razor blade, "mince" the lung, holding one tip with the forceps.
- 6 Pipet the sample into 4ml of the digestion solution / lung (8ml per mouse) and incubate at 37°C for 30min with rocking (in a 15ml falcon tube).
- 7 Use a Pasteur pipet to triturate suspension until tissue is fully dissociated.
- 8 Pass the suspension through a 100  $\mu$ m cell strainer into 50-mL tubes and use 10 mL of DPBS to pass through the strainer to wash it.
- 9 Centrifuge at 1000  $\times$  g for 5 min at 4°C, discard the supernatant.
- 10 Add 1-2 mL of 1 $\times$  RBC lysis buffer to resuspend the pellet and incubate for 5 min at 20°C–26°C.
- 11 Add 10 mL of DPBS to stop the lysing,
- 12 Centrifuge at 1000  $\times$  g for 5 min at 4°C, discard the supernatant.



- 13 Resuspend with 200  $\mu$ L of 0.5% BSA in PBS. Keep on ice.

## Count cells with heamocytometer


- 14 Count two samples each per lung and note down the exact number of cells. This information is important for scRNA library prep. We routinely recover 8-10 million cells / mouse with this protocol.

## Prepare samples for scRNA-seq and scMS

- 15 Split lung samples into 100,000 cell aliquots, according to the cell counts determined above.
- 16 Keep one aliquot on ice for the scRNA-seq sample.
- 17 For the remaining aliquots, spin down at  $1000 \times g$  for 5 min, discard the supernatant and resuspend in 90% FBS, 10% DMSO and freeze at  $-80^{\circ}\text{C}$ . These are the scMS samples.

## scRNA library preparation

- 18 Prepare according to manufacturer's protocol:

 CG000204\_ChromiumNextGEMSin...

## Protocol references

[doi: 10.1084/jem.20190865](https://doi.org/10.1084/jem.20190865)