



Apr 22, 2025

Microbial DNA extraction from the brown seaweed *Ascophyllum nodosum*

 Molecular Ecology Resources

DOI

<https://dx.doi.org/10.17504/protocols.io.8epv52kx6v1b/v1>

Coralie Rousseau¹, Sylvie Rousvoal¹, Catherine Leblanc¹, Simon M. Dittami¹

¹Biology of Marine Models (LBI2M, UMR 8227), Station Biologique de Roscoff, Sorbonne Université, CNRS



Coralie Rousseau

Station Biologique de Roscoff



Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.8epv52kx6v1b/v1>

Protocol Citation: Coralie Rousseau, Sylvie Rousvoal, Catherine Leblanc, Simon M. Dittami 2025. Microbial DNA extraction from the brown seaweed *Ascophyllum nodosum*. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.8epv52kx6v1b/v1>

Manuscript citation:

Rousseau C., Henry N., Rousvoal S., Tanguy G., Legeay E., Leblanc C., Dittami S.M. Comparison of short and long-read metabarcoding sequencing: challenge and solutions for plastid read removal and microbial community exploration of seaweed samples, (2024) 10.22541/au.172675397.78191661/v1



License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: February 11, 2025

Last Modified: April 22, 2025

Protocol Integer ID: 119953

Keywords: DNA extraction, Brown macroalgae, *Ascophyllum nodosum*, Metabarcoding, Polysaccharides, Metagenomics, High molecular weight, Microbial DNA extraction, PCR inhibitors, brown seaweed *ascophyllum nodosum*, brown alga *ascophyllum nodosum*, microbial dna extraction, optimizing microbial dna extraction, other brown seaweed, different dna purification step, polysaccharide contaminant, extraction, purification, dna

Abstract

This protocol is the result of optimizing microbial DNA extraction from the brown alga *Ascophyllum nodosum*. This protocol includes three different DNA purification steps to remove polysaccharide contaminants and other potential PCR inhibitors. This protocol takes 1.5 days for 22 samples

- ✓ Short- and long-read sequencing
- ✓ High molecular weight is obtained (>40 kb)
- ✓ Successfully tested for metagenome sequencing (skip step 1, carry out step 2 in liquid nitrogen, and use 10 mg wet weight of sample in triplicate, and skip the bead beater)

This protocol has been successfully applied to other brown seaweeds:

- ✓ *Sargassum*

Materials

Consumables & Chemicals

Reagents / Buffers:

- Acetone (100%)
- CTAB 10% (Cetyltrimethylammonium bromide)
- PVPP (Polyvinylpolypyrrolidone)
- NaCl 5M (Sodium Chloride)
- Tris HCL pH 8
- EDTA pH 8
- DIECA 0.1M (Sodium Diethyldithiocarbamate, ref. 228680, SigmaAldrich)
- β -mercaptoethanol (ref. M31148, SigmaAldrich)
- Chloroform/isoamyl alcohol (24:1)
- Ethanol 100%
- Ethanol 70% (with Biological grade water)
- Sodium acetate solution 3M
- Biological grade water (ref. 46-000-CV, Corning, USA)
- NucleoSpin Plant II kit (MachereyNagel, Germany)
- Agencourt AMPure XP beads (Beckman Coulter, Brea, CA, USA)

Disposable materials:

- 1.5 mL standart Eppendorf tubes (standard)
- 2 mL standart Eppendorf tubes (standard)
- 2 mL Safe-Lock Eppendorf tubes (ref. 0030120094)
- Parafilm
- Steel beads for TissueLyser

Lab Equipments:

- Freeze-dryer (Iyophilizer) (Qiagen, Hilden, Germany)
- Mortar (for manual grinding)
- Qiagen TissueLyser II (bead beater)
- Vortex mixer
- Centrifuges
- Incubator
- DynaMagTM-2 Magnet (ref. 12321D, ThermoFisher)
- Manual pipettes
- Dry ice
- -80°C freezer
- -20°C freezer



Sample preparation

1m

- 1 Dry matter in freeze-dryer at least  48:00:00
- 2
 - In sterile condition, grind sample by hand with sterile mortar and pestle until a very fine powder is obtained

Note

Clean with **ethanol 70%** mortar and pestle to avoid contaminations between samples

- Take 10 mg of ground algal tissue, put into 2 mL Eppendorf tube with 1 steel bead
- Put in Qiagen TissueLyser II bead beater ( 00:01:00 , 30 Hz)

Sodium diethyldithiocarbamate (CTAB) Buffer preparation (for 22 samples)

10m

- 3 Prepare in a 50 mL Falcon tube:



	Reagent	Initial concentration	Final concentration	Volume to put
	CTAB 10%	10	2	4.6 mL
	PVPP (g)	100	2	0.46 g
	Tris-HCl pH 8	1000	100	2.3 mL
	NaCl 5M	5	1.4	6.44 mL
	EDTA pH 8	500	20	0.92 mL
	Sterile H2O			8.74 mL
				Total: 23 mL

Mix the preparation

Note

Keep CTAB buffer at 21°C or higher to prevent CTAB from precipitating

Use fresh buffer each extraction



Polyphenol Cleaning

16m

- 4
- Add  1 mL **acetone (100%)**
 - Vortex, 10 min
 -  12800 x g, 21°C, 00:01:00
 - Remove acetone by pipetting
 - Repeat the previous steps
- 5 Dry between  00:05:00 to  00:10:00 to remove any remaining acetone at  Room temperature

1m



15m

First DNA purification - Remove polysaccharides

3h 30m

- 6 Remove bead from each tube by reversing the tube to make it roll away (bead added step 2)
- 7 In 2mL Eppendorf Safe-Lock Tubes:
- add  1 mL of **CTAB buffer** to each sample
 - Add  35 μ L of **0.1M DIECA** (Sodium Diethyldithiocarbamate) +  10 μ L of **β -mercaptoethanol**
 - Incubation during  01:00:00 at  65 °C
 - Under Sorbonne, add  1 mL of **chloroform/isoamyl (24:1, v/v)**

1h 30m



Note

Add parafilm around each Eppendorf Safe-Lock Tube to avoid chloroform vapours

- Vortex, 10 min
 -  13200 rpm, 21°C, 00:30:00
 - Collect the **upper phase** ( 750 μ L) in a new 2mL Eppendorf Safe-Lock Tube. Be careful to not collect the interphase
- 8 To the upper phase:
- Add  225 μ L of **ethanol 100%** drop by drop and mix by pipette agitation inside the liquid
 - Add  975 μ L of **chloroform (100%)**

20m



**Note**

Add parafilm around each Eppendorf Safe-Lock Tube to avoid chloroform vapours

- Vortex, 1 min
-  13200 rpm, 21°C, 00:20:00
- Collect the **upper phase** ( 550 µL) in new 2mL Eppendorf Safe-Lock Tubes. Be careful to not collect the interphase

9 To the upper phase:

- Add  110 µL of **acetate sodium 3M** solution
- Add  1.1 mL of **ethanol 100%** to precipitate DNA (mixing by tube inversion)
- Put samples in  -80 °C during  01:00:00

1h 30m

**Note**

Prepare  50 mL of **ethanol 70%** and put the tube on ice

-  13200 rpm, 4°C, 00:30:00
- Carefully remove liquid phase without disturbing DNA pellet

10

- Add  500 µL of **cold ethanol 70%**
- Vortex, 20s
-  13200 rpm, 4°C, 00:10:00
- Carefully remove liquid phase without disturbing DNA pellet
- Repeat step 10 once

10m

**Note**

Pellet could be translucent, white or yellow

11

- Spin down to collect remaining ethanol
- Remove the upper phase by pipetting
- Dry pellet during  00:05:00 at  Room temperature to remove remaining ethanol



- 12 Add  50 μL of **biological grade water** and mix by pipetting 10 times

Note

STOPPING POINT  -20 °C

Second DNA purification (with Nucleospin Plant Kit II)

13m

13

Note

All of the following steps are processed in sterile conditions and performed according to the manufacturer's instructions (except for step 17)

Place **Buffer PE** at  65 °C

- 14
- Add  400 μL of **biological grade water** to the previously extracted DNA and mix by pipetting 10 times
 - Add  400 μL in purple Nucleospin column
 -  11000 x g, 21°C, 00:02:00
 - Keep the flow-through which contains the DNA
 - Add  450 μL of **Buffer PC**
 - Mix by pipetting 5 times

2m



- 15
- Add  450 μL in green NucleoSpin Column
 -  11000 x g, 21°C, 00:01:00
 - Discard the flow-through.
 - Add the remainder in green NucleoSpin Column
 -  11000 x g, 21°C, 00:01:00
 - Discard the flow-through.

2m



- 16
- Add  400 μL of **buffer PW1** to the green column
 -  11000 x g, 21°C, 00:01:00
 - Discard the flow-through.

1m





17 ■ Add  700 μL of **buffer PW2** to the green column

2m

■  11000 x g, 21°C, 00:01:00



■ Discard the flow-through

■ Repeat the previous step

■ Add  200 μL of **buffer PW2** to the green column

■  11000 x g, 21°C, 00:01:00

■ Discard the flow-through

18 ■ Put the green column on a clean 1.5 mL Eppendorf

6m

■ Elution process: Add  25 μL of **Buffer PE** on the green column membrane



■ Incubation at  65 °C during  00:05:00

■  11000 x g, 21°C, 00:01:00

■ Keep the flow-through which contains the DNA

■ Repeat the elution process

Note

STOPPING POINT  -20 °C

Thrid DNA purification (with Protocol Agencourt AMPure XP)

40m 30s

19

30m

Note

All the following steps are processed in sterile conditions and performed according to the manufacturer's instructions

1. Place the bottle of **Agencourt AMPure XP** at  Room temperature  00:30:00 before use and shake it to resuspend any magnetic particles.

2. Prepare fresh 70% ethanol with biological grade water



- 20
- Add  90 μL of **Agencourt AMPure XP** (according to sample reaction volumes:
Volume AMPure = 1.8 x sample volume)
 - Mix by pipetting 10 times. Color must be homogenous
 - Incubation at  Room temperature during  00:05:00
- 21 On DynaMagTM-2 Magnet:
- Put Eppendorf on magnetic plate during  00:02:00 . Wait for the solution to clear
 - Removing slowly and discard the upper phase. Leave a few μL of supernatant behind, otherwise beads are drawn out with the supernatant.
- Add  200 μL of **ethanol 70%**
 - Incubation  00:00:30 at  Room temperature
 - Slowly remove the supernatant by pipetting and discard
 - Repeat the previous steps
- 22 Remove the reaction plate from the magnet plate, dry the pellet to remove any trace of ethanol
Wait until the pellet is a little cracked
- 23
- Add  50 μL of elution **buffer PE** (from kit NucleoSpin)
 - Mix by pipetting 10 times
 - Incubation during  00:02:00 at  Room temperature
- 24 On DynaMagTM-2 Magnet:
- Put Eppendorf on magnetic plate during  00:01:00
 - Transfert the DNA elution (upper phase) into a new 1.5 mL Eppendorf

5m



2m 30s

2m



1m

Note

The end.

Short-reads: Conservation at  -20 °C

Long-reads: Conservation at  -80 °C

Quality control

- 26 DNA must be measured by Qubit Fluorometer (Thermo Fisher Scientific) or equivalent



Protocol references

Panova, M. et al. 2016. DNA Extraction Protocols for Whole-Genome Sequencing in Marine Organisms, pp. 13–44. *In* Bourlat, S.J., ed, Marine Genomics: Methods and Protocols. Springer, New York, NY.

Varela-Alvarez, E., Andreakis, N., Lago-Leston, A., Pearson, G.A., Serrao, E.A., Procaccini, G., Duarte, C.M., and Marba, N. 2006. Genomic data isolation from green and brown brown algae (Caulerpales and Fucales) for microsatellite library construction. *J. Phycol.* 42, 741–745.

Acknowledgements

We would like to thank Gwenn Tanguy (Genomer Platform, FR2424, Station Biologique de Roscoff, Sorbonne Université, CNRS) and Erwan Legeay (Genomer Platform, FR2424, Station Biologique de Roscoff, Sorbonne Université, CNRS, 29680 Roscoff, France & Adaptation and Diversity in the Marine Environment (UMR 7144), Station Biologique de Roscoff, Sorbonne Université, CNRS) for reviewing this protocol.