



Sep 29, 2025

Version 1

## Microglia differentiation V.1

DOI

<https://dx.doi.org/10.17504/protocols.io.4r3l22zbjl1y/v1>

Riana Lo Bu<sup>1,2</sup>, Frank Soldner<sup>1,2</sup>

<sup>1</sup>Albert Einstein College of Medicine, 1301 Morris Park Ave., Bronx, NY 10461, USA.;

<sup>2</sup>Aligning Science Across Parkinson's (ASAP) Collaborative Research Network, Chevy Chase, MD, 20815.



**Oriol Busquets**

Albert Einstein College of Medicine

### Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



**DOI:** <https://dx.doi.org/10.17504/protocols.io.4r3l22zbjl1y/v1>

**Protocol Citation:** Riana Lo Bu, Frank Soldner 2025. Microglia differentiation. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.4r3l22zbjl1y/v1>

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

We use this protocol and it's working



**Created:** January 08, 2024

**Last Modified:** September 29, 2025

**Protocol Integer ID:** 93091

**Keywords:** ASAPCRN, microglia differentiation this protocol, differentiate microglia cells from hesc, differentiate microglia cell, microglia differentiation

**Funders Acknowledgements:**

**Aligning Science Across Parkinson's (ASAP)**

Grant ID: ASAP-000486

**Aligning Science Across Parkinson's (ASAP)**

Grant ID: ASAP-024409

## Abstract

This protocol has been refined to differentiate microglia cells from hESC adapted to feeder free culture systems as described in the following protocol:

### Citation

Douvaras P, Sun B, Wang M, Kruglikov I, Lallois G, Zimmer M, Terrenoire C, Zhang B, Gandy S, Schadt E, Freytes DO, Noggle S, Fossati V  
(2017). Directed Differentiation of Human Pluripotent Stem Cells to Microglia.

<https://doi.org/10.1016/j.stemcr.2017.04.023>

LINK

## Protocol Overview

- A. Flask Preparation
- B. Poly-D-Lysine (PDL) plate coating
- C. Media recipes
- D. Cell passage and differentiation
- E. Microglia precursor isolation and plating for maturation

## General Notes

A list of reagents and relevant vendor information can be found in the table listed under the materials tab.



## Attachments



Microglia differenti...

61KB

## Materials

### Reagent table:

|  | Item                           | Vendor         | Catalog Number |
|--|--------------------------------|----------------|----------------|
|  | DMEM/F12                       | Gibco          | 11320082       |
|  | Reduced growth factor matrigel | Fisher         | CB40230        |
|  | T75 Flask                      | Fisher         | 07202000       |
|  | mTESR                          | Stem Cell Tech | 100-0276       |
|  | STEMPRO-34                     | Gibco          | 10639011       |
|  | BMP4                           | PeProtech      | 120-05         |
|  | VEGF                           | PeProtech      | 100-20         |
|  | FGF                            | PeProtech      | 100-18B        |
|  | SCF                            | PeProtech      | 300-07         |
|  | IL-3                           | PeProtech      | 200-03         |
|  | TPO                            | PeProtech      | 300-18         |
|  | M-CSF                          | PeProtech      | 300-25         |
|  | FLT-3                          | PeProtech      | 300-19         |
|  | GM-CSF                         | PeProtech      | 300-03         |
|  | FLT-3                          | PeProtech      | 300-19         |
|  | M-CSF                          | PeProtech      | 300-25         |
|  | TGF- $\beta$ 1                 | PeProtech      | 100-21         |
|  | IL-34                          | PeProtech      | 200-34         |
|  | M-CSF                          | PeProtech      | 300-25         |
|  | Primaria Plates - 6w           | Corning        | 353846         |
|  | Primaria Plates - 24w          | Corning        | 353847         |
|  | Poly-D-Lysine                  | Sigma          | P6407          |
|  | Neurobasal media               | Gibco          | 21103049       |
|  | N2 Neuroplex                   | Gemini         | 400-163        |
|  | GEM21 Neuroplex                | Gemini         | 400-160        |
|  | Albumax I                      | Gibco          | 11020021       |



|  | Item            | Vendor | Catalog Number |
|--|-----------------|--------|----------------|
|  | Sodium Chloride | Fisher | BP35810        |
|  | Sodium Pyruvate | Gibco  | 11360070a      |
|  | Glutamax I      | Gibco  | 35050061       |
|  | Pen Strep       | Gibco  | 15140122       |



## Flask preparation

1h

- 1 Thaw matrigel Overnight on ice inside a 4 °C fridge. Chill all flasks, pipettes and materials to be used for flask preparation under the same conditions and keep everything as cold as possible throughout the procedure (do not use a freezer as it can cause the matrigel solution to freeze while plating).
- 2 Dilute matrigel in cold DMEM/F12 as described by the manufacturer and add into the T75 flasks (~ 6 mL /flask total volume). Incubate at 37 °C for at least 01:00:00 before use. Flasks can be stored in the incubator for several weeks before use.

1h

## Poly-D-Lysine (PDL) plate coating

6h

- 3 Prepare PDL (0.5 µg/ml) in sterile water and add it to each well and incubate 06:00:00 to Overnight at 37 °C . The day after, wash each well 4 times with sterile water to remove excess PDL (at this point the plates are ready to be used for culture - do not let the wells dry at any step as it can cause PDL crystallization (cytotoxic)).

6h

## Media preparation

- 4 The below media recipes are used in the remainder of the protocol:
  - 4.1 **STEP1 media:** mTeSR media + Penicillin-Streptomycin (100U/ml)+ 80ng/ml BMP4
  - 4.2 **STEP2 media:** StemPro-34 SFM media + 2 mM GlutaMAX + 80 ng/ml VEGF, 25 ng/ml FGF + 100 ng/ml SCF
  - 4.3 **STEP3 media:** StemPro-34 SFM media + 50 ng/ml SCF + 50 ng/ml IL-3 + 5 ng/ml TPO + 50 ng/ul M-CSF + 50 ng/ul Flt3
  - 4.4 **STEP4 media:** StemPro-34 SFM media + 50 ng/ml SCF + 50 ng/ml IL-3 + 5 ng/ml TPO + 50 ng/ul M-CSF + 50 ng/ul Flt3
  - 4.5 **Microglia maturation media:** Neurobasal media + N2 Neuroplex (1x final concentration) + GEM21 Neuroplex (1x final concentration) + 20% AlbuMAX I (0.2% final concentration)



+ NaCl (5M) (50mM final concentration) + sodium pyruvate 100x (1x final concentration)  
+ glutaMAX 100x (1x final concentration) + Penicillin-Streptomycin (100U/ml) + 50 ng/ml  
TGF- $\beta$ 1 + 100 ng/ml IL-34 + 12.5 ng/ml M-CSF

## Cell passage and differentiation

5

### Note

It is important to start the differentiation from pristine, undifferentiated feeder free cultures. For more details consult: <https://doi.org/10.17504/protocols.io.b4mcqu2w>

### Note

When changing media in the flasks, avoid scraping the bottom surface of the flasks to prevent cell and coating loss.

Manually dissect undifferentiated hPSC colonies and transfer 40-60 small aggregates into a matrigel coated flask.

6 Maintain undifferentiated cells in mTeSR media until undifferentiated colonies reach an average colony size of 1 mm. Usually,  7 mL of media/flask and media change every other day works well to maintain the cells until they reach the right size.

7 STEP1 (days 0-3): Once the hPSC colonies reach the correct size, microglia differentiation is initiated (day 0) by culturing the undifferentiated cells in  6 mL of STEP1 media. STEP1 media should be replaced with fresh media on day 2.

8 STEP2 (days 4-5): Culture cells in  6 mL of STEP2 media. No media change until the next step.

9 STEP3 (days 6-13): Culture cells in  6 mL of STEP3 media. Media should be changed on day 10.

10 STEP4 (days 14-28): Culture cells in  6 mL of STEP4 media. Media should be replaced every  96:00:00 .

4d

## Microglia precursor isolation and plating for maturation

4m

11

4m

### Note

The timepoint of microglial precursor generation as indicated by the appearance of a large number of floating cells with microglial precursor typical morphology in the supernatant is variable between hPSC cell lines and starts usually around day 28 of differentiation. Microglia precursor differentiation can be validated by FACS analysis (described in another protocol in this same collection). We normally only initiate microglial differentiation from cultures with more than 80% CD11b/CD45 and CD14/CD16 double positive microglial precursors.

### Note

If properly maintained with STEP4 media changes every 4 days, cultures in STEP4 are capable of generating microglia precursors for several months (as monitored by FACS). Microglia precursors can be collected every 4 days for microglial maturation as described below.

### Note

For media changes in STEP4 on days without collecting microglial precursors for maturation (as described below), we usually try to return microglial precursors after media change. Because microglia precursor cells grow mainly in suspension, we collect ~  5 mL of microglial containing supernatant in a conical tube, centrifuged at  150 rcf for  00:04:00, and subsequently resuspend the cells in  5 mL of fresh STEP4 media to return to the same flask.

### Note

Microglial precursors can be plated on either Primaria plates or Poly-D-Lysine (PDL) coated cell culture plates or glass cover slips. However, based on our experience, plating on primaria plates results in more consistent differentiation to mature microglia.



Once microglia precursor differentiation is validated on day 28 or later, microglia progenitors are collected from STEP4 flasks by transferring the supernatant to a conical tube followed by centrifugation at  150 rcf, Room temperature for  00:04:00 .

The cell pellet is resuspended in microglia maturation media and plated at 75k cells/cm<sup>2</sup> (Primaria or PDL-coated plates) in microglia maturation media. Media should be changed every 4 days for at least 2 weeks to allow the cells to mature.

## Citations

Douvaras P, Sun B, Wang M, Kruglikov I, Lallo G, Zimmer M, Terrenoire C, Zhang B, Gandy S, Schadt E, Freytes DO, Noggle S, Fossati V. Directed Differentiation of Human Pluripotent Stem Cells to Microglia.

<https://doi.org/10.1016/j.stemcr.2017.04.023>