

Sep 02, 2024

## Microplastics extraction from oyster tissue

DOI

<https://dx.doi.org/10.17504/protocols.io.j8nlk81dxl5r/v1>

Benjamin-Rafael Mingoa<sup>1</sup>, Archana Anand<sup>1</sup>

<sup>1</sup>San Francisco State University

Anand Lab



Archana Anand

San Francisco State University

### Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



**DOI:** <https://dx.doi.org/10.17504/protocols.io.j8nlk81dxl5r/v1>

**Protocol Citation:** Benjamin-Rafael Mingoa, Archana Anand 2024. Microplastics extraction from oyster tissue . protocols.io <https://dx.doi.org/10.17504/protocols.io.j8nlk81dxl5r/v1>

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



**Protocol status:** Working

**We use this protocol and it's working**

**Created:** September 01, 2024

**Last Modified:** September 02, 2024

**Protocol Integer ID:** 106797

**Keywords:** Oyster, Microplastics, Tissue, Olympia oyster, microplastics extraction from oyster tissue, microplastics in oyster tissue, microplastics extraction, microplastic, microplastic particle, oyster tissue, size of microplastic, microplastic particles every time, plastic particle, such as detergent, detergent, many water filtration system, such as the san francisco bay, tissue, extraction, such as acrylic, nylon, polyester, acrylic

**Funders Acknowledgements:**

College of Science & Engineering, San Francisco State University

## Abstract

Microplastics are defined as any plastic particles that are smaller than 5mm. They have been found in many of our household products, such as detergents, toothpaste, and other beauty products such as exfoliants. It has also been found that clothes with synthesized fibers, such as acrylics, nylon, and polyester shed hundreds of thousands of microplastic particles every time you wear or wash them. Due to the size of microplastics, they easily pass through many water filtration systems and end up in bodies of water, such as the San Francisco Bay or the Pacific Ocean. This protocol describes the procedure to extract and quantify microplastics in oyster tissue.

## Preparation and extraction of microplastics from Oyster Gastrointestinal Tract/Tissues

- 1 Take samples of the gastrointestinal tract (GIT) /biofilm and gill tissues (fish gills have the same biofilm membrane)
- 2 Weigh out 1-2 g of the tissues. Use a glass cover for weighing.
- 3 Cut the sample with forceps to increase surface area
- 4 10% KOH, 50 mL,  
(as a guideline 1:50 ratio of tissue to KOH works well)
- 5 incubate at 60 C for 24 hours
- 6 Filter through coarse mesh
- 7 Set up a filtration system (Whatman 1 filter, flask)
- 8 Soak the Whatman 1 filter in distilled water to prevent it from floating
- 9 USE VACUUM (or samples will be compromised due to floating)
- 10 Pour the KOH and tissue solution into the Whatman 1 filter. Vacuum.
- 11 Rinse sides of filter with distilled water
- 12 Carefully take the filter and place it into a glass plate with forceps



- 13 Wash with 50 mL boiling water
- 14 Incubate for 10 min with 10 mL acetone
- 15 Wash with another 10 mL of acetone
- 16 Staining with 1mL of diluted Nile Red
- 17 Rinse with 10 mL of water (wash the sides of the filter)
- 18 Place into a Petri plate
- 19 Observe, use a fluorescent scope