

Jul 31, 2023

Microscopy-based bead protein-protein interaction assay

DOI

<https://dx.doi.org/10.17504/protocols.io.dm6gp3erpvzp/v1>

Elisabeth Holzer^{1,2}

¹Max Perutz Labs, University of Vienna, Vienna Biocenter, Vienna, Austria;

²Vienna Biocenter PhD Program, Doctoral School of the University of Vienna and Medical University of Vienna, Vienna, Austria



OLIVIA HARDING

Children's Hospital of Philadelphia

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.dm6gp3erpvzp/v1>

Protocol Citation: Elisabeth Holzer 2023. Microscopy-based bead protein-protein interaction assay. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.dm6gp3erpvzp/v1>

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this protocol and it's working

Created: June 07, 2023

Last Modified: May 31, 2024

Protocol Integer ID: 83013

Keywords: ASAPCRN, based bead protein, bead protein, proteins as prey, protein interaction, proteins as bait, tagged protein, purified protein, protein interaction in an equilibrium state, protein, assay this protocol, assay, microscopy

Funders Acknowledgements:

Aligning Science Across Parkinson's

Grant ID: : Mechanisms of mitochondrial damage control by PINK1 and Parkin (ASAP-000350)

Abstract

This protocol describes how to perform microscopy-based bead protein-protein interaction assay with GST- or mCherry-tagged proteins as baits and fluorescently-tagged proteins as preys. The protocol requires to have purified proteins and allows to monitor protein-protein interaction in an equilibrium state. The fluorescent signal can be quantified.

Guidelines

Experiment should be repeated at least three times for statistical analysis.

Materials

Glutathione Sepharose 4B (Cytiva)

RFP-Trap Agarose beads (ChromoTek)

384-well glass-bottom microplate (Greiner Bio-One)

confocal microscope

ImageJ software

Before start

SEC Buffer:

25 mM HEPES pH 7.5

150 mM NaCl

Freshly added: 1 mM DTT

Purify tagged bait and prey proteins

Protocol

NAME

EXPRESSION AND PURIFICATION OF HUMAN NEMO (GST-GFP-NEMO)

CREATED BY

OLIVIA HARDING

Preview



Prepare bait-coated beads

1h 2m

- 1 Equilibrate  20 μL Glutathione Sepharose 4B or RFP-Trap Agarose beads with  200 μL SEC buffer

Note

Prepare these bead equilibrations for each condition.

- 2 Incubate equilibrated beads with GST- or mCherry-tagged bait protein for a final concentration of  5 micromolar (μM) in SEC buffer for  01:00:00 at  4 $^{\circ}\text{C}$ with gentle rotation.

1h

- 3 Centrifuge beads at 3000 rcf for  00:02:00 at  4 $^{\circ}\text{C}$

2m

- 4 Remove the supernatant and wash beads with  200 μL SEC buffer

- 4.1 Repeat for a total of 2 washes, then discard buffer

- 5 Add  20 μL SEC buffer to achieve a beads:buffer ratio of 1:1

Interaction assay set-up

30m

- 6 Pipette prey proteins into the wells of a 384-well glass-bottom microplate (Greiner Bio-One)

Note

A volume of at least  20 μL should be used to cover the bottom of the well



Note

Prey concentration should be [M] 0.1 micromolar (μM) to [M] 1 micromolar (μM) , but should be adjusted depending on the strength of the interaction

Note

Different preys should be conjugated to different fluorophores

7 Pipette  1 μL of bait-coated beads into each well

8 Incubate the plate for  00:30:00 in the dark at  Room temperature

30m

Signal detection

9 Use a microscope configured to detect fluorescent signal (e.g. Zeiss LSM 700 confocal microscope equipped with Plan-Apochromat 20X/0.8 objective)

10 Acquire fluorescent images in the middle section of the beads and collect more than one image for each well

11 Also acquire bright field images for each field

Quantification using ImageJ

12 Draw several lines across each bead in the fluorescent channel and measure the intensity along the lines

13 Record the maximum intensity for each bead

14 For background correction, measure the average intensity of a rectangular ROI that covers an area of each field of view with no beads



- 15 Subtract the average intensity of the background ROI from each bead maximum in that field
- 16 Calculate the average of the background-corrected maximum intensities of beads for each sample