

Nov 20, 2020

Microtome Dissection of Beetles

DOI

dx.doi.org/10.17504/protocols.io.bnuqmevw



Jiri Hulcr¹

¹University of Florida

Protocols Bark Beetle M...



Bark Beetle Mycobiome Research Coordination Network Bark

Beetle

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.bnuqmevw>

Document Citation: Jiri Hulcr 2020. Microtome Dissection of Beetles. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.bnuqmevw>



License: This is an open access document distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Created: October 23, 2020

Last Modified: November 20, 2020

Document Integer ID: 43632

Keywords: microtome dissection of beetle, part of the bark beetle mycobiome, bark beetle mycobiome, microtome dissection, fungus symbiosis, beetle, symbiosis, widespread insect, fungus, more information on the bbm, mandibular mycangia, bbm

Abstract

This protocol describes how visualize internal structures, from mandibular mycangia, to beetle heads or bodies.

This protocol is part of the Bark Beetle Mycobiome (BBM) Research Coordination Network. For more information on the BBM international network: Hulcr J, Barnes I, De Beer ZW, Duong TA, Gazis R, Johnson AJ, Jusino MA, Kasson MT, Li Y, Lynch S, Mayers C, Musvuugwa T, Roets F, Selmann KC, Six D, Vanderpool D, & Villari C. 2020. Bark beetle mycobiome: collaboratively defined research priorities on a widespread insect-fungus symbiosis. Symbiosis 81: 101–113 <https://doi.org/10.1007/s13199-020-00686-9>.

Troubleshooting



Fixating

For visualizing mandibular mycangia, beetle heads or bodies should be fixed in 96% ethanol, immersed in 30% hydrogen peroxide for 24 h, and embedded in paraffin.

Microtome

Preparates can be sectioned on a microtome at 5 μm and stained with hematoxylin (to emphasize nuclei in fungi) and eosin (to visualize proteins, muscles, and other beetle tissue).

Important points

1. The immersion in hydrogen peroxide is critical for softening beetle exoskeleton and to avoid fracturing. The angle of beetle immersion in the wax is also critical. It depends on what mycangia you are interested in, what angle you want to slice them, and what is the angle of attachment of the wax sample to your microtome. You may need to play with the beetle position as the wax is hardening.