

Apr 26, 2019

Microwell-based Single-Cell RNA-seq

DOI

<https://dx.doi.org/10.17504/protocols.io.yqmfvu6>



Peter Sims¹, Jinzhou Yuan¹, Yim L. Cheng¹

¹Columbia University Irving Medical Center

Human Cell Atlas Metho...



Peter Sims

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.yqmfvu6>

Protocol Citation: Peter Sims, Jinzhou Yuan, Yim L. Cheng 2019. Microwell-based Single-Cell RNA-seq. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.yqmfvu6>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol in our group and it is working.

Created: March 01, 2019

Last Modified: April 26, 2019

Protocol Integer ID: 20973

Keywords: single-cell, RNA-seq, single cell rna, cell rna, chip reverse transcription version, rna, reverse transcription version, seq, cell

Abstract

(On-chip reverse transcription version, Sims Lab – Jinzhou Yuan and Yim L. Cheng) 7/6/2017

Attachments



Microfluidic Single-...

38KB

Guidelines

Our single-cell RNA-seq technology uses a custom microfluidic device that we produce in our lab.

The CAD drawings of these devices that one can submit to a foundry for photolithographic production of molds from which these devices can be routinely generated:

http://www.columbia.edu/~pas2182/localhost/protocols/First_Mold_JY_CU_6inchwafer_microwellarrays_50umthick.dwg

http://www.columbia.edu/~pas2182/localhost/protocols/Second_Mold_JY_CU_6inchwafer_flowcell_100umthick.dwg

Safety warnings

 Please refer to the SDS (Safety Data Sheet) for information about hazards, and to obtain advice on safety precautions.



Before start

Day before experiment:

Device Preparation

- Fill a new device with wash buffer (20mM Tris-HCl pH7.9, 50mM NaCl, 0.1% Tween 20). Cover both the inlet and outlet of the device with a puddle of wash buffer.
- Store the wash buffer filled device in a humid chamber (P1000 pipette tip container filled with water) at room temperature.

Experiment day 1

1 Live stain cells

Harvest and resuspend 500,000 cells in  1 mL TBS **(500 cells/uL)**.

2 Incubate the cells in 4 uM CalceinAM (4ul 1mM CalceinAM per ml cells) for  00:15:00 on ice.

3 Preparation of lysis buffer, and wash buffer, and fluorinated oil

While the cells are being live stained, make

 1 mL lysis buffer ( 990 µL buffer TCL +  10 µL 2-Mercaptoethanol)

 2 mL RNase inhibitor-doped wash buffer ( 1998 µL wash buffer +  2 µL SUPERaseIN).

4 Load buffers and fluorinated oil to reagent reservoirs

- reagent reservoir 1:  1 mL lysis buffer
- reagent reservoir 4:  1.8 mL fluorinated oil
- reagent reservoir 5:  2 mL RNase inhibitor-doped wash buffer

5 Prepare mRNA-capture-beads

Gently vortex the beads, transfer  500 µL of beads into a microcentrifuge tube.

6 Wash the beads with  1 mL TBS buffer and resuspend the beads in  400 µL TBS buffer . Keep the beads on ice before use.

7 mRNA-capture-bead and cell loadings

Flush  1 mL TBS through the wash buffer-filled device. Exercise caution to avoid flowing air into the device.



- 8 Load  500 μ L live-stained cells into the TBS-flushed device. Fine tune the size of liquid puddles on both inlet and outlet of the device until the cells stop moving in the device. Let cells settle for  00:03:00 .
- 9 Flush  1 mL TBS through the cell-loaded device. Flush gently so that the loaded cells don't get washed out.
- 10 Confirm 5-10% of wells contain single cell under microscope.
- 11 Load  400 μ L TBS-washed beads into the device. Reuse the beads to repeat the loading until >60% wells contain beads. Load gently so that the bead flow does not wash cells out.
- 12 Flush  1 mL TBS through the cell-bead-loaded device. Flush gently so that the loaded cells and beads don't get washed out.
- 13 Confirm >60% of wells contain mRNA-capture-bead under microscope.
- 14 **Pre-run fluidics wash**

Open channel 1 for  00:00:10 by executing the following command in the MinGW terminal: **"Valvesonly_3_dual_3.exe 10000 1"**.
- 15 Open channel 4 for  00:00:10 by executing the following command in the MinGW terminal: **"Valvesonly_3_dual_3.exe 10000 4"**.
- 16 Open channel 5 for  00:00:30 by executing the following command in the MinGW terminal: **"Valvesonly_3_dual_3.exe 30000 5"**.
- 17 **Connect the device to the fluidics system and start the run**

Plug the inflow (red) and outflow (green) tubing to the inlet and outlet of the device respectively. **Make sure that the tubings are NOT inserted all the way to the bottom of the inlet or outlet as this can clog the device**, which will prevent liquid from flowing through or flow at a significantly reduced flow rate.

- 18 Start the run by executing the following command in the MinGW terminal:
"**SingleCellIRNACapture_TCL_dual_large_array_3.exe**". Lysis buffer will be pushed through the device immediately followed by fluorinated oil, which will seal the wells, physically isolating the lysate of each individual cell.
- 19 The device will be heated to  50 °C to promote cell lysis and mRNA capture, and then cooled back to room temperature.

20 Cell lysis QC and mRNA capture

Once the device has been cooled to room temperature, disconnect the device from the fluidics system. Seal both inlet and outlet of the device with disc tapes.

- 21 Imaging the device under fluorescence microscope with a 10X objective. Check the distribution of lysate (traced by the live stain dye) in the device. You should see the fluorescent dye fill up the wells that had cells, and the absence of fluorescent dye in wells that did not have cell. If this is the case, unseal the inlet and outlet of the device, reconnect the device to the fluidics system, and proceed to the next step. Otherwise, abort the experiment.
- 22 Reverse transcription

Prepare the following reverse transcription mixture

	Final Concentration	Volume	Stock Concentration
	1X Maxima RT buffer	100 µL	5X Maxima RT buffer
	1 mM dNTPs	50 µL	10 mM dNTP
	1 U/ul SupraseIN	25 µL	20 U/ul SupraseIN
	2.5 uM SMRT_TSO	12.5 µL	100 uM SMRT_TSO
	10 U/ul Maxima H-RT	25 µL	200 U/ul Maxima H- RT
	0.1% Tween 20	5 µL	10% Tween20
	Nuclease free water	282.5 µL	
	Total	500 µL	

23 Load all  500 μ L reverse transcription mixture to reagent reservoir 2 before mRNA capture step is completed.

24 **Bead extraction**

Assemble a syringe-tubing set: 2.5 cm green PEEK tubing and 1ml syringe connected with 12 cm tygon tubing.

25 Prepare  3 mL wash buffer , keep on ice

26 “Massage” the device to dislodge the beads from the device, then flow  1 mL wash buffer through the device with syringe-tubing set, **COLLECT the overflow.**

Note

cDNA-coated beads will be washed out from the device with the buffer flow, so you want to collect the overflow.

27 Repeat step 26 until >95% of the beads have been removed from the device.

 [go to step #26](#) Repeat step 26

28 Wash the cDNA-beads with

 1 mL TE/SDS buffer (10 mM Tris pH 8.0, 1 mM EDTA, 0.5% SDS)

29 Wash the cDNA-beads with

 1 mL TE/TW buffer (10 mM Tris pH 8.0, 1 mM EDTA, 0.01% Tween-20) (1/2)

30 Wash the cDNA-beads with

 1 mL TE/TW buffer (10 mM Tris pH 8.0, 1 mM EDTA, 0.01% Tween-20) (2/2)

31 Resuspend the cDNA-beads in  1 mL TE/TW buffer .

32 Store cDNA-beads at  4 °C .

33 Wash fluidics system

Rinse and fill ALL reagent reservoirs (except reservoir 4) with DI water.

34 Execute the following command in the MinGW terminal:
"SingleCellRNACapture_wash_dual_3.exe"

Experiment day 2

35 Exo-I reaction

Prepare the following Exo-I reaction mixture

	Final Concentration	Volume	Stock Concentration
	1X Exo I buffer	5 µL	10X Exo I buffer
	1 U/uL Exo I	2.5 µL	20 U/uL
	Nuclease free water	42.5 µL	
	Total	50 µL	

36 Wash the cDNA-beads with  1 mL nuclease-free water . (1/2)

37 Wash the cDNA-beads with  1 mL nuclease-free water . (2/2)

38 Resuspend the beads in  50 µL Exo-I reaction mixture .

39 Run the following program on a thermocycler:

	Temperature	Time	Number of cycles
--	-------------	------	------------------

	37 °C	30 minutes	1
	4 °C	Hold	

40 Wash the Exol-treated cDNA-beads with  1 mL TE/SDS buffer .

41 Wash the Exol-treated cDNA-beads with  1 mL TE/TW buffer . (1/2)

42 Wash the Exol-treated cDNA-beads with  1 mL TE/TW buffer . (2/2)

43 Wash the Exol-treated cDNA-beads with  1 mL nuclease-free water . (1/2)

44 Wash the Exol-treated cDNA-beads with  1 mL nuclease-free water . (2/2)

45 **SMRT PCR reaction (cDNA PCR amplification)**

Prepare the following SMRT PCR reaction mixture

	Final Concentration	Volume	Stock Concentration
	1X PCR Ready Mix	500 µL	2X PCR Ready Mix
	1 uM SMRT PCR primer	10 µL	100 uM SMRT PCR primer
	Nuclease free water	490 µL	
	Total	1000 µL	

46 Resuspend cDNA-beads in  1000 µL SMRT PCR reaction mixture . Aliquot  50 µL per PCR tube.

47 Run the following PCR program on a thermocycler:

Temperature	Time	Number of cycles
95 °C	3 minutes	1
98 °C	20 seconds	4
65 °C	45 seconds	
72 °C	3 minutes	
98 °C	20 seconds	8
67 °C	20 seconds	
72 °C	3 minutes	
72 °C	5 minutes	1
4 °C	hold	

48 **SMRT PCR product purification with 0.6X Ampure beads**

Combine all SMRT PCR products into a single microcentrifuge tube.

49 Spin for  00:01:00 , transfer  900 µL supernatant to a new microcentrifuge tube.

50 Add  540 µL Ampure beads (0.6X volume) to tube. Incubate at room temperature for  00:08:00 .

51 Place tube on a magnetic stand for  00:05:00 . Discard the supernatant.

52 Wash Ampure beads with 80% ethanol twice

- Add  1 mL freshly-made 80% ethanol in nuclease-free water to tube, wait for  00:00:30 , remove 80% ethanol.
- Repeat 80% ethanol wash.
- Quick spin for  00:00:30 , place it back on a magnetic stand for  00:00:30 , and remove residual 80% ethanol.



- 53 Cover the tube with a KimWipe, air dry the Ampure beads on the magnetic stand for  00:05:00 .
- 54 Elute cDNA off Ampure beads with 20 ul nuclease-free water
- Add  20 µL nuclease-free water directly on the Ampure beads, remove the tube from the magnetic stand, mix by pipetting.
 - Incubate at room temperature for  00:05:00 .
 - Place tube on a magnetic stand for  00:03:00 .
 - Keep supernatant (amplified cDNA).
- 55 Store cDNA at  -20 °C .
- 56 cDNA QC with Qubit and Bioanalyzer (follow vendor's exact instructions).
- 57 **Nextera tagmentation reaction (Nextera XT kit)**
- Add  10 µL TD Buffer to a PCR tube.
- 58 Add  5 µL of 0.6 ng cDNA plus nuclease free water.
- 59 Add  5 µL ATM to the PCR tube. Mix by pipetting.
- 60 Run the following program on a thermocycler:
- | | Temperat
ure | Time | Number of
cycles |
|--|-----------------|--------------|---------------------|
| | 55 °C | 5
minutes | 1 |
| | 10 °C | hold | |
- 61 Immediately after the temperature reaches 10 °C, add  5 µL NT Buffer to the PCR tube and mix by pipetting.
- 62 Incubate at room temperature for  00:05:00 .

**63 Selective amplification of the 3' end fragment of cDNA (Nextera XT kit)**

Add the following reagents to the end product in step 62 and mix by pipetting:

Reagent	Volume
NPM	15 µL
N7 PCR primer	5 µL
2 uM Custom P5 primer (not from Nextera XT kit)	5 µL

64 Run the following program on a thermocycler:

Temperature	Time	Number of cycles
72 °C	3 minutes	1
95 °C	30 seconds	1
95 °C	10 seconds	12
55 °C	30 seconds	
72 °C	30 seconds	
72 °C	5 minutes	1
4 °C	hold	

65 Nextera PCR product purification with 0.6X + 1X Ampure bead

Quick spin the Nextera PCR product  00:01:00 , transfer  50 µL PCR product to a new microcentrifuge tube.

66 Add  30 µL Ampure beads (0.6X volume) to tube. Incubate at room temperature for  00:08:00 .

67 Place tube on a magnetic stand for  00:05:00 . Discard the supernatant.



- 68 Wash Ampure beads with 80% ethanol twice
- Add  200 μ L freshly-made 80% ethanol in nuclease-free water to tube, wait for  00:00:30 , remove 80% ethanol.
 - Repeat 80% ethanol wash.
 - Quick spin for  00:00:30 , place it back on a magnetic stand for  00:00:30 , and remove residual 80% ethanol.
- 69 Cover the tube with a KimWipe, air dry the Ampure beads on the magnetic stand for  00:05:00 .
- 70 Elute the Nextera library off Ampure beads with 50 μ L nuclease-free water
- Add  50 μ L nuclease-free water directly on the Ampure beads, remove the tube from the magnetic stand, mix the beads and the water by pipetting.
 - Incubate at room temperature for  00:05:00 .
 - Place tube on a magnetic stand for  00:03:00 .
 - Keep supernatant (amplified cDNA).
- 71 Repeat Ampure bead purification steps (steps 65-70) with 1x volume Ampure bead (50 μ L) and elute the Nextera library off Ampure beads with  20 μ L nuclease-free water .
-  [go to step #65](#)
- 72 Store the Nextera library at  -20 $^{\circ}$ C .
- 73 Nextera library QC with Qubit and Bioanalyzer (follow vendor's exact instructions).
- 74 **Sequencing**
- Use custom read 1 sequencing primer.
Use 20% phix.
21 cycles on read 1.
63 cycles on read 2.
8 cycles on index read 1 (if multiple samples are pooled together).